HEALICOIL[®] REGENESORB Suture Anchors

10601068-H

Smith & Nephew Orthopaedics GmbH
Alemannenstraße 14
78532 Tuttlingen
Germany
ec.rep@smith-nephew.com

Endoscopy
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810
USA

10/2018 10601068 Rev. H

CE mark on the package insert (IFU) is not valid unless there is a CE mark on the product (description) label.

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HEALICOIL[®] REGENESORB Suture Anchors

Device Description

Smith & Nephew HEALICOIL REGENESORB Suture Anchors consist of an absorbable suture anchor with attached non-absorbable sutures, assembled with an insertion device. The suture anchor is intended to provide secure reattachment of soft tissue to bone. Attachment of the soft tissue is performed by the surgeon's preferred technique.

This device is provided sterile, for single use only.

Contents

The suture anchors are packaged in the following configurations. Refer to individual suture anchor product labels for suture size, type and quantity.

Each HEALICOIL REGENESORB Suture Anchor contains:

- 1 ea. Absorbable, biocomposite suture anchor – poly L-lactide-co-glycolide, calcium sulfate and β -tricalcium phosphate
- 1 ea. Insertion device – stainless steel shaft with polycarbonate and ABS handle
- 2 or 3 ea. Sutures – braided, uncoated, UHMW polyethylene, blue/cobraid blue, non-absorbable and/or
UHMW polyethylene with monofilament nylon cobraid black, non-absorbable

Indications for Use

HEALICOIL REGENESORB Suture Anchors are intended for use only for the reattachment of soft tissue to bone for the following indications:

Shoulder

- Bankart lesion repairs
- SLAP lesion repairs
- Acromioclavicular separation repairs
- Rotator cuff tear repairs
- Capsular shift or capsulolabral reconstructions
- Biceps tenodesis
- Deltoid repairs

Foot and Ankle

- Hallux valgus repairs
- Medial or lateral instability repairs/reconstructions
- Achilles tendon repairs/reconstructions
- Midfoot reconstructions
- Metatarsal ligament/tendon repairs/reconstructions

Elbow

- Ulnar or radial collateral ligament reconstructions
- Lateral epicondylitis repair
- Biceps tendon reattachment

Knee

- Extra-capsular repairs
 - Medial collateral ligament
 - Lateral collateral ligament
 - Posterior oblique ligament
- Iliotibial band tenodesis
- Patellar realignment and tendon repairs
 - Vastus medialis obliquus advancement

Contraindications

- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be conducted and sensitivity ruled out prior to implantation.
- Pathological conditions of bone, such as cystic changes or severe osteopenia, which would compromise secure anchor fixation.
- Pathological conditions in the soft tissues to be attached that would impair secure fixation by suture.
- Comminuted bone surface, which would compromise secure anchor fixation.
- Physical conditions which would eliminate, or tend to eliminate, adequate anchor support or impair healing, i.e., limitation of blood supply, infection, etc.

Warnings

 Do not use if package is damaged. Do not use if the product sterilization barrier or its packaging is compromised.

- Contents are sterile unless package is opened or damaged. **DO NOT RESTERILIZE.** For single use only. Discard any open, unused product. Do not use after the expiration date.
- It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.
- Read these instructions completely prior to use.
- Product must be stored in the original sealed pouch.
- Incomplete anchor insertion may result in poor anchor performance.
- Breakage of the suture anchors can occur if insertion sites are not prepared with appropriate instrumentation prior to implantation.
- Do not attempt to implant this device within cartilage epiphyseal growth plates or non-osseous tissue.
- Do not resterilize or reuse anchors, sutures, insertion devices, or disposable hole preparation instrumentation.

Precautions

R_{ONLY} U.S. Federal law restricts this device to sale by or on the order of a physician.

- Hazards associated with reuse of this device include, but are not limited to, patient infection and/or device malfunction.
- Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.
- The suture anchor and suture are not intended to provide indefinite biomechanical integrity. As in all suture anchor or suturing techniques, until biological attachment of tissue to bone is complete the fixation should be considered temporary and may not withstand weight bearing or other unsupported stresses.
- Use of excessive force during insertion can cause failure of the suture anchor or insertion device.
- Bone quality must be adequate to allow proper placement of the suture anchor.
- Inadequate bone quality could result in loss of fixation or pullout of the suture anchor.
- Do not alter the implant or instrumentation or performance may be compromised.
- Do not use sharp instruments to manage or control the suture.
- Breakage of the suture can occur if there is damage caused by sharp instruments.
- After use, this device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.

Adverse Reactions

- Mild inflammatory reaction
- Foreign body reaction
- Infection, both deep and superficial
- Allergic reaction

Instructions for Use

1. Using routine sterile technique, remove the HEALICOIL[®] REGENESORB Suture Anchor from the package.
2. Use the appropriate Smith & Nephew reusable or disposable threaded dilator or drill to prepare the insertion site. Refer to the table below for the recommended threaded dilator or drill.

Notes:

Smith & Nephew reusable or disposable threaded dilators specific to the suture anchor size are sold separately.

When a larger hole is required, as in the case with hard bone, use a drill to create a pilot hole for the anchor, then insert the threaded dilator to better prepare the bone for the anchor.

Suture Anchor	For General Use	For Hard Bone
HEALICOIL REGENESORB 4.75 mm Suture Anchor	REF 72203709 4.75 mm Reusable Threaded Dilator or REF 72203951 4.75 mm Disposable Threaded Dilator	REF 72203482 3.5 mm Spade Tip Drill REF 72203709 4.75 mm Reusable Threaded Dilator or REF 72203951 4.75 mm Disposable Threaded Dilator
HEALICOIL REGENESORB 5.5 mm Suture Anchor	REF 72203710 5.5 mm Reusable Threaded Dilator or REF 72203952 5.5 mm Disposable Threaded Dilator	REF 72203483 4.5 mm Spade Tip Drill REF 72203710 5.5 mm Reusable Threaded Dilator or REF 72203952 5.5 mm Disposable Threaded Dilator

3. Insert the threaded dilator fully into the bone by rotating the device clockwise until the indicator line on the dilator is flush with the surface of the bone.
or
If drilling is required, place the distal tip of the drill onto the bone surface at the desired location. Advance the drill until the indicator line on the drill bit is flush with the bone surface.
 - a. Remove the drill.
 - b. Insert the threaded dilator fully into the bone by rotating the device clockwise until the indicator line on the dilator is flush with the surface of the bone.
4. Remove the threaded dilator from the insertion site by rotating it counterclockwise.
5. Establish axial alignment of the suture anchor to the insertion site. While supporting the insertion site, rotate the anchor clockwise into the prepared bone hole. Continue rotating the anchor until the indicator line on the insertion device is flush with the surface of the bone.
6. Visually inspect the site to confirm that desired placement is achieved.

Notes:

Insert the insertion device to the top of the horizontal indicator line to place the anchor approximately 1 mm below the bone surface.

Axial lines on the insertion device indicate suture orientation within the anchor.

CAUTION: Use of excessive force during insertion can cause failure of the suture anchor or insertion device.

7. If more torque is required to insert the anchor, stop and ensure that the hole size and depth are correct for the bone conditions encountered. It is the responsibility of the surgeon to determine the patient's bone condition, appropriately prepare the insertion site, and determine the suitability of the implant for the procedure.
8. Disengage the insertion device from the suture anchor:
 - a. Hold the handle securely in the palm and pull back on the finger grips with two fingers to retract the suture tensioning mechanism.
 - b. Hold the tensioning mechanism in the back position with the finger grips and slowly remove the inserter from the insertion site until the suture tips are visible. The suture releases and feeds through the insertion device as it is removed.
 - c. If necessary, tap lightly on the underside of the handle to remove the insertion device from the insertion site.

Alternatively:

- a. Pull back on the finger grips and unwind the suture with the free hand.
- b. Slowly remove the insertion device.
- c. As the insertion device is removed, pinch the suture below the insertion device; use the free hand to verify suture release from the inserter.

CAUTION: Do not use sharp instruments to manage or control the suture.

9. Discard the insertion device.

Warranty

For single use only. This product is warranted to be free from defects in material and workmanship. Do not reuse.

For Further Information

If further information on this product is needed, contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized representative.

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©2014, 2018 Smith & Nephew, Inc. Alle Rechte vorbehalten.