



Smith+Nephew

VISIONAIRE ♦
Adaptive Guides

Surgical Technique

For Total Knee Systems:

GENESIS ♦ II Total Knee System
LEGION ♦ Total Knee System
JOURNEY ♦ II Total Knee System
ANTHEM ♦ Total Knee Systems

Table of Contents

Introduction	3
Femur Guide	4
Positioning and Exposure	4
Pinning and Drilling	4
Resection	5
Additional ANTHEM [◇] Total Knee System instructions for femur guide	6
Tibia Guide	7
Positioning and Exposure	7
Placement and Pinning	7
Resection for Cruciate-Retaining (CR), Posterior-Stabilized (PS) or Bi-Cruciate Stabilized (BCS) Tibia.....	8
Pin Positioning for Bi-Cruciate Retaining (XR) Tibia.....	9
Appendix A: Standard Distal Femur and Proximal Tibia Resection Guides	11
Appendix B: ANTHEM Total Knee System Standard Distal Femur and Proximal Tibia Resection Guides	12
Appendix C: Adaptive Femur and Tibia Guide Surgeon Design References	13
Catalog Items.....	14

Introduction

Smith+Nephew VISIONAIRE[®] Adaptive Guides are used as a patient-specific surgical instrument to assist in the positioning of total knee replacement components intraoperatively and in guiding the marking of bone before cutting provided that anatomic landmarks necessary for alignment and positioning of the implant are identifiable on patient imaging scans.

The VISIONAIRE Adaptive Guides are indicated for use with the following existing Smith+Nephew knee systems and their cleared indications for use: GENESIS[®] II Total Knee System, LEGION[®] Total Knee System, and JOURNEY[®] II Total Knee Systems. VISIONAIRE Adaptive Guides are available for use with ANTHEM[®] Total Knee System in select markets outside the United States.

The VISIONAIRE sample guides and bone models are intended to inform the surgeon of the expected fit of the VISIONAIRE guide on the patient's anatomy. The bone models represent the intended contact areas of the VISIONAIRE guide and are not intended to inform clinical investigation of the knee anatomy.

Nota Bene

The following technique is for informational and educational purposes only. It is not intended to serve as medical advice. It is the responsibility of treating physicians to determine and utilize the appropriate products and techniques according to their own clinical judgment for each of their patients. For more information on the VISIONAIRE Adaptive Guides, including its indications for use, contraindications, and product safety information, please refer to the product's label and the Instructions for Use.

Note: Please verify the case information and labeling matches the information inscribed on the part.

Note: It is recommended that a sterile back-up instrument set be made available nearby in the event that an Adaptive Guide is intraoperatively determined to be unsuitable for intended use.


Tip: Do not remove osteophytes before evaluating guide fit. If an osteophyte prevents proper fit, remove only the interfering osteophyte before reevaluating fit.

Tip: Avoid over-driving rimmed speed pins. If over-driven into the guide, rimmed speed pins can be extracted by removing the pin capture using a rongeur and backing them out with the driver.


Tip: Paddle extensions are perforated to allow optional removal prior to placement.

For a full list of guide preferences available refer to Appendix A.


Removable Pin Capture for Embedded Pins



Perforated Extension



Removable Tab Extensions



Femur Guide

Positioning and Exposure

1. Flex the knee.
 2. Carefully remove only soft tissue from the anterior femur cortex which may prevent proper guide fit.
-

Pinning and Drilling

3. Place the femur guide on the distal femur by pushing the guide into the trochlear groove and down on the distal condyles.

Note: The proximal contact area of the femur guide should contact the anterior cortex to prevent unintended flexion.

Tip: The alignment rod may be used to verify alignment before pinning the femur guide (see Step 6).




Step 3

4. While the femur guide is firmly held in place, secure the femur guide distally by inserting pins into the distal holes.

Note: The distal pin holes of the femur guide correspond to the spikes of the femur A/P cutting guide associated with the implant. Insert pins deeper than the resection plane.


Tip: It is best for the surgeon to hold the femur guide while an assistant pins.



Step 4

5. Secure the femur guide anteriorly by inserting headless pins into the appropriate anterior parallel pin holes.

Tip: Use of the anterior oblique pin is optional for additional security if desired.



Step 5

Resection

6. Use the external alignment rod to verify proper alignment prior to making the distal resection.


Note: The A/P etch line points to the Rod when assembled correctly and is oriented neutral to the planned alignment of the femur implant. The M/L etch line length is equal to the M/L width of the planned femur implant. The rotation of the M/L etch line is parallel to the preoperatively determined transepicondylar axis.

Note: If interference occurs between soft tissue and the alignment rod, the Quick Connect Handle (71440044) can be mated with the VISIONAIRE[®] alignment connector for additional clearance.

Alignment Rod Assembly



Insert alignment guide



Use light pressure when inserting connector into guide

Step 6


-
7. Remove one distal pin prior to beginning the distal resection. To complete the resection, move the remaining distal pin to the opposite side of the femur guide to preserve three pin fixation.



Step 7

-
8. After completing the resection, remove the femur guide and complete the procedure per the surgical technique recommended for the implant, taking care to ensure resection alignment is acceptable.

Note: The anterior parallel pin holes created through the femur guide correspond to the standard distal cutting block associated with the implant. See Appendix A for the list of available standard guides.




Step 8


Femur Guide (continued)

Additional ANTHEM[◇] Total Knee System Instructions for Femur Guide

9. For ANTHEM Total Knee System, first re-insert non-headed pins into the pre-drilled distal pin holes, then slide the Orthomatch 4-in-1 cut block on to the distal femur.




Step 9a



Step 9b

10. Remove the anterior pins and complete the procedure per the recommended ANTHEM Total Knee System Surgical Technique.



Step 10

See Appendix B for a pin placement guide.

Positioning and Exposure

1. Sublux the tibia.
2. Remove the meniscus, including the posterior medial meniscus.
3. Carefully remove only soft tissue from the anterior tibia cortex which may prevent proper guide fit.

Placement and Pinning

4. Place the tibia guide on the proximal tibia. The primary key contact area for the tibia guide is the anterior medial tibial cortex, the secondary contact area are the medial and lateral plateau.

Note: The key contact areas of the tibia guide should be mated flush to the corresponding anatomy. If not, remove only osteophytes or soft tissue which may be preventing proper fit.

Tip: The alignment rod may be used to verify alignment before pinning the tibia guide (see Step 7).

5. While the tibia guide is firmly held in place, secure the guide proximally by inserting pins into the proximal holes.

Note: With the exception of the TCF Guide, the proximal pin holes correspond to the pin holes of the tibial baseplate trial associated with the implant. Insert pins deeper than the resection plane.

Tip: It is best for the surgeon to hold the tibia guide while an assistant pins.

6. Secure the tibia guide anteriorly by inserting headless pins into the appropriate anterior parallel pin holes.

Tip: Use of one or more oblique anterior pin holes with short headed pins is optional for added security if desired.

TCF Approach




DCF Approach



Step 4



Step 5



Step 6

Tibia Guide (continued)

7. Use the external alignment rod to verify proper alignment prior to making the proximal resection.

Note: The anterior line corresponds to the sagittal mid-plane of the tibial implant and the intended M/L position of the alignment rod.

Tip: The medial side of the guide includes a symbol which indicates that the Alignment Rod will align either Parallel \parallel to the tibia mechanical axis or Perpendicular \perp to the cutting slot.



Step 7

Resection for Cruciate-Retaining (CR), Posterior-Stabilized (PS) or Bi-Cruciate-Stabilized (BCS) Tibia

8. Remove one distal pin prior to beginning the distal resection. To complete the resection, move the remaining distal pin to the opposite side of the femur guide to preserve three pin fixation.



Step 8

9. After completing the resection, remove the tibia guide and complete the procedure per the surgical technique recommended for the implant, taking care to ensure resection alignment is acceptable.

Note: For ANTHEM[®] CR or PS, the proximal pin holes align with the anterior pin hole location on the Orthomatch Stemless Tibia Trials.

Tip: The anterior parallel pin holes created through the tibia guide correspond to the standard proximal tibia cutting guide specified for the case. See Appendix A for the list of available standard guides.



Step 9


Pin Positioning for Bi-Cruciate Retaining (XR) Tibia

10. Mark the tibia with the M/L and A/P position of the tibial eminence to be preserved using the 'goal post' feature of the guide. Mark the internal/external rotation of the tibia eminence to be preserved using the straight inside edge of the proximal medial paddle. Pin through the XR anterior parallel pin holes labelled 'XR'.

Note: Verify that the ACL is viable before proceeding. If not, continue with the JOURNEY[®] II Cruciate Retaining Knee System or JOURNEY II Bi-Cruciate Stabilized Knee System technique using the non-XR anterior parallel pin holes.



Step 10a



Step 10b

11. If satisfied with the position and rotation of the tibia eminence to be preserved relative to the ACL footprint, predrill through the pin hole lateral to the cut slot using a 1/8" (3.2mm) drill or 110mm Headless Speed Pin.

Note: The ACL may be obscured in the MRI, therefore the VISIONAIRE[®] Engineer sets the rotation and placement as a result of best fit for the asymmetric base plate.



Step 11a

Potential damage: Verify guide placement does not interfere with ACL.



Step 11b

Incorrect Position and Rotation – Do NOT Drill


Tibia Guide (continued)

- 12.** Remove the guide and mate the 3° Datum Block to the anterior parallel pins. Center the XR Orientation Stylus on the 'goal post' mark and orientate the arms parallel with the medial line. Complete the procedure using the JOURNEY[®] II XR[®] Bi-Cruciate Retaining Knee System surgical technique taking care to ensure resection alignment is acceptable.


Note: The Adaptive Guide XR pin guide is designed to match the 3° XR Datum Block. If additional slope is needed, the 5° XR Datum Block should be used as a recut guide.

Tip: If predrilled, the pin hole in-line with the cut slot can be used to position and rotate the Orientation Stylus.

Tip: Later in the XR Technique, use the 'goal post' marking to position the XR Anterior Eminence Chisel.



Step 12a



Step 12b

Appendix A: Standard Distal Femur and Proximal Tibia Resection Guides

VISIONAIRE[◇] Adaptive Guides provide anterior parallel pin holes which correspond to Smith+Nephew conventional distal femur and proximal tibia resection guides. On the anterior face of each Adaptive Guide, the conventional resection guide intended for use is labeled. Below are pictured the standard resection guides which correspond to each text label.



Corresponding Standard Femur Guides



Corresponding Standard Tibia Guides

SYNC





7401-6204

SYNC



L 7401-6216
R 7401-6217

SYNC SYM



7401-6227

GII





7144-1147

GII MIS



L 7144-1136
R 7144-1137

GII TAA




L 7144-0209
R 7144-0205

JNY




7401-2233

JNY MIS




L 7401-8411
R 7401-7411

JII MIS



L 7401-4473
R 7401-4474


JII XR[◇]



L 7401-4463
R 7401-4464


Appendix B: ANTHEM^o Standard Distal Femur and Proximal Tibia Resection Guides

ANTHEM Femoral



7144-2109

ANTHEM Tibial



7144-2109



Use Posterior Referencing Holes

Appendix C: Adaptive Femur and Tibia Guide Surgeon Design Preferences

Femur	Options	Tibia	
	 <p data-bbox="427 617 477 642">DCF</p>	<p data-bbox="599 331 878 600">To accommodate techniques with tighter joint spaces, a lower profile tibia guide with angled proximal pins is available. These pin holes do not set rotation.</p>  <p data-bbox="1036 617 1086 642">DCF</p>	 <p data-bbox="1360 617 1411 642">TCF</p>
 <p data-bbox="180 884 224 909">Tab</p>	 <p data-bbox="410 884 485 909">No Tab</p> <p data-bbox="599 674 873 806">Guides can be designed with extended contact tabs or without.</p>	 <p data-bbox="1040 884 1084 909">Tab</p>	 <p data-bbox="1344 884 1419 909">No Tab</p>
	 <p data-bbox="427 1163 477 1188">MIS</p>	<p data-bbox="599 947 850 1146">Along with the standard MIS cutting slot, the tibia guide is offered with an elongated TAA cutting slot.</p>  <p data-bbox="1040 1163 1084 1188">MIS</p>	 <p data-bbox="1360 1163 1411 1188">TAA</p>
	 <p data-bbox="375 1451 532 1476">Perpendicular</p>	<p data-bbox="599 1220 873 1419">The guide can be designed to position the drop rod parallel to the mechanical axis or perpendicular to the cutting slot.</p>  <p data-bbox="906 1356 997 1430">Compare to VISIONAIRE legacy alignment</p> <p data-bbox="987 1451 1144 1476">Perpendicular</p>	 <p data-bbox="1446 1356 1537 1430">Compare to conventional instrument alignment</p> <p data-bbox="1360 1451 1451 1476">Parallel</p>
	 <p data-bbox="367 1703 537 1728">JOURNEY II CR</p>	<p data-bbox="599 1514 873 1608">The Adaptive Guide for JOURNEY[◇] II supports BCS, CR and XR[◇].</p>  <p data-bbox="1040 1703 1084 1728">XR</p>	 <p data-bbox="1365 1703 1409 1728">CR</p>
 <p data-bbox="154 1976 245 2001">Rimmed</p>	 <p data-bbox="380 1976 526 2001">Non Rimmed</p> <p data-bbox="599 1766 857 1860">Rimmed pin holes are offered for use with headed speed pins.</p>	 <p data-bbox="1013 1976 1104 2001">Rimmed</p>	 <p data-bbox="1312 1976 1458 2001">Non Rimmed</p>

Catalog Items

VISIONAIRE[◊] Adaptive Guide Alignment Verification Tools

71440302 VISIONAIRE Lightweight Alignment Rod

71440303 VISIONAIRE Alignment Connector

VISIONAIRE Adaptive Guide Case Specific Samples

V0100103 VISIONAIRE Adaptive Guide Sample Bone And Cutting Guide

VISIONAIRE Adaptive Guide Sterile Part Numbers

V0100104 VISIONAIRE Adaptive Guide Distal Femur GENESIS[◊] II

V0100105 VISIONAIRE Adaptive Guide Proximal Tibia GENESIS II

V0100106 VISIONAIRE Adaptive Guide Kit GENESIS II

V0100107 VISIONAIRE Adaptive Guide Distal Femur LEGION[◊] Primary

V0100108 VISIONAIREz Adaptive Guide Proximal Tibia LEGION Primary

V0100109 VISIONAIRE Adaptive Guide Kit LEGION Primary

V0100110 VISIONAIRE Adaptive Guide Distal Femur JOURNEY[◊] II

V0100111 VISIONAIRE Adaptive Guide Proximal Tibial JOURNEY II

V0100112 VISIONAIRE Adaptive Guide Kit JOURNEY II

V0100125 VISIONAIRE Adaptive Guide Distal Femur ANTHEM[◊]

V0100126 VISIONAIRE Adaptive Guide Proximal Tibia ANTHEM

V0100127 VISIONAIRE Adaptive Guide Kit ANTHEM

VISIONAIRE[◊] Adaptive Guide Non-Sterile Part Numbers

V0200104	Nonsterile VISIONAIRE Adaptive Guide Distal Femur GENESIS [◊] II
V0200105	Nonsterile VISIONAIRE Adaptive Guide Proximal Tibia GENESIS II
V0200106	Nonsterile VISIONAIRE Adaptive Guide Kit GENESIS II
V0200107	Nonsterile VISIONAIRE Adaptive Guide Distal Femur LEGION [◊] Primary
V0200108	Nonsterile VISIONAIRE Adaptive Guide Proximal Tibia LEGION Primary
V0200109	Nonsterile VISIONAIRE Adaptive Guide Kit LEGION Primary
V0200110	Nonsterile VISIONAIRE Adaptive Guide Distal Femur JOURNEY [◊] II
V0200111	Nonsterile VISIONAIRE Adaptive Guide Proximal Tibial JOURNEY II
V0200112	Nonsterile VISIONAIRE Adaptive Guide Kit JOURNEY II
V0200125	Nonsterile VISIONAIRE Adaptive Guide Distal Femur ANTHEM [◊]
V0200126	Nonsterile VISIONAIRE Adaptive Guide Proximal Tibia ANTHEM
V0200127	Nonsterile VISIONAIRE Adaptive Guide Kit ANTHEM

SPEED PIN Adaptors

74013489	SPEED PIN Quick Connect
74016461	SPEED PIN Quick Connect

Headed SPEED PIN

74013471	45mm
74013472	65mm
74016465	MIS 45mm
74016466	MIS 65mm

Non-Headed SPEED PIN

74013480	65mm
----------	------

Drill Bit

114969	3.2mm X 76mm
--------	--------------

Trocar Pin

71210002	GEN TROCAR Pin 1/8 X 3
74016461	3.2mm X 75mm

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area.

Smith & Nephew, Inc.
1450 Brooks Road
Memphis, Tennessee 38116
USA

www.smith-nephew.com

Telephone: 1-901-396-2121
Information: 1-800-821-5700
Orders/Inquiries: 1-800-238-7538

◇ Trademark of Smith+Nephew
All Trademarks acknowledged
©2022 Smith+Nephew
All rights reserved.
07684 V5 71282151 REVD 06/22