

## Tibiototalcalcaneal arthrodesis with the ILIZAROV method



# Tibiotalocalcaneal arthrodesis with the ILIZAROV™ method

## Surgical Technique

### Table of contents

<b>Introduction</b> .....	2
Ring sizing .....	3
Pre-building the frame .....	3
Safe zones in the tibia and foot .....	4
Wire insertion technique .....	5
Half pin insertion technique .....	6
<b>Surgical technique</b> .....	8
Operative dressings .....	14
Postoperative pin care .....	14
<b>Catalog information</b> .....	16

### **Nota Bene**

The technique description herein is made available to the healthcare professional to illustrate the author's suggested treatment for the uncomplicated procedure. In the final analysis, the preferred treatment is that which addresses the needs of the specific patient. For more information on the products in this surgical technique, including indications for use, contraindications, effects, precautions and warnings, please consult the products' Instructions for Use (IFU).

## Introduction

In patients requiring ankle arthrodesis, the ILIZAROV™ method allows a percutaneous approach and early weight-bearing which has been shown to lead to high fusion rates and improved patient outcomes when compared to traditional internal fixation treatment methods.<sup>7-9</sup>



## Ring sizing

Place a full ring over the affected limb. Make sure there is a two finger clearance between the inner diameter of the ring and the soft tissues of the leg and the foot. See **Figure 1**.




Figure 1

## Pre-building the frame

A pre-assembled ILIZAROV™ foot frame may be selected (see **Figure 2**) or a frame may be built by the physician. When building the frame, for the tibia ring fixation block choose two full pre-sized circular rings, one pre-sized foot ring and one half ring, all of the same diameter. Using 10mm nuts, place four 150mm threaded rods equidistant from each other between the two full rings, allowing for maximum separation between the rings. Place two 120mm threaded rods between the foot ring and the tibia fixation block in the posterior aspect of the frame. Leave approximately 15mm to 20mm of thread length remaining proximally to the distal tibia ring to allow for compression or distraction at the ankle joint and for application versatility. Using two male hinges, two bolts and four 10mm nuts, close the open end of the foot ring with the chosen half ring.




Figure 2

## Safe zones in the tibia and foot

Care must be taken to avoid neurovascular structures and intraarticular penetration. **Figure 3** and **Figure 4** show suggested safe pin and wire pathways that allow for minimal transfixion of muscle compartments and neurovascular structures of the lower leg, ankle and foot.




Figure 3




Figure 4

## Wire insertion technique

When inserting thin wires it is important to follow the guidelines below.

### Step 1

Manually insert the wire to the bone to minimize soft tissue trauma.

### Step 2

Only pass wires through muscle compartments in their elongated state to avoid excessive tension. See **Figures 5.1** and **5.2**.




Figure 5.1




Figure 5.2

### Step 3

Consider holding the wire with a wet sponge to dissipate the heat generated from drilling. See **Figure 6**. After drilling the wire through the bone use a mallet to pass wires through the remaining soft tissues to prevent twisting.




Figure 6

#### Step 4

Pause frequently when drilling to allow heat to dissipate.

#### Step 5

Fixate the wire to the frame in its natural resting place to prevent pin site irritation.

### Half pin insertion technique

When inserting half pins it is important to follow the guidelines below.

#### Step 1

Attach a rancho cube to the ring and place a drill sleeve with a trocar through the furthest hole from the ring. See **Figure 7**.

#### Step 2

Make a skin incision and push the trocar and sleeve down to bone. See **Figure 7**.




Figure 7



**Step 3**

Remove the trocar and place the drill through the sleeve. Drill through the bone taking frequent pauses to allow for heat dissipation. See **Figure 8**.




Figure 8

**Step 4**

Select half pin length and insert through the drill sleeve until bicortical purchase is obtained. See **Figure 9**.




Figure 9

**Step 5**

Remove the drill sleeve. Insert the proper diameter centering sleeve and lock with a short bolt or set screw. See **Figure 10**.




Figure 10

# Surgical technique

## Step 1

The patient is placed supine on a radiolucent table with the affected side hip bolstered up. The patella should remain pointing to the ceiling without any additional support. See **Figure 11**.




Figure 11

## Step 2

A direct lateral approach is made centered on the fibula. See **Figure 12**.




Figure 12

## Step 3

Cut the fibula above the level of the syndesmosis and rotate the fibula posteriorly by releasing the soft tissue attachments anteriorly. This approach allows direct access to the lateral ankle joint and the subtalar joint. The distal fibula can be discarded if the bone is infected or atrophic. See **Figure 13**.




Figure 13

## Step 4

Prepare the subtalar and ankle joint surfaces for arthrodesis by removing residual cartilage and preparing subchondral bone. Any multiplanar alignment correction should be done at this time. See **Figure 13**.

### Step 5

Provisionally secure the ankle and the hindfoot in neutral position with a Steinman pin.

### Step 6

Close incision sites. See **Figure 14**.




Figure 14

### Step 7

Place a bolster under the knee. Then place the pre-assembled frame over the foot and tibia. Place a stack of towels between the posterior aspect of the frame and the posterior soft tissues of the tibia and foot until the threaded rods of the frame are parallel to the posterior cortex of the tibia. Adjust the proximal/distal positioning of the frame until the foot ring is at the midpoint of the calcaneus while the ankle is in neutral position. See **Figure 15**.



Figure 15

### Step 8

Follow the wire and pin safe zone guidelines and insertion techniques section on pages 5 through 7 for wire and half pin insertion locations. Insert a bicortical reference wire from medial to lateral on the inferior side of the proximal tibial ring that is perpendicular to the axis of the shaft of the tibia. Attach the wire directly to the ring with wire fixation bolts and tension to 130 ft lbs. Verify that the axis of the frame is in line with the axis of the tibia in both the AP and lateral radiograph views. Insert a bicortical wire from medial to lateral on the inferior side of the distal tibia ring. Attach the wire directly to the ring with wire fixation bolts and tension to 130 ft lbs. If a wire is not directly abutting the ring, do not pull it to the ring. Washers, hinges or supports can be used in conjunction with the wire fixation bolts to build the frame to the wire. See **Figure 16**.



Figure 16

### Step 9

Place two wires on the superior aspects of each tibial ring from anterolateral to posteromedial. Attach the wire directly to the ring with wire fixation bolts and tension to 130 ft lbs. See **Figure 17**.



Figure 17

### Step 10

Place a rancho cube on the anteromedial superior side of the distal and proximal tibia ring. Using the appropriate diameter drill bit for the chosen diameter half pin, drill a bicortical hole for half pin insertion and place the half pin through the centering sleeve until one or two threads have cleared the far cortex. Verify bicortical thread purchase under fluoroscopic radiograph. Lock the half pin into the rancho cube with a short bolt or a set screw. See **Figure 18**.



Figure 18

### Step 11

Attention is directed to the application of the foot ring. If necessary, Achilles tendon or gastrocnemius lengthening may be performed at this time. Ensure the hindfoot and ankle are in neutral position to the tibia and insert an olive wire in the calcaneus from anteromedial to posterolateral on the inferior side of the foot ring. Attach and tension the wire to 90-110 ft lbs. Insert an olive wire in the calcaneus from anterolateral to posteromedial on the superior side of the foot ring. Attach and tension the wire to 90-110 ft lbs. See **Figure 19**.




Figure 19

### Step 12

Insert an olive wire in the first and second metatarsals from medial to lateral and a second olive wire in the fifth, fourth and third metatarsals from lateral to medial. Attach them to the ring and tension to 90 ft lbs. See **Figure 20**.




Figure 20

### Step 13

Attach two 120mm threaded rods to empty holes in the distal portion of the foot ring. Connect the proximal portion of these to the distal tibia ring with the appropriate length connection plate to stabilize the foot ring to the tibial fixation block. See **Figure 21**.



Figure 21

### Optional step

If compression of the subtalar joint is desired, wires may be inserted into the talus. A bayonet wire is inserted from medial to lateral through the talar neck and/or a bayonet wire is inserted from anterolateral to posteromedial through the talar body. One or each of these is arched slightly toward the foot plate and attached to the foot ring. Tension is applied to one or both of the talar wires achieving compression of the subtalar joint as the wire(s) translate inferiorly. See **Figure 22**.



Figure 22

### Step 14

Loosen the 10mm nuts on the inferior side of the distal tibia ring approximately 1cm below the ring. Remove the provisional Steinman pin from the hindfoot and manually compress and hold the foot ring toward the distal tibia ring. Tighten the 10mm nuts on the superior side of the distal tibia ring to hold the compression across the arthrodesis site. Tighten the 10mm nuts on the inferior side to the distal tibia ring. See **Figure 23**.



Figure 23



**Step 15**

Cut half pins with a rod cutter and place protective pin caps on the ends. Check and retighten all nuts and bolts on the frame to ensure stability. See **Figures 24** and **25**.



Figure 24



Figure 25

## Operative dressings

Apply sterile dressings over surgical incisions. Apply ACTICOAT<sup>®</sup> EXFIX or ILIZAROV<sup>™</sup> sponges with disc clips around the pin sites and compressive dressings completely around the leg within the frame. Wrap the exterior of the frame with a long 6-inch Ace<sup>®</sup> bandage. Leave in place for five to seven days before beginning pin site care. See **Figure 26**.




Figure 26

## Postoperative pin care

On postoperative day four initiate routine daily pin site care. Using sterile cotton swabs, clean each pin and wire site with a 1:1 ratio mixture of saline and hydrogen peroxide. After the surgical wounds have healed and the sutures have been removed (two-four weeks) the patient can shower daily with antibacterial soap and should continue routine daily pin site care. See **Figure 27**.



Figure 27

**Note** Do not cross contaminate pin sites with the cotton swabs.



# Catalog information



## ILIZAROV™ Foot and Ankle Tray

Set No. 7107-0605

## Instruments

Cat. No.	Description	Tray Qty	Cat. No.	Description	Tray Qty
102905	Combination Wrench, 10mm	2	112716	Hex Fix Pin Driver/Extractor	1
102911	Fixation Bolt Wrench	1	7107-0344	Wire Cutter	1
103101	Dynamometric Wire Tensioner	2	102910	10mm Slotted Hexdriver	1
L202002	Combination Wrench, 10mm x 13mm	2			

## Hardware

Cat. No.	Description	Tray Qty	Cat. No.	Description	Tray Qty
101800	Short Connector Plate, 3 Hole, 45mm	4	102302	Threaded Rod, 100mm	4
102550	Universal Joint	4	102303	Threaded Rod, 120mm	4
71011000	Adjustable Pin Clamp/Locking Hinge	2	102304	Threaded Rod, 150mm	4
103300	Nut, 10mm	60	102305	Threaded Rod, 200mm	4
103301	Nut, 10mm, Nylon Insert	10	201700	Conical Washer Couple	16
103302	DC Counter, 4 Point	4	100600	Wire Fixation Bolt, Cannulated	12
101400	Male Support, 2 Hole	4	100700	Wire Fixation Bolt, Slotted	12
101401	Male Support, 3 Hole	2	103405	5mm Centering Sleeve	6
101402	Male Support, 4 Hole	2	7192-9154	Russian Nuts	6
102301	Threaded Rod, 80mm	4			

## Hardware – Mayo Tray

<b>Cat. No.</b>	<b>Description</b>	<b>Tray Qty</b>	<b>Cat. No.</b>	<b>Description</b>	<b>Tray Qty</b>
71070394	Ilizarov Foot and Ankle Mayo Stand Caddy	1	103200	Bolt, 10mm	20
103452	Rancho Cube, 2 Hole	4	103201	Bolt, 16mm	20
103453	Rancho Cube, 3 Hole	4	103406	6mm Centering Sleeve	6
103454	Rancho Cube, 4 Hole	2	100600	Wire Fixation Bolt, Cannulated	12
101600	Male Hinge	6	100700	Wire Fixation Bolt, Slotted	12
101602	Low Profile Male Hinge	6	102706	Washer, 2.0mm, Fixation Bolt	20
103300	Nut, 10mm	40	102707	Washer, 4.0mm, Fixation Bolt	20

## Wires

<b>Cat. No.</b>	<b>Description</b>	<b>Tray Qty</b>
7107-0390	Wire Caddy	1
102107	Wire with Stopper, 400mm	10
102102	Wire Bay Point Cortical, 370mm	10

# References

- 1 Pounder NM, Harrison AJ. Low intensity pulsed ultrasound for fracture healing: a review of the clinical evidence and the associated biological mechanism of action. *Ultrasonics*. 2008;48(4):330-338
- 2 Sant'Anna EF, Leven RM, Viridi AS, Sumner DR. Effect of low intensity pulsed ultrasound and BMP-2 on rat bone marrow stromal cell gene expression. *J Orthop Res*. 2005;23(3):646-652.
- 3 Naruse K, Miyauchi A, Itoman M, Mikuni-Takagaki Y. Distinct anabolic response of osteoblast to low-intensity pulsed ultrasound. *J Bone Miner Res*. 2003;18(2):360-369.
- 4 Chen YJ, Wang CJ, Yang KD, et al. Pertussis toxin-sensitive Galphai protein and ERK-dependent pathways mediate ultrasound promotion of osteogenic transcription in human osteoblasts. *FEBS Lett*. 2003;554(1-2):154-158.
- 5 Wang FS, Kuo YR, Wang CJ, et al. Nitric oxide mediates ultrasound-induced hypoxia-inducible factor-1 alpha activation and vascular endothelial growth factor-A expression in human osteoblasts. *Bone*. 2004;35:114-123.
- 6 Zhou S, Schmelz A, Seufferlein T, Li Y, Zhao J, Bachem MG. Molecular mechanisms of low intensity pulsed ultrasound on human skin fibroblasts. *J Biol Chem*. 2004;279(52):54463-54469
- 7 Morasiewicz P, Dejnek M, Urbanski W, Dragan S, Kulej M, Dragan S. Radiological evaluation of ankle arthrodesis with Ilizarov fixation compared to internal fixation. *Injury*. 2017;48(7):1678-1683.
- 8 Ankle Arthrodesis using Ilizarov Ring Fixator: A Primary or Salvage Procedure? An Analysis of Twenty Cases Hasan O, FCPS, Fahad S, MBBS, Sattar S, MBBS, Umer M, FCPS, Rashid H, FCPS Malaysian Orthopaedic Journal 2018 Vol 12 No 3
- 9 Complex Ankle Arthrodesis Using the Ilizarov Method Yields High Rate of Fusion Austin T. Fragomen MD, Eugene Borst BA, Lindsay Schachter BS, Stephen Lyman PhD, S. Robert Rozbruch MD *Clin Orthop Relat Res* (2012) 470:2864–2873

**Orthopaedics**

Smith & Nephew, Inc.  
1450 Brooks Road  
Memphis, TN 38116  
USA

[www.smith-nephew.com](http://www.smith-nephew.com)

Telephone: 1-901-396-2121  
Information: 1-800-821-5700  
Orders/Inquiries: 1-800-238-7538