MICRORAPTOR[®] REGENESORB Suture Anchor

10601382-D

Smith & Nephew Orthopaedics GmbH
Alemannenstraße 14
78532 Tuttlingen
Germany
ec.rep@smith-nephew.com

Endoscopy
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810
USA

04/2019 10601382 Rev. D

MICRORAPTOR REGENESORB Suture Anchor

MICRORAPTOR REGENESORB-
Fadenanker

Anclaje para sutura

MICRORAPTOR REGENESORB

Ancrage de suture

MICRORAPTOR REGENESORB

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MICRORAPTOR REGENESORB

MICRORAPTOR REGENESORB

suturankare

MICRORAPTOR REGENESORB

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Âncora de sutura

MICRORAPTOR REGENESORB

MICRORAPTOR REGENESORB

suturanker

MICRORAPTOR REGENESORB

suturanker

MICRORAPTOR REGENESORB

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MICRORAPTOR REGENESORB

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Instructions for Use

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Mode d'emploi

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Bruksanvisning

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Brugsanvisning

Bruksanvisning

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MICRORAPTOR[®] REGENESORB Suture Anchor

CE mark on the package insert (IFU) is not valid unless there is a CE mark on the product (description) label.

Device Description

Smith & Nephew MICRORAPTOR REGENESORB Suture Anchor is a fixation device intended to provide secure attachment of soft tissue to bone. It consists of a suture anchor with attached non-absorbable suture preassembled to a flexible insertion device. The anchor is preloaded with one #1 suture. This device is provided sterile, for single use only.

Material

The suture anchor is made from PLGA (poly(lacticoglycolic acid)), β -TCP (beta tricalcium phosphate), and calcium sulfate.

Contents

MICRORAPTOR REGENESORB Suture Anchor contains:

- 1 ea. Absorbable biocomposite suture anchor
- 1ea. Flexible insertion device – flexible stainless steel shaft with a polycarbonate and ABS handle
- 1 ea. Sutures – Braided, uncoated, UHMW polyethylene, non-absorbable

Refer to individual suture anchor product labels for suture size, type, and quantities.

Intended Use

The Smith & Nephew MICRORAPTOR REGENESORB Suture Anchor is intended for the reattachment of soft tissue to bone.

Indications for Use

The Smith & Nephew MICRORAPTOR REGENESORB Suture Anchor is intended for the reattachment of soft tissue to bone for the following indications:

Hip

Acetabular labrum repair/reconstruction

Shoulder

- Capsular stabilization
 - Bankart repair
 - Anterior shoulder instability
 - SLAP lesion repairs
 - Capsular shift or capsulolabral reconstructions
- Biceps tenodesis

Contraindications

- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Pathological conditions of bone, such as cystic changes or severe osteopenia, which would compromise secure anchor fixation.
- Pathological conditions in the soft tissues to be attached that would impair secure fixation by suture.
- Comminuted bone surface, which would compromise secure anchor fixation.
- Physical conditions which would eliminate, or tend to eliminate, adequate anchor support or retard healing.
- The presence of infection.
- Conditions which would limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- Contraindications may be relative or absolute and must be carefully weighed against the patient’s entire evaluation.

Warnings

Do not use if package is damaged. Do not use if the product sterilization barrier or its packaging is compromised.

- Contents are sterile unless package is opened or damaged. **DO NOT RESTERILIZE.** For single use only. Discard any open, unused product. Do not use after the expiration date.
- It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.
- Read these instructions completely prior to use.
- Product must be stored in the original sealed pouch.
- Incomplete anchor insertion may result in the anchor not functioning as intended.
- Breakage of the inserter can occur if the insertion site is not prepared with recommended Smith & Nephew, drill guides & drills prior to insertion.
- Only use the recommended drills, drill guides, and obturators intended for use with the MICRORAPTOR® REGENESORB Suture Anchor. Use of other instruments may injure the patient, damage the instruments, or compromise fixation.
- Maintaining drill guide alignment throughout drilling is required to ensure drill integrity.
- Do not attempt to implant this device within cartilage epiphyseal growth plates or non-osseous tissue.
- Do not resterilize or reuse anchors, sutures and insertion devices packaged with the MICRORAPTOR REGENESORB Suture Anchor.

Precautions

R_{ONLY} U.S. Federal law restricts this device to sale by or on the order of a physician.

- Hazards associated with reuse of this device include, but are not limited to, patient infection and/or device malfunction.
- Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.
- Postoperative care is important. Instruct the patient on the limitations of the implant and caution them regarding weight bearing and body stresses on the appliance prior to secure bone healing.
- Do not use sharp instruments to manage or control the suture.
- As in all suture anchor or suturing techniques, the fixation given should be considered as only temporary, until biological attachment of tissue to bone is completed, and may not withstand weight bearing or other unsupported stresses. The suture anchor and suture are not intended to provide indefinite biomechanical integrity.
- Implantation of the suture anchor requires preparation of the insertion site. Predrilling with the appropriate Smith & Nephew MICRORAPTOR drill is the preferred method of site preparation.
- Ensure the anchor placement is aligned with the drilled hole. Proper alignment is essential for successful anchor implantation.
- Use of excessive force during insertion can cause failure of the suture anchor or insertion device.
- Bone quality must be adequate to allow proper placement of suture anchor.
- Do not alter implant or instrumentation or performance may be compromised.
- Once seated, do not rotate the suture anchor device in the bone as this may cause device failure.
- Postoperative range of motion is to be determined by the physician.
- After use, this device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.
- Do not remove the drill guide from the anchor insertion site
- For dense bone, it may be required to cycle the drill in and out of the implantation site to ensure complete removal of bone.

Adverse Reactions

- Breakage of the suture can occur
- Loss of fixation or pullout of suture anchors can occur
- Mild inflammatory reaction
- Foreign body reaction
- Infection, both deep and superficial
- Allergic reaction

Instructions for Use

 **WARNING:** Use of the appropriate Smith & Nephew instrumentation is necessary to prepare the insertion site and maintain axial alignment between insertion site and the MICRORAPTOR® REGENESORB Suture Anchor.

When using a MICRORAPTOR REGENESORB Suture Anchor, use a Smith & Nephew drill guide(s), obturator(s), and drill(s) specific to the MICRORAPTOR REGENESORB suture anchor to prepare the insertion site (each sold separately).

1. Place the distal tip of the drill guide onto the bone at the desired implantation site. Drill guide placement may be aided using an optional obturator.
2. Hold the drill guide firmly in place, use the following recommended drills to prepare the insertion site:

| Shoulder General Use | Shoulder Hard Bone |
|-------------------------------|-------------------------------|
| REF 72204988, 1.6 mm Drill | REF 72205267, 1.8 mm Drill |

| Hip General Use | Hip Soft Bone |
|-------------------------------|-------------------------------|
| REF 72205267, 1.8 mm Drill | REF 72204988, 1.6 mm Drill |

- a. Verify that the drill is rotating forward (clockwise) while drilling. The depth stop on the drill bottoms out on the handle of the drill guide when proper hole depth is reached.
- b. Please note 1.8 mm Hard bone drill will have a laser marked band on shaft near hub.

CAUTION: For dense bone, it may be required to cycle the drill in and out of the implantation site to ensure complete removal of bone.

3. While holding the drill guide steady and maintaining alignment, remove the drill from the insertion site.
CAUTION: Do not remove the drill guide from the anchor insertion site.
4. Holding the drill guide firmly in place, insert the preloaded suture anchor into the proximal end of the drill guide. Advance the insertion device down the drill guide until resistance is met and anchor is visible in the windows of the drill guide.
5. While maintaining distal pressure on the insertion device, push and/or tap on the proximal end of the insertion device to insert the anchor.

Note: It is recommended to insert with several light taps. Proper insertion depth is achieved when the insertion device handle makes contact with the proximal end of the drill guide handle. This positions the suture anchor at the appropriate depth below the surface of the bone.

6. Disengage the suture anchor from the insertion device:
 - a. Hold the handle securely in the palm of the hand and pull back on the finger grips using two fingers to retract the suture tensioning mechanism.
 - b. While holding the suture tensioning mechanism in the back position with the finger grips, slowly remove the insertion device from the insertion site until the suture tips are visible. The suture releases and feeds through the insertion device as it is removed.

Alternatively:

- a. Pull back on the finger grips and unwind the suture with the free hand, then release the finger grips.
- b. Slowly remove the insertion device.
- c. As the insertion device is removed, grasp the suture below the insertion device using the free hand to verify suture release from the insertion device.

CAUTION: Do not use sharp instruments to manage or control the suture.

7. Discard the insertion device.

CAUTION: *After use, this device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.*

8. Pass the repair sutures through or around the soft tissue to be refixated using the desired suture passing techniques and devices.
9. Secure the repair of the tissue with standard knot tying techniques. When complete, cut and remove any excess suture.
10. Place additional suture anchors as required to complete the procedure.

Warranty

For single use only. This product is warranted to be free from defects in material and workmanship. Do not reuse.

For Further Information

If further information on this product is needed, contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized representative.

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©2018, 2019 Smith & Nephew, Inc. Med ensamrätt.
°Varumärke som tillhör Smith & Nephew, och har registrerats vid U.S. Patent and Trademark Office (USA:s Patent- och registreringsverk).
°Varemærke, der tilhører Smith & Nephew og er registreret i det amerikanske patent- og varemærkedirektorat (U.S. Patent and Trademark Office).