
# HEALICOIL<sup>◊</sup> Knotless PK Suture Anchor

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Brugsanvisning

Bruksanvisning

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Reference product (description) label for MR compatibility.

## Device Description

The Smith & Nephew HEALICOIL Knotless PK Suture Anchor is a fixation device intended to provide secure reattachment of soft tissue to bone.

The device consists of a knotless suture anchor assembled on to an insertion device.

The device is to be used with Smith & Nephew ULTRABRAID, ULTRATAPE, and MINITAPE Sutures.

This device is provided sterile, for single use only.

## Material

	HEALICOIL Knotless PK Self-Tapping	HEALICOIL Knotless PK
Distal Implant	Titanium	PEEK
Internal Plug	PEEK	
Proximal implant	PEEK	

## Contents

The Smith & Nephew HEALICOIL Knotless PK Suture Anchor contains:

1 ea. Non-absorbable PEEK suture anchor with non-absorbable PEEK distal implant and non-absorbable PEEK internal plug

or

1 ea. Non-absorbable PEEK suture anchor with non-absorbable titanium distal implant and non-absorbable PEEK internal plug

with

1 ea. Insertion device – stainless steel shaft with ABS/polycarbonate handle

1 ea. Suture threader – nitinol wire crimped into a stainless steel tube with ABS handle.

## Intended Use

The Smith & Nephew HEALICOIL Knotless PK Suture Anchor is intended for the reattachment of soft tissue to bone.

## Indication for Use

Smith & Nephew HEALICOIL Knotless PK Suture Anchor is intended for use only for the reattachment of soft tissue to bone for the following indications:

### Shoulder

- Biceps tenodesis
- Rotator cuff tear repair

## Contraindications

- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Pathological conditions of bone, such as cystic changes or severe osteopenia, which would compromise secure anchor fixation.
- Pathological conditions in the soft tissues to be attached that would impair secure fixation by suture.
- Comminuted bone surface, which would compromise secure anchor fixation.
- Physical conditions which would eliminate, or tend to eliminate, adequate anchor support or retard healing.
- The presence of infection.
- Conditions which would limit the patient's ability or willingness to restrict activities or follow direction during the healing period.
- Contraindications may be relative or absolute and must be carefully weighed against the patient's entire evaluation.

## Warnings

-  Do not use if package is damaged. Do not use if the product sterilization barrier or its packaging is compromised.
-  Contents are sterile unless package is opened or damaged. **DO NOT RESTERILIZE.** For single use only. Discard any open, unused product. Do not use after the expiration date.
- It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.
- Read these instructions completely prior to use.
- Product must be stored in the original sealed pouch.
- For HEALICOIL Knotless PK Self-Tapping devices; titanium alloys contain elements that may stimulate allergic hypersensitive responses by the immune system. These elements are titanium, aluminum, and vanadium (Ti, Al, and V). When sensitivity is anticipated, appropriate preoperative testing should be conducted.
- Incomplete anchor insertion may result in the anchor not functioning as intended.
- Breakage of the suture anchor can occur if insertion sites are not prepared with the appropriate instrumentation prior to implantation.
- Only use the recommended hole preparation and drills for use with the HEALICOIL Knotless PK Self-Tapping Suture Anchor. Use of other instruments may injure the patient, damage the instruments, or compromise fixation.
- Maintaining inserter alignment throughout insertion is required to ensure implant integrity.
- Do not attempt to implant this device within cartilage epiphyseal growth plates or non-osseous tissue.
- Incomplete anchor insertion may result in poor anchor performance.

## Precautions

 U.S. Federal law restricts this device to sale by or on the order of a physician.

- Hazards associated with reuse of this device include, but are not limited to, patient infection and/or device malfunction.
- Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.
- Postoperative care is important. A patient should be instructed on the limitations of the implant and should be cautioned regarding weight bearing and body stresses on the appliance prior to secure bone healing. Postoperative range of motion is to be determined by the physician.
- Do not use sharp instruments to manage or control the suture.
- As in all suture anchor or suturing techniques, until biological attachment of tissue to bone is complete the fixation should be considered temporary and may not withstand weight bearing or other unsupported stresses. The suture anchor and suture are not intended to provide indefinite biomechanical integrity.
- For HEALICOIL Knotless PK devices, implantation of the HEALICOIL Knotless PK Suture Anchor requires preparation of the insertion site, hole preparation with the appropriate Smith & Nephew instruments is the recommended method of site preparation.
- For HEALICOIL Knotless PK devices, ensure the anchor placement is aligned with the hole. Proper alignment is essential for successful repair.

- Use of excessive force during insertion can cause failure of the suture anchor or insertion device.
- Bone quality must be adequate to allow proper placement of suture anchor. Inadequate bone quality could result in loss of fixation or pullout of the suture anchor.
- Do not alter implant or instrumentation or performance may be compromised.
- For HEALICOIL Knotless PK Self-Tapping devices, interact with distal implant with care as it may puncture personal protective equipment and/or cause injury.
- A minimum of 2 suture limbs must be used or insertion and/or fixation performance may be compromised
- Use of suture quantities above the maximum specification (6 suture limbs or 4 tape limbs) or non-Smith & Nephew products may compromise insertion and/or fixation performance.
- Advancing the proximal implant before locking the sutures in the distal implant may compromise the inserter integrity, compromise insertion and/or fixation performance.
- After use, this device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.

## Adverse Reactions

- Loss of fixation or pullout of suture anchors can occur
- Mild inflammatory reaction
- Foreign body reaction
- Infection, both deep and superficial
- Allergic reaction

## MAGNETIC RESONANCE (MR) IMAGING

**MR** The HEALICOIL Knotless PK Suture Anchor includes a PEEK distal implant that is composed of non-metallic material. Non-clinical testing has demonstrated that the HEALICOIL Knotless PK Suture Anchor is MR safe.

**MR** The HEALICOIL Knotless PK Self-Tapping Suture Anchor includes a titanium distal implant that is composed of non-ferromagnetic titanium alloy (Ti-6Al-4V). Non-clinical testing has demonstrated that the HEALICOIL Knotless PK Self-Tapping Suture Anchor is MR Conditional; MR Conditional refers to an item that has been demonstrated to pose no known hazards in a specified MR environment under specified conditions of use. This implant can be scanned safely under the following conditions:

- Static magnetic field of 1.5-Tesla (1.5 T) or 3 Tesla (3.0 T).
- Spatial gradient field of up to:
- 2,000 G/cm (20 T/m) for 3 T systems.
- Maximum whole body averaged specific absorption rate (SAR) of:
- 2.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 1.5 T.
- 2.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 3 T.

### 1.5 T RF Heating

In non-clinical testing with body coil excitation, the HEALICOIL Knotless PK Self-Tapping Suture Anchor produced a temperature rise of 2.1°C or less at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg, as assessed by calorimetry for 15 minutes of continuous scanning in a 1.5 T Siemens Espree (MRC30732) MR scanner with SYNGO MR B19 software.

### 3 T RF Heating

In non-clinical testing with body coil excitation, the HEALICOIL Knotless PK Self-Tapping Suture Anchor produced a temperature rise of less than 3.2°C at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg, as assessed by calorimetry for 15 minutes of continuous scanning in a 3 T Siemens Trio (MRC20587) MR scanner with SYNGO MR A35 4VA35A software.

**Caution:** The RF heating behavior does not scale with static field strength. Devices which do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

### 3 T MR Artifact

The image artifact extends approximately 1.1 cm from the HEALICOIL Knotless PK Self-Tapping Suture Anchor when scanned in nonclinical testing using a spin-echo sequence in a 3 T Siemens Trio (MRC20587) MR scanner with SYNGO MR A35 4VA35A software and body coil.

### Intended User Groups

- Orthopedic surgeons

### Target Patient Population

- Patients with soft tissue damage that requires reattachment to the bone that do not fall into any of the contraindicated groups.

### Instructions for Use

 **WARNING:** Use of the appropriate Smith & Nephew instrumentation is recommended to prepare the insertion site for the HEALICOIL Knotless PK Suture Anchor.

- Using routine sterile technique, remove the suture anchor from the packaging.  
**Note:** Both knobs located at the proximal end of the handle are preset and should not be turned until appropriate step is reached in the instructions for use.
- Use the recommended Smith & Nephew hole preparation device to prepare the insertion site.
  - Smith & Nephew disposable and reusable awls and dilators specific to the suture anchor are sold separately.
  - When a larger hole is required, as in the case of hard bone, reusable drills specific to the suture anchor are sold separately.
  - Refer to table below for appropriate awls, dilators, and drills.

Suture Anchor	For Soft Bone	For General Use	For Hard Bone
REF 72205137 HEALICOIL Knotless PK Suture Anchor, 5.0 mm	REF 72202621 3.8 mm Tapered Awl, Disposable or REF 72201915 3.8 mm Tapered Awl, Reusable	REF 72203710 5.5 mm Reusable Threaded Dilator or REF 72203952 5.5 mm Disposable Threaded Dilator	REF 72205308 4.75 mm Spade Tip Drill with REF 72203710 5.5 mm Reusable Threaded Dilator or REF 72203952 5.5 mm Disposable Threaded Dilator
REF 72205138 HEALICOIL Knotless PK Suture Anchor, 5.0 mm, Self-Tapping	No hole Preparation Required		REF 72203952 5.5 mm Disposable Threaded Dilator

3. For HEALICOIL Knotless PK only,  
For hard bone, place the distal tip of the spade drill onto the bone at the desired implantation site. With the drill powered on, advance the spade drill until the shoulder on the drill bit reaches the bone surface and then remove.
  - a. Place the distal tip of the awl or dilator onto the bone at the desired implantation site.
  - b. Hold the hole preparation instrument in place and use a mallet to lightly strike the proximal end of the awl or thread the dilator to prepare the insertion site. The proper hole depth is achieved when the laser depth mark is fully inserted into the bone.
  - c. Hold the awl/dilator in axial alignment and remove from the insertion site.
4. Thread sutures or tape through the suture threader loop, remove the suture threader from the shaft of the insertion device, and pull to feed the sutures through the distal implant.

Maximum suture limb load:

Suture type	Maximum per pass		Total
	1st pass	2nd pass	
ULTRABRAID #2	4 limbs	2 limbs	6 limbs
ULTRATAPE/ MINITAPE	2 limbs	2 limbs	4 limbs

**Note:** While passing suture it is recommended to support the proximal end of the distal implant.

5. Slide the loaded suture anchor either percutaneously or through a cannula, depending on technique preference. Minimize suture entanglement by removing excess suture slack and advance the anchor to the insertion site.
 

**Note:** Do not attempt to tension the tissue at this time.

**Note:** Ensure the sutures are not twisted around the suture anchor or inserter.
6. For HEALICOIL Knotless PK, establish and maintain axial alignment of the suture anchor to the prepared insertion site and place the tip of the distal implant into the prepared hole. Push or lightly tap the inserter until the distal portion of the proximal implant's first thread is positioned partially within the prepared hole.
 

For HEALICOIL Knotless PK Self-Tapping Suture Anchor, establish and maintain perpendicular alignment of the suture anchor to the insertion site. Using a mallet, lightly tap the inserter handle until the distal portion of the proximal implant's first thread is positioned partially within the bone.
7. Apply tension and lock the suture.
  - a. Maintain slight downward pressure on the inserter handle. Manually pull sutures individually to apply the desired tension. Tensioned sutures may be locked in the suture cleats.
  - b. Maintain slight downward pressure on the inserter and locate the smaller, most proximal knob on the inserter handle. Rotate knob clockwise until a minimum of one click is heard.
  - c. Check the repair. If the final tension on the suture is not adequate, rotate the smaller, most proximal knob counterclockwise until the suture can slide easily. Repeat steps "a" and "b" to apply appropriate tension and lock the suture. Repeat as necessary to achieve the desired repair tension.

**Note:** Verify the proximal anchor's first thread is positioned partially within the bone. If not, lightly tap the inserter handle until starting thread is seated in bone.

- Continue to maintain axial alignment and turn the larger, more distal knob clockwise to advance the proximal implant into the bone until the laser mark is flush with the surface of the cortical bone.
- If suture limbs have been cleated, uncleat them. Disengage the suture anchor from the inserter by pulling straight back or tapping the distal end of the handle with a mallet. Do not wiggle the handle to disengage it from the suture anchor. Discard the inserter.

## Warranty

For single use only. This product is warranted to be free from defects in material and workmanship. Do not reuse.

## For Further Information

If further information on this product is needed, contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized representative. This product may be covered by one or more US patents. See '[smith-nephew.com/patents](http://smith-nephew.com/patents)' for details.

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©2020 Smith & Nephew, Inc. Reservados todos los derechos.
Smith & Nephew Orthopaedics GmbH
Smith & Nephew, Inc.