Smith+Nephew

LEGION[◇] CONCELOC[◇]

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Bibliografia

1. Watson J, Jordan J. LEGION™ Primary Knee System for total knee arthroplasty: Design rationale and early results. *Bone & Joint Science*. 2015;5(1):1-8.
2. Nadaud M, Snyder MA, Brown D, McPherron A. New trends and early clinical outcomes with a modern knee revision system. *Orthopedics*. 2007;30(8):97.
3. Bieger R, Kappe T, Fraitzl CR, Reichel H. The aetiology of total knee arthroplasty failure influences the improvement in knee function. *Arch Orthop Trauma Surg*. 2013;133(2):237-241.
4. Christensen DD, Klement MR, Moschetti WE, Fillingham YA. Current Evidence-based Indications for Modern Noncemented Total Knee Arthroplasty. *J Am Acad Orthop Surg*. 2020;28(20):823-829.
5. Yayac M, Harrer S, Hozack WJ, Parvizi J, Courtney M. The use of cementless components does not significantly increase procedural costs in total knee arthroplasty. *J Arthroplasty*. 2020;35:407-712.
6. Behery OA, Clair AJ, Long WJ, Deshmukh A, Schwarzkopf R. Cementless primary total knee arthroplasty will this be the future? *Bull Hosp Jt Dis*. 2013;79(1):6-10.
7. Donaldson AJ, Harper NJ, Kenny NW, Thomson HE. Bone cement implantation syndrome. *Br J Anaesth*. 2009;102:12-22.
8. Bieger R, Kappe T, Fraitzl CR, Reichel H. The aetiology of total knee arthroplasty failure influences the improvement in knee function. *Arch Orthop Trauma Surg*. 2013;133(2):237-241.
9. Mufty S, Vandenuecker H, Bellemans J. The influence of leg length difference on clinical outcome after revision TKA. *Knee*. 2014;21(2):424-427.
10. National Joint Registry for England, Wales, Northern Ireland and the Isle of Man. Herefordshire, UK 2020. 17th Annual Report.
11. Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR) 2020. Hip, Knee & Shoulder Arthroplasty Annual Report.
12. McCalden RW, Hart GP, MacDonald SJ, Naudie DD, Howard JH, Bourne RB. Clinical Results and Survivorship of the GENESIS II Total Knee Arthroplasty at a Minimum of 15 Years. *J Arthroplasty*. 2017;32(7):2161-2166.
13. Bhandari M, Pascale W, Sprague S, Pascale V. The Genesis II in primary total knee replacement: a systematic literature review of clinical outcomes. *Knee*. 2012;19(1):8-13.
14. ODEP. ODEP Product Rating: GENESIS II OXINIUM PS GenII HF Resurfacing. Available at: <https://www.odep.org.uk/product.aspx?pid=9715>. Accessed 09 June, 2021.
15. ODEP. ODEP Product Rating: GENESIS II CR Cemented OXINIUM. Available at: <https://www.odep.org.uk/product.aspx?pid=9532>. Accessed 09 June, 2021.
16. Papisoulis E, Karachalios T. A 13- to 16-year clinical and radiological outcome study of the genesis II cruciate retaining total knee arthroplasty with an oxidised zirconium femoral component. *Knee*. 2019;26(2):492-499.
17. Lee J, Wang S, Kim K. Is there a difference in joint line restoration in revision Total knee arthroplasty according to prosthesis type? *BMC Musculoskelet Disord*. 2018;19(1):382.
18. Innocenti M, Matassi F, Carulli C, Soderi S, Villano M, Civinini R. Joint line position in revision total knee arthroplasty: the role of posterior femoral off-set stems. *The Knee*. 2013;20(6):447-450.
19. Gil-Martínez P, Sanz P, López-Torres I, Arnal-Burró J, Chana F, Vaquero J. Influence of the cause of the revision on the outcome after revision knee arthroplasty with condylar constrained implant. *Rev Esp Cir Ortop Traumatol (English Edition)*. 2016;60(3):184-191.
20. Stockwell KD, Malleck S, Gascoyne TC, Turgeon TR. Clinical and radiographic outcomes of a hybrid fixation revision total knee arthroplasty system at short to mid-term follow-up. *Knee*. 2019;26(1):240-249.
21. Smith + Nephew 2015. Orthopaedic Research Report OR-14-091A.
22. S&N Material Specification 0071120.
23. Kienapfel H, Sprey C, Wilke A, Griss P. Implant fixation by bone ingrowth. *J Arthroplasty*. 1999;14(3):355-368.
24. Bobyn JD, Pilliar RM, Cameron HU, Weatherly GC. The optimum pore size for the fixation of porous-surfaced metal implants by the ingrowth of bone. *Clin Orthop Relat Res*. 1980;150:263-270.
25. Sidambe A. Biocompatibility of Advanced Manufactured Titanium Implants-A Review. *Materials*. 2014;7(12):8168-8188.
26. Heuer DA. Structural comparison of advanced and conventional bone ingrowth technologies. Poster presented at Annual Meeting of the ORS, Long Beach CA, 2011. Poster No. 1034.
27. Smith+Nephew 2019. TM-19-081.
28. Stryker. Tritanium primary acetabular shells. SODTR-SS. 2008.
29. Bertollo N, Matsubara M, Shinoda T, et al. Effect of surgical fit on integration of cancellous bone and implant cortical bone shear strength for a porous titanium. *J Arthroplasty*. 2011;26(7):1000-1007.
30. Tong W, Smith B, Henninger D, Kane R, Muench T. Bone ingrowth of RTT porous structure. Paper presented at Annual Meeting of the ORS, 2020. Paper No.0337.
31. Morrison M, Woodard E. Comparison of Friction between Traditional Porous Coatings and an Additive Manufactured Porous Structure. Poster presented at Annual Meeting of the ORS, 2017. Poster No. 1095.
32. Bobyn JD, Hacking SA, Chan SP, et al. Characterization of a new porous tantalum biomaterial for reconstructive orthopaedics. Exhibit presented at 1999 Annual Meeting of the AAOS, Anaheim, CA.
33. Morrison M, Taylor ME, Weaver C. Characterization of a Reticulated, Randomized Porous Structure Produced by Additive Manufacturing. Poster presented at Annual Meeting of the ORS, 2016. Poster No. 2004.
34. Morgan EF, Keaveny TM. Dependence of yield strain of human trabecular bone on anatomic site. *J Biomech*. 2001;34(5):569-577.
35. Smith+Nephew 2015. OR-15-114.
36. Zardiackas LD, Parsell DE, Dillon LD, Mitchell DW, Nunnery LA, Poggie R. Structure, metallurgy, and mechanical properties of a porous tantalum foam. *J Biomed Mater Res B*. 2001;58(2):180-187.
37. Smith+Nephew 2015. OR-15-022A.
38. Fenwick S, Wilson D, Williams M, et al. A load bearing model to assess osseointegration of novel surfaces - A pilot study. *World Biomater Cong*, Amsterdam, NL, May 28-Jun 1, 2008, 233.
39. Smith & Nephew 2015. NCS249.
40. Khodadadyan-Klostermann C, Seebach M, Taylor MR, Duda GN, Hass NP. Distribution of bone mineral density with age and gender in the proximal tibia. *Clinical Biomechanics*. 2004;19:370-376.
41. Abdel MP, Gesell MW, Hoedt CW, Meyers KN, Wright TM, Haas SB. Polished Trays Reduce Backside Wear Independent of Post Location in Posterior-stabilized TKAs. *Clinical Orthopaedics and Related Research*. 2014;472(8):2477-2482.
42. Dolan MM, Kelly NH, Nguyen JT, Wright TM, Haas SB. Implant design influences tibial post wear damage in posterior-stabilized knees. *Clin Orthop Relat Res*. 2011;469(1):160-167.
43. Sisko ZW, Teeter MG, Lanting BA, et al. Current Total Knee Designs: Does Baseplate Roughness or Locking Mechanism Design Affect Polyethylene Backside Wear? *Clin Orthop Relat Res*. 2017;475(12):2970-2980.
44. Teeter MG, Thoren J, Yuan X, et al. Migration of a cemented fixed-bearing, polished titanium tibial baseplate (Genesis II) at ten years. *Bone Joint J*. 2016;98-B:616-621.
45. Hicks CA, Noble P, Tullos H. The Anatomy of the Tibial Intramedullary Canal. *Clin Orthop Relat Res*. 1995;321:111-116.
46. Westich, GH, Laskin RS, Haas SB, Sculco TD. Resection Specimen Analysis of Tibial Coverage in Total Knee Arthroplasty. *Clin Orthop Rel Res*. 1994;309:163-175.
47. Shah S, Agarwal N, Jain A. MRI Based Comparison of Tibial Bone Coverage by Five Knee Prosthesis: Anthropometric Study in Indians. *J Arthroplasty*. 2015;30:1643-1646.
48. Martin S, Saurez A, Ismaili S, et al. Maximizing Tibial Coverage Is Detrimental to Proper Rotational Alignment. *Clin Orthop Relat Res*. 2014;472:121-125.
49. Westrich GH, Haas SB, Insall JN, Frachie A. Resection specimen analysis of proximal tibial anatomy based on 100 total knee arthroplasty specimens. *J Arthroplasty*. 1995;10:47-51.
50. Hartel MJ, Loosli Y, Gralla J, et al. The mean anatomical shape of the tibial plateau at the knee arthroplasty resection level: an investigation using MRI. *The Knee*. 2009;16:452-457.
51. Westich GH, Agulnick MA, Laskin RS, Haas SB, Sculco TP. Current analysis of tibial coverage in total knee arthroplasty. *The Knee*. 1997;4:87-91.
52. Smith+Nephew 2020. OR-20-149.
53. ASTM F1800-19.
54. Smith+Nephew 2020. OR-20-157.
55. Smith+Nephew 2020. OR-20-148.
56. Smith+Nephew 2020. OR-20-153.
57. Grau LC, Ong AC, Restrepo S, et al. Survivorship, Clinical and Radiographic Outcomes of a Novel Cementless Metal-Backed Patella Design. *J Arthroplasty*. 2021;1-6.
58. Smith+Nephew 2020. OR-20-165.
59. Smith+Nephew 2020. OR-20-169.
60. Hardeman F, Vandenuecker H, Van Lauwe J, Bellemans J. Cementless Total Knee Arthroplasty with Profix: A 8 to 10 Year Follow-up Study. *The Knee*. 2006;13(6): 419-421.
61. Whiteside LA, Vigano R. Young and heavy patients with a cementless TKA do as well as older and lightweight patients. *Clin Orthop Relat Res*. 2007;464:93-98.
62. Kaper BP, Woolfrey M, Bourne RB. The Effect of Built-In External Rotation on Patellofemoral Tracking in the GENESIS II Total Knee Arthroplasty. *J Arthroplasty*. 2000; 15:964-969.
63. Demcoe AR, Bohm ER, Hedden DR, et al. Does oxidized zirconium make a difference? Midterm cohort survivorship of symmetric posterior condyle posterior-stabilized total knee arthroplasty. *Can J Surg*. 2019;62:118-122.
64. ODEP. Orthopaedic Data Evaluation Panel. Available at <http://www.odep.org.uk>. Accessed 30 June 2021.
65. Stockwell KD, Gascoyne TC, Singh M, et al. Survivorship of constrained polyethylene inserts in primary total knee replacements. *Knee*. 2020;27:1343-1348.

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©2021 Smith & Nephew, Inc.