

## Flexion/Extension Gap Balancing



# Balancer components

Balancer Shaft and Handle  
Shaft: Part #71933037  
Handle: Part #71674076



Balancer Body  
Part #71933036



Balancer Tibial Paddle  
Part #71933038



Tibial Paddle assembled onto the body



Shaft assembled into the body and tibial paddle



Femoral Paddle  
Part #71933039




Balancer in Extension with the Femoral Paddle



JOURNEY® Sizing Guide and Stylus  
Sizing Guide: Part #71933040  
Stylus: Part #71933041



Balancer in Flexion with Sizing Guide and Stylus



## Introduction

Correct mechanical alignment, sizing, and soft tissue balancing of the knee may produce a more functional knee replacement. Improper alignment, sizing, or balancing may cause flexion instability, flexion/extension mismatch, or residual deformity. These are common causes of patient pain, complaints, and revisions.

One method to improve alignment, sizing, and balancing is a Gap Balancing technique using a tensioning device. It's important to consider that TKA is a bone AND soft tissue operation. Often soft tissues are pathologic and respond variably to releases. Using a tensioner to set the size and rotation of the femoral component aids in balancing the flexion gap by adjusting bone resections rather than potentially requiring additional soft tissue releases.


A common misconception is that the utilization of a gap balancer may be cumbersome and less efficient to traditional monolithic balancing instruments. The goal of this surgical technique is to demonstrate the robustness, efficiency, and accuracy of utilizing a soft-tissue gap balancer to achieve a well-aligned TKA.

## 1. Perform distal femoral and proximal tibia resection

per recommended surgical technique

JOURNEY® II TKA - #00829

**Note:** Leave pins in order to recut if needed



## 2. Balance Extension Gap with Spacer Blocks/Tensioner

### a. Check Extension Gap w/ Spacer Block

Extension Gap	Next Step
Good	Remove all pins
Loose	Insert Larger Spacer Block until good
Tight	Balance through initial Varus/Valgus Releases
Tight after V/V Releases	Recut Bone

### b. Once gross extension balance obtained by spacer block, Insert Tensioner

#### i. Tension knee to match the insert thickness found with the Spacer Block

**Note:** Note tension on torque wrench required to obtain correct insert thickness

#### ii. Fine tune extension balance with needling techniques in order to read a "0" setting on guide<sup>1</sup>

1. Bellemans J. Multiple Needle Puncturing: Balancing the Varus Knee. ORTHOPEDICS. 1; 32: 693. doi: 10.3928/10477447-20110714-48



**Nota Bene:** The technique description herein is made available to the healthcare professional to illustrate the authors' suggested treatment for the uncomplicated procedure. In the final analysis, the preferred treatment is that which addresses the needs of the patient.

### 3. Balance Flexion Gap

- Change out Tensioner device for use in flexion
- Tension knee in 90° of flexion to same torque wrench setting recorded in extension (make sure tensioner is flush on Distal cut)
- Size femoral component with anterior stylus on lateral trochlea
- Loosen set screw and shift mechanism anterior to posterior as needed to match flexion gap to extension gap and finalize size of femoral component

#### Femoral Rotation Guidance:



**Varus Knee:** Balancer should automatically externally rotate femoral component block

**Valgus Knee:** CAUTION! Avoid internal rotation of femoral component cutting block by considering releasing popliteus insertion on lateral femur or manually overriding rotation prior to drilling holes



- Lock the Sizing Guide and drill Holes for 5-in-1 Cutting Block

**Note:** Drill the 3-10 Holes for all sizes (1-10)

The drill holes for 1-2 are not used for JOURNEY® II TKA.



### 4. Insert 5-in-1 Cutting Block and finish femoral cuts



## Contributing Surgeons


Mark Freeman, MD (Chattanooga, TN)

G van Hellemond, MD (Nijmegen, Netherlands)

Matt Nadaud, MD (Knoxville, TN)

## Reference JOURNEY II TKA Surgical Technique 00829

The following technique guide was prepared under the guidance of Dr. Mark Freeman, Prof. G van Hellemond and Dr. Matt Nadaud under close collaboration with each physician. It contains a summary of medical techniques and opinions based upon their training and expertise in the field, along with their knowledge of Smith & Nephew products. It is provided for educational and informational purposes only. Smith & Nephew does not provide medical advice and it is not intended to serve as such. It is the responsibility of the treating physician to determine and utilize the appropriate products and techniques according to their own clinical judgment for each of their patients. For more information on the products in this surgical technique, including indications for use, contraindications, effects, precautions and warnings, please consult the products' Instructions for Use (IFU)





Supporting healthcare professionals for over 150 years

**Smith & Nephew, Inc.**  
1450 Brooks Road  
Memphis, TN 38116  
USA

[www.smith-nephew.com](http://www.smith-nephew.com)

Telephone: 1-901-396-2121  
Information: 1-800-821-5700  
Orders/inquiries: 1-800-238-7538