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
A comparison of rollback ratio between bicruciate substituting total knee arthroplasty and Oxford unicompartmental knee arthroplasty¹²

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



Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics¹³

Lutes W and Fitch D. Presented at: 39th SICOT Orthopaedic World Congress; 2018; Montreal, Canada

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
*Fabbricato da DePuy Synthes Orthopaedics, Warsaw, IN, USA

The bicruciate substituting knee design and initial experience¹⁴

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
A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY II bi-cruciate stabilizing knee system and total hip arthroplasty patients¹⁵

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

[Versione integrale](#)

Midterm performance of a guided-motion bicruciate-stabilized total knee system: results from the international study of over 2000 consecutive primary total knee arthroplasties¹⁶

Harris AI, Christen B, Malcorps JJ, O'Grady CP, Kopjar B, Sensiba PR, Vandenneucker H, Huang BK, Cates HE, Hur J, Marra DA. *J Arthroplasty*. 2019; 34:S201-S208

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



In vivo knee kinematics: how important are the roles of femoral geometry and the cruciate ligaments?¹⁷

Smith LA, Nachtrab J, LaCour M, Cates H, Freeman MG, Komistek RD. *J Arthroplasty*. 2020: doi: <https://doi.org/10.1016/j.arth.2020.10.020>

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Boese K, et al.	Early clinical and patient-reported results of a bi-cruciate retaining total knee implant: six-month results of a prospective multicentre study of 149 primary TKAs	Harris AI, et al.	Guided motion total knee arthroplasty (TKA) in patients with BMI of 40kg/m ² or more: results from the international multicentre study of 2,059 primary TKAs with up to 6 years follow-up
Christen B, Kojar B.	Second-generation bi-cruciate stabilized total knee system has a lower reoperation and revision rate than its predecessor	Heir S, et al.	Clinical and functional outcomes of a second generation guided motion total knee arthroplasty system: two-year results of a prospective multicentre study
Di Benedetto P, et al.	Pre-operative and post-operative kinematic analysis in total knee arthroplasty	Hommel H, Wilke K.	Good Early Results Obtained with a Guided-Motion Implant for Total Knee Arthroplasty
Grieco TF, et al.	In vivo kinematic comparison of a bicruciate stabilized total knee arthroplasty and the normal knee using fluoroscopy	Hino K, et al.	Bi-cruciate substituting total knee arthroplasty provides varus-valgus stability throughout the midflexion range
Iriuchishima T, Ryu K.	A comparison of rollback ratio between bicruciate substituting total knee arthroplasty and Oxford unicompartmental knee arthroplasty	Hyodo K, et al.	Systematic review & meta-analysis of total knee arthroplasty kinematics ⁹
Harris AI, et al.	Short-term safety and effectiveness of a second generation motion-guided total knee system	Inui H, et al.	The relationship between soft-tissue balance and intraoperative kinematics of guided motion total knee arthroplasty
Harris AI, et al.	Guided motion total knee arthroplasty system: five-year outcomes of the prospective multicentre US study	Inui H, et al.	The relationship between anteroposterior stability and medial-lateral stability of the bi-cruciate stabilized total knee arthroplasty
Harris AI, et al.	Midterm performance of a guided-motion bicruciate-stabilized total knee system: results from the international study of over 2000 consecutive primary total knee arthroplasties	Iriuchishima T and Ryu K.	Bicruciate substituting total knee arthroplasty improves stair climbing ability when compared with cruciate-retain or posterior stabilizing total knee arthroplasty
Harris AI, et al.	Guided motion total knee arthroplasty (TKA) system in younger patients has a lower revision rate than registry controls: results from the international multicenter study with up to 6 years follow-up	Ishibashi T, et al.	Kinematics of bicruciate and posterior stabilized total knee arthroplasty during deep knee flexion and stair climbing

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Ishida K, Shibayama N, Toda A, et al.	Comparison of intra-operative navigation-based kinematics between bi-cruciate-stabilised total knee arthroplasty (TKA) and conventional posterior stabilised TKA	Kosse NM, et al.	Minor adaptations in implant design bicruciate-substituted total knee system improve maximal flexion
Kaneko T, et al.	The influence of compressive forces across the patellofemoral joint on patient-reported outcome after bi-cruciate stabilized total knee arthroplasty	Lutes W, Fitch D.	Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics
Kaneko T, et al.	Bi-cruciate substituting total knee arthroplasty improved medio-lateral instability in mid-flexion range	Mayman DJ, et al.	Hospital related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patients
Kaneko T, et al.	The influence of compressive forces across the patellofemoral joint on patient-reported outcome after bi-cruciate stabilized total knee arthroplasty	Moewis P, et al.	Retention of posterior cruciate ligament alone may not achieve physiological knee joint kinematics after total knee arthroplasty: a retrospective study
Kono K, et al.	In vivo kinematics of bicruciate-retaining total knee arthroplasty with anatomical articular surface under high-flexion conditions	Murakami, et al.	In vivo kinematics of gait in posterior-stabilized and bicruciate-stabilized total knee arthroplasties using image matching techniques
Kono K, et al.	Bicruciate-stabilised total knee arthroplasty provides good functional stability during high-flexion weight-bearing activities	Murakami K, et al.	Knee kinematics in bicruciate stabilized total knee arthroplasty during squatting and stair-climbing activities
Kopjar B, et al.	Clinical and functional outcomes of JOURNEY II CR total knee system. Interim results of an ongoing, prospective, multicenter study	Murakami K, et al.	Preoperative tibial mechanical axis orientation and articular surface design influence on the coronal joint line orientation relative to the ground during gait after total knee arthroplasties
		Nodzo, SR; et al.	The bicruciate substituting knee design and initial experience
		Oikonomidis L, et al.	The JOURNEY bicruciate knee replacement: design modifications yield better early functional results and reduce complications

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Patel AR, Delhougne G.	Hospital related clinical and economical outcomes of two premium knee systems in total knee arthroplasty (TKA) patients
Salzmann M, et al.	Does postoperative mechanical axis alignment have an effect on clinical outcome of primary total knee arthroplasty?
Smith LA, et al.	In vivo knee kinematics: how important are the roles of femoral geometry and the cruciate ligaments?
Snyder MA, et al.	A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY[®] II bi-cruciate stabilizing knee system and total hip arthroplasty patients
Takubo A, et al.	Comparison of muscle recovery following bicruciate substituting versus posterior stabilized total knee arthroplasty in an Asian population
Tomite T, et al.	Gait analysis of conventional total knee arthroplasty and bicruciate stabilized total knee arthroplasty using a triaxial accelerometer
West JA, et al.	Clinical outcomes and patient satisfaction after total knee arthroplasty; a follow-up of the first 50 cases by a single surgeon
Zambianchi F, et al.	Changes in total knee arthroplasty design affect in-vivo kinematics in a redesigned total knee system: A fluoroscopy study

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
Hyodo K, et al. 2020
Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty²¹

Boese K, et al. 2019
Early clinical and patient-reported results of a bi-cruciate retaining total knee implant: six-month results of a prospective multicentre study of 149 primary TKAs²⁴

Di Benedetto P, et al. 2019
Pre-operative and post-operative kinematic analysis in total knee arthroplasty. A pilot study²²

Heir S, et al. 2019
Clinical and functional outcomes of a second-generation guided motion total knee arthroplasty system: two-year results of a prospective multicentre study²⁵

Patel AR, et al. 2019
Hospital related clinical and economical outcomes of two premium knee systems in total knee arthroplasty (TKA) patients²⁶

West JA, et al. 2019
Clinical outcomes and patient satisfaction after total knee arthroplasty: a follow-up of the first 50 cases by a single surgeon²⁷

Christen B, et al. 2018
Second-generation bi-cruciate stabilized total knee system has a lower reoperation and revision rate than its predecessor²⁸

Harris AI, et al. 2018
Short-term safety and effectiveness of a second-generation motion-guided total knee system²⁹

Lutes W, et al. 2018
Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics¹³

Mayman DJ, et al. 2018
Hospital related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patients²³

Nodzo SR, et al. 2018
The Bicruciate Substituting Knee Design and Initial Experience¹⁴

Snyder MA, et al. 2018
A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY^o II bi-cruciate stabilizing knee system and total hip arthroplasty patients¹⁵

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Hyodo K, et al. 2020	Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty ²¹	Kono K, et al. 2019	Bicruciate-stabilised total knee arthroplasty provides good functional stability during high-flexion weight-bearing activities ⁴¹
Kaneko T, et al. 2020	The influence of tibiofemoral joint forces on patient-reported outcome measurements after bicruciate stabilized total knee arthroplasty ³⁵	Kono K, et al. 2019	In vivo kinematics of bicruciate-retaining total knee arthroplasty with anatomical articular surface under high-flexion conditions ⁴⁰
Moewis P, et al. 2020	Retention of posterior cruciate ligament alone may not achieve physiological knee joint kinematics after total knee arthroplasty: a retrospective study ³⁶	Kopjar B, et al. 2019	Clinical and functional outcomes of JOURNEY ^o II CR total knee system. Interim results of an ongoing, prospective, multicenter study ⁴²
Oikonomidis L, et al. 2020	The Journey bicruciate knee replacement: design modifications yield better early functional results and reduce complications ³⁷	West JA, et al. 2019	Clinical outcomes and patient satisfaction after total knee arthroplasty: a follow-up of the first 50 cases by a single surgeon ²⁷
Di Benedetto P, et al. 2019	Pre-operative and post-operative kinematic analysis in total knee arthroplasty. A pilot study ²²	Grieco TF, et al. 2018	In vivo kinematic comparison of a bicruciate stabilized total knee arthroplasty and the normal knee using fluoroscopy¹¹
Heir S, et al. 2019	Clinical and functional outcomes of a second-generation guided motion total knee arthroplasty system: two-year results of a prospective multicentre study ²⁵	Harris AI, et al. 2018	Short-term safety and effectiveness of a second-generation motion-guided total knee system ²⁹
Inui H, et al. 2019	The relationship between soft-tissue balance and intraoperative kinematics of guided motion total knee arthroplasty ³⁸	Hino K, et al. 2018	Bi-cruciate substituting total knee arthroplasty provides varus- valgus stability throughout the midflexion range ⁴³
Iriuchishima T, et al. 2019	Bicruciate substituting total knee arthroplasty improves stair climbing ability when compared with cruciate-retain or posterior stabilizing total knee arthroplasty ³⁹	Inui H, et al. 2018.	The relationship between anteroposterior stability and medial-lateral stability of the bi-cruciate stabilized total knee arthroplasty ⁴⁴

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Iriuchishima T, et al. 2018	A comparison of rollback ratio between bicruciate substituting total knee arthroplasty and Oxford unicompartmental knee arthroplasty¹²	Zambianchi F, et al. 2018	Changes in total knee arthroplasty design affect in-vivo kinematics in a redesigned total knee system: A fluoroscopy study ⁴⁷
Kosse NM, et al. 2018	Improved maximal flexion after minor adaptations in implant design bicruciate-substituted total knee arthroplasty ³²	Hommel H, et al. 2017	Good Early Results Obtained with a Guided-Motion Implant for Total Knee Arthroplasty: A Consecutive Case Series ⁴⁸
Lutes W, et al. 2018	Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics¹³	Ishida K, et al. 2017	Comparison of intra-operative navigation-based kinematics between bi-cruciate-stabilised total knee arthroplasty (TKA) and conventional posterior-stabilised TKA ⁴⁹
Murakami K, et al. 2018	In vivo kinematics of gait in posterior-stabilized and bicruciate-stabilized total knee arthroplasties using image-matching techniques ³⁴	Kaneko T, et al. 2017	Bi-cruciate substituting total knee arthroplasty improved medio-lateral instability in midflexion range ⁵⁰
Murakami K, et al. 2018	Knee kinematics in bi-cruciate stabilized total knee arthroplasty during squatting and stair-climbing activities ⁴⁵	Salzman M, et al. 2017	Does postoperative mechanical axis alignment have an effect on clinical outcome of primary total knee arthroplasty? A retrospective cohort study ⁵¹
Murakami K, et al. 2018	Preoperative tibial mechanical axis orientation and articular surface design influence on the coronal joint line orientation relative to the ground during gait after total knee arthroplasties ⁴⁶	Takubo A, et al. 2017	Comparison of muscle recovery following bicruciate substituting versus posterior stabilized total knee arthroplasty in an Asian population ³³
Nodzo SR, et al. 2018	The Bicruciate Substituting Knee Design and Initial Experience¹⁴	Tomite T, et al. 2017	Gait analysis of conventional total knee arthroplasty and bicruciate stabilized total knee arthroplasty using a triaxial accelerometer ⁵²
Snyder MA, et al. 2018	A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY^o II bi-cruciate stabilizing knee system and total hip arthroplasty patients¹⁵		

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Moewis P, et al. 2020 Retention of posterior cruciate ligament alone may not achieve physiological knee joint kinematics after total knee arthroplasty: a retrospective study³⁶

Oikonomidis L, et al. 2020 The Journey bicruciate knee replacement: design modifications yield better early functional results and reduce complications³⁷

Boese K, et al. 2019 Early clinical and patient-reported results of a bi-cruciate retaining total knee implant: six-month results of a prospective multicentre study of 149 primary TKAs²⁴

Heir S, et al. 2019 Clinical and functional outcomes of a second-generation guided motion total knee arthroplasty system: two-year results of a prospective multicentre study²⁵

Kopjar B, et al. 2019 Clinical and functional outcomes of JOURNEY^o II CR total knee system. Interim results of an ongoing, prospective, multicenter study⁴²

West JA, et al. 2019 Clinical outcomes and patient satisfaction after total knee arthroplasty: a follow-up of the first 50 cases by a single surgeon²⁷

Harris AI, et al. 2018 Short-term safety and effectiveness of a second-generation motion-guided total knee system²⁹

Kaneko T, et al. 2018 The influence of compressive forces across the patellofemoral joint on patient-reported outcome after bi-cruciate stabilized total knee arthroplasty⁵⁶

Lutes W, et al. 2018 Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics¹³

Nodzo SR, et al. 2018 The Bicruciate Substituting Knee Design and Initial Experience¹⁴

Snyder MA, et al. 2018 A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY II bi-cruciate stabilizing knee system and total hip arthroplasty patients¹⁵

Kaneko T, et al. 2017 Bi-cruciate substituting total knee arthroplasty improved medio-lateral instability in midflexion range⁵⁰

Salzman M, et al. 2017 Does postoperative mechanical axis alignment have an effect on clinical outcome of primary total knee arthroplasty? A retrospective cohort study⁵¹

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Boese K, et al.
2019 Early clinical and patient-reported results of a bi-cruciate retaining total knee implant: six-month results of a prospective multicentre study of 149 primary TKAs²⁴

Harris AI, et al.
2019 **Midterm performance of a guided-motion bicruciate-stabilized total knee system: results from the international study of over 2000 consecutive primary total knee arthroplasties¹⁶**

Harris AI, et al.
2019 Guided motion total knee arthroplasty system: five-year outcomes of the prospective multicenter US study⁶¹

Harris AI, et al.
2019 Guided motion total knee arthroplasty (TKA) system in younger patients has a lower revision rate than registry controls: results from the international multicenter study with up to 6 Years follow-up⁶²

Harris AI, et al.
2019 Guided motion total knee arthroplasty (TKA) in patients with BMI of 40kg/m² or more: results from the international multicentre study of 2,059 primary TKAs with up to 6 years follow-up⁶³

Kopjar B, et al.
2019 Clinical and functional outcomes of JOURNEY^o II CR total knee system. Interim results of an ongoing, prospective, multicenter study⁴²

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
Patel AR et al.
2019

Hospital related clinical and economical outcomes of two premium knee systems in total knee arthroplasty (TKA) patients²⁶

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Guided motion total knee arthroplasty (TKA) in patients with BMI of 40kg/m² or more: results from the international multicentre study of 2,059 primary TKAs with up to 6 years follow-up⁶³

Harris A, O'Grady C, Sensiba PR, Cates HE, Christen B, Malcorps J, Kopjar B. Abstract presented at: EKS;2019; Valencia, Spain

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Guided motion total knee arthroplasty (TKA) system in younger patients has a lower revision rate than registry controls: results from the international multicenter study with up to 6 years follow-up⁶²

Harris A. O'Grady C, Sensiba PR, Cates HE, Christen B, Malcorps J, Kopjar B. Abstract presented at: EKS;2019; Valencia, Spain

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
Early clinical and patient-reported results of a bi-cruciate retaining total knee implant: six-month results of a prospective multicentre study of 149 primary TKAs²⁴

Boese K, MacDonald J, Huang W, Schwarzkopf R, Gerlinger T, Swank M, Huff T, Schinsky M, Amin N, Ast M, Ries M, Roche M, Jones J, Cooper H. Poster presented at: European Orthopaedic Research Society (EORS); 2019; Maastricht, the Netherlands

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Short-term safety and effectiveness of a second-generation motion-guided total knee system²⁹

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Guided motion total knee arthroplasty system: five-year outcomes of the prospective multicenter US study⁶¹

Harris AI, Luo TD, Lang JE, Snyder M, Haas S, Shields J, Cates H, Swank M, Kopjar B. Presented at: EFORT 2019; Lisbon, Portugal

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Clinical and functional outcomes of a second-generation guided motion total knee arthroplasty system: two-year results of a prospective multicentre study²⁵

Heir S, Catani F, van Hellemond G, Victor J, Jansen E, Zambianchi F, Kopjar B. number P37 presented at: European Knee Society; 2019; Valencia, Spain

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The relationship between soft-tissue balance and intraoperative kinematics of guided motion total knee arthroplasty³⁸

Inui H, Taketomi S, Yamagami R, Shirakawa N, Kawaguchi K, Tanaka S. *J Knee Surg.* 2019;32:91–96

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Bicruciate substituting total knee arthroplasty improves stair climbing ability when compared with cruciate-retain or posterior stabilizing total knee arthroplasty³⁹

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Clinical and functional outcomes of JOURNEY II CR total knee system. Interim results of an ongoing, prospective, multicenter study⁴²

Kopjar B, Archan P, MacDonald J, Gerlinger T, Chow J, Swank M, Geller J, Pifarre PT, Miles J, Shields J. presented at: ISTA; 2019; Toronto, Canada

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In vivo kinematics of bicruciate-retaining total knee arthroplasty with anatomical articular surface under high-flexion conditions⁴⁰

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Bicruciate-stabilised total knee arthroplasty provides good functional stability during high-flexion weight-bearing activities⁴¹

Kono K, Inui H, Tomita T, Yamazaki T, Taketomi S, Sugamoto K, Tanaka S. *Knee Surg Sports Traumatol Arthrosc.* 2019;27:2096-2103

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*Fabbricato da Stryker Orthopaedics, Mahwah, NJ USA

Clinical outcomes and patient satisfaction after total knee arthroplasty: a follow-up of the first 50 cases by a single surgeon²⁷

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Second-generation bi-cruciate stabilized total knee system has a lower reoperation and revision rate than its predecessor²⁸

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Bi-cruciate substituting total knee arthroplasty provides varus-valgus stability throughout the midflexion range⁴³

Hino K, Kutsuna T, Watamori K, Ishimaru Y, Kiyomatsu H, Shiraishi Y, Miura H. *Knee*. 2018;25:897-902

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
The relationship between anteroposterior stability and medial-lateral stability of the bi-cruciate stabilized total knee arthroplasty⁴⁴

Inui H, Taketomi S, Yamagami R, Kawaguchi K, Nakazaro K, Tanaka S. *Knee*. 2018;25:1247- 1253

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The influence of compressive forces across the patellofemoral joint on patient-reported outcome after bi-cruciate stabilized total knee arthroplasty⁵⁶

Kaneko T, Kono N, Mochizuki Y, Hada M, Sunakawa T, Ikegami H, Musha Y. *Bone Joint J.* 2018;100-B:1585-1591

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Improved maximal flexion after minor adaptations in implant design bicruciate-substituted total knee arthroplasty³²

Kosse NM, Heesterbeek PJC, Defoort KC, Wymenga AB, van Hellemond GG. presented at 19th EFORT Congress; 2018; Barcelona, Spain

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
Hospital related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patients²³

Mayman DJ, Patel AR, Carroll KM. Poster presented at: ISPOR Symposium;2018; Baltimore, Maryland, USA

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Knee kinematics in bi-cruciate stabilized total knee arthroplasty during squatting and stair-climbing activities⁴⁵

Murakami K, Hamai S, Okazaki K, Gondo H, Wang Y, Ikebe S, Higaki H, Shimoto T, Mizu-uchi H, Akasaki Y, Nakashima Y. *J Orthop.* 2018;15:650-654

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
In vivo kinematics of gait in posterior-stabilized and bicruciate-stabilized total knee arthroplasties using image-matching techniques³⁴

Murakami K, Hamai S, Okazaki K, Wang Y, Ikebe S, Higaki H, Shimoto T, Mizu-uchi H, Akasaki Y, Nakashima Y. *Int Orthop.* 2018;42:2753-2581

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
*Fabbricato da Zimmer Biomet, Warsaw, Indiana, USA

Preoperative tibial mechanical axis orientation and articular surface design influence on the coronal joint line orientation relative to the ground during gait after total knee arthroplasties⁴⁶

Murakami K, Hamai S, Okazaki K, Ikebe S, Higaki H, Shimoto T, Nakashima Y. *Knee Surg Sports Traumatol Arthrosc.* 2018;26:3368–3376

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
Guarigione

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Survivorship

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Changes in total knee arthroplasty design affect in-vivo kinematics in a redesigned total knee system: A fluoroscopy study⁴⁷

Zambianchi F, Fiacchi F, Lombardi V, Venturelli L, Marcobigi A, Giorgini A, Catani F. *Clin Biomech.* 2018;54:92-102

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Good early results obtained with a guided-motion implant for total knee arthroplasty: a consecutive case series⁴⁸

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Comparison of intra-operative navigation-based kinematics between bi-cruciate-stabilised total knee arthroplasty (TKA) and conventional posterior-stabilised TKA⁴⁹

Ishida K, Shibamura N, Toda A, Kodato K, Inokuchi T, Matsumoto T, Takayama K, Kuroda R, Kurosaka M. *Orthop Proceedings*. 2017;99B:Suppl 4

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*Fabbricato da Stryker Orthopaedics, Mahwah, NJ, USA

**Fabbricato da Zimmer Biomet, Warsaw, Indiana, USA

Bi-cruciate substituting total knee arthroplasty improved medio-lateral instability in midflexion range⁵⁰

Kaneko T, Kono N, Mochizuki Y, Hada M, Toyoda S, Musha Y. *J Orthop.* 2017;14:201-206

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
Does postoperative mechanical axis alignment have an effect on clinical outcome of primary total knee arthroplasty? A retrospective cohort study⁵¹

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Comparison of muscle recovery following bicruciate substituting versus posterior stabilized total knee arthroplasty in an Asian population³³

Takubo A, Ryu K, Iriuchishima T, Tokuhashi Y. *J Knee Surg.* 2017;30:725–729

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
*Fabbricato da Nakashima Medical, Japan

Gait analysis of conventional total knee arthroplasty and bicruciate stabilized total knee arthroplasty using a triaxial accelerometer⁵²

Tomite T, Saito H, Aizawa T, Kijima H, Miyakoshi N, Shimada Y. *Case Report Orthop.* 2016;2016:6875821

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Pre-operative and post-operative kinematic analysis in total knee arthroplasty. A pilot study²²

Di Benedetto P, Vidi D, Colombo, Buttironi MM, Cainero V, Causero A. *Acta Biomed.* 2019;90:91-97

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*Fabbricato da DePuy Synthes Orthopaedics, Warsaw, IN, USA

Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty²¹

Hyodo K, Kanamori A, Kadone H, Takahashi T, Kajiwara M, Yamazaki M. *Arthroplast Today*. 2020;6:338-342

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
The influence of tibiofemoral joint forces on patient-reported outcome measurements after bicruciate stabilized total knee arthroplasty³⁵

Kaneko T, Kono N, Mochizuki Y, Hada M, Toyoda S, Ikegami H, Musha Y. *J Orthop Surg.* 2020;28:1-2

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BCS

CR

Retention of posterior cruciate ligament alone may not achieve physiological knee joint kinematics after total knee arthroplasty: a retrospective study³⁶

Moewis P, Duda GN, Trepczynski A, Krahl L, Boese CK, Hommel H. *J Bone Joint Surg Am.* 2020; doi: 10.2106/JBJS.20.00024

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The JOURNEY bicruciate knee replacement: design modifications yield better early functional results and reduce complications³⁷

Oikonomidis L, Santini AJA, Davidson JS, Banks JV, Phillipson A, Pope J. *J Knee Surg.* 2020; doi: 10.1055/s-0040-1718599

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
Kinematics of bicruciate and posterior stabilized total knee arthroplasty during deep knee flexion and stair climbing⁶⁴

Ishibashi T, Tomita T, Yamazaki T, Tsuji S, Yoshikawa H, Sugamoto K. *J Orthop Res.* 2020: doi: 10.1002/jor.24773

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Smith+Nephew

BIBLIOGRAFIA

- Collins M, Lavigne M, Girard J, Vendittoli PA. Joint perception after hip or knee replacement surgery. *Orthop Traumatol Surg Res.* 2012;98:275–280.
- Noble PC, Gordon MJ, Weiss JM, et al. Does Total Knee Replacement Restore Normal Knee Function? *Clin Orthop Relat Res.* 2005;431:157–165.
- Scott CEH, Howie CR, MacDonald D, Biant LC. Predicting dissatisfaction following total knee replacement. *J Bone Joint Surg Am.* 2010;92-B(9):1253–1258.
- Lee A, Park J, Lee S. Gait analysis of elderly women after total knee arthroplasty. *J Phys Ther Sci.* 2015;27:591–595.
- Hyodo K, Kanamori A, Kadone H, Takahashi T, Kajiwara M, Yamazaki M. Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty. *Arthroplast Today.* 2020;6:338–342.
- Thomas AC, Stevens-Lapsley JE. Importance of attenuating quadriceps activation deficits after total knee arthroplasty. *Exerc Sport Sci Rev.* 2012;40:95–101.
- Birdsall PD, Hayes JH, Cleary R, Pinder IM, Moran CG, Sher JL. Health outcome after total knee replacement in the very elderly. *J Bone Joint Surg Br.* 1999;81:660–662.
- Evans JT, Evans JP, Walker RB, Blom AW, Whitehouse MR, Sayers A. How long does a hip replacement last? A systematic review and meta-analysis of case series and national registry reports with more than 15 years of follow-up. *Lancet.* 2019;393:647–654.
- Angerame MR, Holst DC, Jennings JM, Komistek RD, Dennis DA. Total knee arthroplasty kinematics. *J Arthroplasty.* 2019;34:2502–2510.
- Van Onsem S, Verstraete M, Van Eenoo W, Van Der Straeten C, Victor J. Are TKA Kinematics During Closed Kinetic Chain Exercises Associated With Patient-reported Outcomes? A Preliminary Analysis. *Clin Orthop Relat Res.* 2020;478:255–263.
- Grieco TF, Sharma A, Dessinger GM, et al. In vivo kinematic comparison of a bicruciate stabilized total knee arthroplasty and the normal knee using fluoroscopy. *J Arthroplasty.* 2018;33:565–571.
- Iriuchishima T, Ryu K. A Comparison of rollback ratio between bicruciate substituting total knee arthroplasty and Oxford unicompartmental knee arthroplasty. *J Knee Surg.* 2018;31:568–572.
- Lutes W and Fitch D. Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics. Presented at: 39th SICOT Orthopaedic World Congress, October 10–13, 2018, Montreal, Canada.
- Nodzo SR, Carroll KM, Mayman DJ. The Bicruciate Substituting Knee Design and Initial Experience. *Tech Orthop.* 2018;33:37–41.
- Snyder MA, Sympon A, Gregg J, Levit A. A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY II bi-cruciate stabilizing knee system and total hip arthroplasty patients. *Orthop Trauma Prosth.* 2018; 2018; 3:5–10.
- Harris AI, Christen B, Malcorps JJ, et al. Mid-term performance of a guided motion bicruciate stabilized total knee system; results from the international study of over 2,000 consecutive primary total knee arthroplasties. *J Arthroplasty.* 2019;34:5201–5208.
- Smith LA, Nachtrab J, LaCour M, Cates H, Freeman MG, Komistek RD. In vivo knee kinematics: how important are the roles of femoral geometry and the cruciate ligaments? *J Arthroplasty.* 2020; doi: <https://doi.org/10.1016/j.arth.2020.10.020>.
- Husain A, Lee GC. Establishing realistic patient expectations following total knee arthroplasty. *J Am Acad Orthop Surg.* 2015;23:707–713.
- Mutsuzaki H, Takeuchi R, Mataka Y, Wadano Y. Target range of motion for rehabilitation after total knee arthroplasty. *J Rural Med.* 2017;12:33–37.
- Hiyama Y, Wada O, Nakakita S, Mizuno K. Joint awareness after total knee arthroplasty is affected by pain and quadriceps strength. *Orthop Traumatol Surg Res.* 2016;102:435–439.
- Hyodo K, Kanamori A, Kadone H, Takahashi T, Kajiwara M, Yamazaki M. Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty. *Arthroplast Today.* 2020;6:338–342.
- Di Benedetto P, Vidi D, Colombo, Buttironi MM, Cainero V, Causero A. Pre-operative and post-operative kinematic analysis in total knee arthroplasty. A pilot study. *Acta Biomed.* 2019;90:91–97.
- Mayman DJ, Patel AR, Carroll KM. Hospital Related Clinical and Economic Outcomes of a Bicruciate Knee System in Total Knee Arthroplasty Patients. Poster presented at: ISFOR Symposium; May 19–23, 2018; Baltimore, Maryland, USA.
- Boese K, MacDonald J, Huang W, et al. Early clinical and patient-reported results of a bi-cruciate retaining total knee implant: six-month results of a prospective multicentre study of 149 primary TKAs. Poster presented at: European Orthopaedic Research Society (EORS); October 2–4, 2019; Maastricht, the Netherlands.
- Heir S, Catani F, van Hellemond G, et al. Clinical and functional outcomes of a second generation guided motion total knee arthroplasty system: two-year results of a prospective multicentre study. Abstract number P37 presented at: European Knee Society; May 2–3, 2019; Valencia, Spain.
- Patel AR, Delhougne G. Hospital related clinical and economical outcomes of two premium knee systems in total knee arthroplasty (TKA) patients. Poster presented at: ISPOR Symposium; May 18–22, 2019; New Orleans, Louisiana, USA.
- West JA, Scudday T, Anderson S, Amin NH. Clinical outcomes and patient satisfaction after total knee arthroplasty; a follow-up of the first 50 cases by a single surgeon. *J Int Med Res.* 2019;47:1667–1676.
- Christen B, Kopjar B. Second-generation bi-cruciate stabilized total knee system has a lower reoperation and revision rate than its predecessor. *Arch Orthop Trauma Surg.* 2018;138:1591–1599.
- Harris AI, Luo TD, Lang JE, Kopjar B. Short-term safety and effectiveness of a second generation motion-guided total knee system. *Arthroplast Today.* 2018;4:240–243.
- Matsuda S, Kawahara S, Okazaki K, Tashiro Y, Iwamoto Y. Postoperative alignment and ROM affect patient satisfaction after TKA. *Clin Orthop Relat Res.* 2013;471:127–133.
- Parcells BW, Tria AJ Jr. The Cruciate Ligaments in Total Knee Arthroplasty. *Am J Orthop.* 2016;45:153–160.
- Kosse NM, Heesterbeek PJ, Defoort KC, et al. Minor adaptations in implant design bicruciate-substituted total knee system improve maximal flexion. Poster presented at: 2nd World Arthroplasty Congress; 19–21 April, 2018; Rome, Italy.
- Takubo A, Ryu K, Iriuchishima T, Tokuhashi Y. Comparison of muscle recovery following bicruciate substituting versus posterior stabilized total knee arthroplasty in an Asian population. *J Knee Surg.* 2017;30:725–729.
- Murakami K, Hamai S, Okazaki K, et al. In vivo kinematics of gait in posterior-stabilized and bicruciate-stabilized total knee arthroplasties using image matching techniques. *International orthopaedics. Int Orthop.* 2018;42:2573–2581.
- Kaneko T, Kono N, Mochizuki Y, Hada M, Toyoda S, Ikegami H, Musha Y. The influence of tibiofemoral joint forces on patient-reported outcome measurements after bicruciate stabilized total knee arthroplasty. *J Orthop Surg.* 2020;28:1–2.
- Moewis P, Duda GN, Treczynski A, Krahl L, Boese CK, Hommel H. Retention of posterior cruciate ligament alone may not achieve physiological knee joint kinematics after total knee arthroplasty: a retrospective study. *J Bone Joint Surg Am.* 2020; doi: 10.2106/JBJS.20.00024.
- Oikonomidis L, Santini AJA, Davidson JS, Banks JV, Phillipson A, Pope J. The JOURNEY bicruciate knee replacement: design modifications yield better early functional results and reduce complications. *J Knee Surg.* 2020; doi: 10.1055/s-0040-1718599.
- Inui H, Taketomi S, Yamagami R, et al. The Relationship between Soft-Tissue Balance and Intraoperative Kinematics of Guided Motion Total Knee Arthroplasty. *J Knee Surg.* 2019;32:91–96.
- Iriuchishima T and Ryu K. Bicruciate substituting total knee arthroplasty improves stair climbing ability when compared with cruciate-retain or posterior stabilizing total knee arthroplasty. *Indian J Orthop.* 2019. DOI:10.4103/ortho.IJOrtho_392_18.
- Kono K, Inui H, Tomita T, Yamazaki T, Taketomi S, Tanaka S. In vivo kinematics of bicruciate-retaining total knee arthroplasty with anatomical articular surface under high-flexion conditions. *J Knee Surg.* 2019 Sep 9. [Epub ahead of print].
- Kono K, Inui H, Tomita T, et al. Bicruciate-stabilised total knee arthroplasty provides good functional stability during high-flexion weight-bearing activities. *Knee Surg Sports Traumatol Arthrosc.* 2019;27:2096–2103.
- Kopjar B, Archan P, MacDonald J, et al. Clinical and functional outcomes of JOURNEY II CR total knee system. Interim results of an ongoing, prospective, multicenter study. Abstract presented at: International Society for Technology in Arthroplasty Congress (ISTA); October 2–5, 2019; Toronto, Canada.
- Hino K, Kutsuna T, Watamori K, et al. Bi-cruciate substituting total knee arthroplasty provides varus–valgus stability throughout the midflexion range. *Knee.* 2018;25:897–902.
- Inui H, Taketomi S, Yamagami R, et al. The relationship between anteroposterior stability and medial-lateral stability of the bi-cruciate stabilized total knee arthroplasty. *Knee.* 2018;25:1247–1253.
- Murakami K, Hamai S, Okazaki K, et al. Knee kinematics in bicruciate stabilized total knee arthroplasty during squatting and stair-climbing activities. *J Orthop.* 2018;15:650–654.
- Murakami K, Hamai S, Okazaki K, et al. Preoperative tibial mechanical axis orientation and articular surface design influence on the coronal joint line orientation relative to the ground during gait after total knee arthroplasties. *Knee Surg Sports Traumatol Arthrosc.* 2018;26:3368–3376.
- Zambianchi F, Fiacchi F, Lombardi V et al. Changes in total knee arthroplasty design affect in-vivo kinematics in a redesigned total knee system: A fluoroscopy study. *Clin Biomech.* 2018;54:92–102.
- Hommel H, Wilke K. Good Early Results Obtained with a Guided-Motion Implant for Total Knee Arthroplasty: A Consecutive Case Series. *Open Orthop J.* 2017;11:51–56.
- Ishida K, Shibanuma N, Toda A, et al. Comparison of intra-operative navigation-based kinematics between bi-cruciate-stabilised total knee arthroplasty (TKA) and conventional posteriorstabilised TKA. *Orthopaedic proceedings.* 2017;99B:Suppl 4.
- Kaneko T, Kono N, Mochizuki Y, et al. Bi-cruciate substituting total knee arthroplasty improved medio-lateral instability in mid-flexion range. *J Orthop.* 2017;14:201–206.
- Salzmann M, Fennema P, Becker R, Hommel H. Does Postoperative Mechanical Axis Alignment Have an Effect on Clinical Outcome of Primary Total Knee Arthroplasty? A Retrospective Cohort Study. *Open Orthop J.* 2017;11:1330–1336.
- Tomite T, Saito H, Aizawa T, Kijima H, Miyakoshi N, Shimada Y. Gait Analysis of Conventional Total Knee Arthroplasty and Bicruciate Stabilized Total Knee Arthroplasty Using a Triaxial Accelerometer. *Case Report Orthop.* 2016;2016:6875821.
- Lyu H, Wick E, Housman M, Freischlag J, Makary M. Patient satisfaction as a possible indicator of quality surgical care. *JAMA Surg.* 2013;148:362–367.
- Otani K, Waterman B, Faulkner K, Boslaugh S, Burroughs T, Dunagan W. Patient satisfaction: focusing on “excellent”. *J Healthc Manag.* 2009;54:93–102.
- Walker LC, Clement ND, Bardgett M, et al. The WOMAC score can be reliably used to classify patient satisfaction after total knee arthroplasty. *Knee Surg Sports Traumatol Arthrosc.* 2018;26:3333–3341.
- Kaneko T, Kono N, Mochizuki Y, et al. The influence of compressive forces across the patellofemoral joint on patient-reported outcome after bi-cruciate stabilized total knee arthroplasty. *Bone Joint J.* 2018;100-B:1585–1591.
- Inacio MCS, Paxton EW, Graves SE, Namba RS, Nemes S. Projected increase in total knee arthroplasty in the United States – an alternative projection model. *Osteoarthritis Cartilage.* 2017;25:1797–1803.
- Molloy IB, Brook IM, Wayne EM, Jevesevar DS. Effects of the length of stay on the cost of total knee and total hip arthroplasty from 2002 to 2013. *J Bone Joint Surg Am.* 2017;99:402–407.
- Hamilton DF, Howie CR, Burnett R, Simpson AH, Patton JT. Dealing with the predicted increase in demand for revision total knee arthroplasty; challenges, risks and opportunities. *Bone Joint J.* 2015;97-b:723–728.
- Delanois RE, Mistry JB, Gwam CU, Mohamed NS, Choksi US, Mont MA. Current epidemiology of revision total knee arthroplasty in the United States. *J Arthroplasty.* 2017;32:2663–2668.
- Harris A, Luo TD, Lang JE, et al. Guided motion total knee arthroplasty system: five-year outcomes of the prospective multicentre US study. Presentation number 2761 presented at: EFORT; June 5–7, 2019; Lisbon, Portugal.
- Harris A, O’Grady C, Sensiba PR, et al. Guided motion total knee arthroplasty (TKA) system in younger patients has a lower revision rate than registry controls: results from the international multicenter study with up to 6 years follow-up. Abstract number P35 presented at: European Knee Society; May 2–3, 2019; Valencia, Spain.
- Harris A, O’Grady C, Sensiba PR, et al. Guided motion total knee arthroplasty (TKA) in patients with BMI of 40kg/m2 or more: results from the international multicentre study of 2,059 primary TKAs with up to 6 years follow-up. Abstract number P36 presented at: European Knee Society; May 2–3, 2019; Valencia, Spain.
- Ishibashi T, Tomita T, Yamazaki T, Tsuji S, Yoshikawa H, Sugamoto K. Kinematics of bicruciate and posterior stabilized total knee arthroplasty during deep knee flexion and stair climbing. *J Orthop Res.* 2020; doi: 10.1002/jor.24773.

In vivo kinematic comparison of a bicruciate stabilized total knee arthroplasty and the normal knee using fluoroscopy¹¹

Grieco TF, Sharma A, Dessinger GM, Cates HE, Komistek RD. *J Arthroplasty*. 2018;33:565-571

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[Riassunto integrale](#)

+ Evidence in focus




Study summary: Harris AI, et al. *J Arthroplasty* (2019)*

Smith+Nephew

JOURNEY[◇] II BCS total knee arthroplasty (TKA) shows favourable mid-term revision rates compared with registry-reported cemented posterior stabilised (PS) implants

+ Plus points

At 5 years post-TKA, JOURNEY II BCS demonstrated:


Overview

- International, multicentre, retrospective case series analysis of 2,059 JOURNEY II BCS in 1,648 consecutive TKA patients (mean age, 64.3 years; females, 58.5%; patella resurfaced, 95.9%)
- Procedures were performed at seven US and three European sites, with up to 6.1 years of follow up (median, 4.2 years)
- Cumulative incidence of implant revision was compared to data for primary TKAs with cemented PS implants from the Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR)

Results

Compared to the AOANJRR registry control, JOURNEY II BCS demonstrated:

- A comparable, but numerically lower, 5-year revision rate for JOURNEY II BCS (3.6 vs 4.1%; Figure)
- Lower revision rate in patients under 55 (males, 3.0 vs 7.0%; females, 3.1 vs 6.9%)
- Lower percentage of major revisions, defined as all revisions involving the removal of the tibial and/or femoral component (percentage of total revisions, 29.9 vs 41.6%)
 - Major revisions accounted for less than 1% of all JOURNEY II BCS TKAs
 - Total of 67 revisions
 - Most common reason for revision was infection (33%), followed by mechanical loosening (21%), fracture of bone around the joint (16%) and implant instability (15%)

Conclusions

The JOURNEY II BCS knee system performs favourably when compared to PS TKA control from the AOANJRR, particularly in patients <55 years of age. At up to 6.1 years post TKA, less than 1% of all JOURNEY II BCS TKAs required major revision.

Citation

*Harris AI, Christen B, Malcorps JJ, et al. Mid-term performance of a guided motion bicruciate stabilized total knee system; results from the international study of over 2,000 consecutive primary total knee arthroplasties. *J Arthroplasty*. 2019;34:S201–208.
Available at: [The Journal of Arthroplasty](https://doi.org/10.1054/j.1083-7539.2019.34.s201.x)

+ Evidence in focus

Study summary: Grieco TF, et al. *J Arthroplasty* (2018)*

JOURNEY[◇] II BCS exhibits normal-like knee kinematic patterns

+ Plus points

JOURNEY II BCS subjects exhibited similar patterns of femoral rollback and axial rotation to the normal knee in early and late flexion

The dual cam/post mechanism and asymmetric articulate geometries may replicate ACL and PCL function

Overview

- A retrospective, comparative, single-surgeon analysis of the in vivo kinematics of 50 knees through a full weight-bearing range of motion to determine whether the dual cam/post mechanism is able to replicate the cruciate ligament function:
 - 40 implanted with JOURNEY II BCS (average follow up, 14.3 months; mean age, 69.8 years)
 - 10 normal asymptomatic knees (mean age, 57.4 years)
- All TKA subjects were selected as a result of having a well-functioning TKA, with passive flexion of at least 100°, and Knee Society Score (KSS) ≥ 90
- Each patient was asked to perform a weight-bearing deep knee bend (full extension to full flexion)
- Kinematics were measured at full extension and at 30° increments to full flexion

Results

Early flexion (0-30°; ACL function)

- JOURNEY II BCS subjects exhibited similar patterns of femoral rollback and axial rotation compared with normal knee subjects

Mid flexion (30-60°; ACL/PCL translation)


- JOURNEY II BCS subjects experienced minimal anterior-posterior motions and axial rotation, whereas normal knees continued to rollback and externally rotate

Deeper flexion (60-90°; PCL dominant)

- JOURNEY II BCS resumed posterior motion

Late flexion (90°+)

- After 90° axial rotation increased in a normal-like fashion

Conclusions

JOURNEY II BCS exhibited normal-like kinematic patterns and moved as designed under in vivo observation. Similarities in early and late kinematic patterns between the two groups suggest the dual cam-post design and asymmetric articular geometries of the JOURNEY II BCS adequately replicate ACL and PCL function.

Citation

*Grieco TF, Sharma A, Dessinger GM, Cates HE, Komistek RD. In vivo kinematic comparison of a bicruciate stabilized total knee arthroplasty and the normal knee using fluoroscopy. *J Arthroplasty*. 2018;33:565-571.

Available at: [The Journal of Arthroplasty](https://www.jarthroplasty.com)

+ Evidence in focus


Study summary: Nodzo SR, et al. *Tech Orthop* (2018)*

Smith+Nephew

JOURNEY[◇] II BCS is associated with significantly improved flexion and patient-reported outcomes compared with posterior stabilized (PS) total knee arthroplasty (TKA)

+ Plus points

Compared with patients receiving PS TKA, patients with JOURNEY II BCS showed:


Overview

- A retrospective analysis of 200 patients (mean age, 51 years; males, 87; females, 113) who underwent primary TKA by a single surgeon between 2013 and 2014
 - First 100 patients received JOURNEY II BCS
 - Next 100 patients received LEGION[°] PS TKA
- Computer-navigated guides were used to minimize surgical alignment error
- Range of motion (ROM) and Knee Society Scores (KSS) were recorded pre-TKA and at 6 weeks and 1 year post-TKA

Results

Compared with patients receiving LEGION PS TKA, patients with JOURNEY II BCS showed:

- Numerically improved mean KSS score at 6 weeks post-TKA (88 vs 73)
- Significantly improved mean ROM at 1 year post-TKA (119° vs 96°; p<0.0001; Figure)
- Significantly improved mean KSS scores at 1 year post-TKA (89 vs 81; p<0.001)

Conclusions

JOURNEY II BCS led to significant improvements in ROM and patient-reported outcomes at 1 year post-TKA, compared with standard PS TKA. The authors noted that the results suggest that the more anatomic design of the implant, which is intended to replicate a more normal knee position and kinematic patterns, may be responsible for the improved flexion and patient satisfaction, compared with PS TKA.

Citation

*Nodzo, SR; Carroll KM, Mayman DJ. The Bicruciate Substituting Knee Design and Initial Experience. *Tech Orthop*. 2018;33:37-41. Available at: [Techniques in Orthopaedics](#)

+ Evidence in focus

Study summary: Snyder MA, et al. *Orthop Trauma Prosth* (2018)*

Smith+Nephew

JOURNEY[◇] II BCS total knee arthroplasty (TKA) patients experience comparable levels of satisfaction and activity to total hip arthroplasty (THA) patients in short-term follow-up

+ Plus points

Compared to THA patients, JOURNEY II BCS patients reported:

Overview

- Retrospective review of data from a total joint registry in Cincinnati, Ohio, comparing patient outcomes between clinically matched JOURNEY II BCS TKA and THA patients
 - 48 JOURNEY II BCS patients (mean age, 58.3 years; male, 54.2%)
 - 48 THA patients (mean age, 55.9 years; male, 64.6%)
- Overall patient satisfaction, University of California Los Angeles (UCLA) activity scores and EuroQol five-dimension scores (EQ-5D) were compared at 3 months and 1 year post-op

Results

Compared to THA patients, JOURNEY II BCS patients reported

- No significant difference in overall satisfaction at 3 months ($p = 0.398$) and 1 year ($p = 0.590$; Figure)
- Significantly improved median UCLA activity scores at 3 months (8 vs 7; $p = 0.028$) and 1 year (8 vs 7; $p < 0.001$)
- Significantly improved median EQ-5D scores at 3 months (90 vs 80; $p < 0.001$), no significant difference in EQ-5D scores at 1 year
- No significant difference in patient quality of life measures (time to return to work, activities of daily living or sports activities) at 3 months and 1 year
- Similar clinical improvements over time in pain and joint function
 - 43.8% patients scored ≥ 95 on the Knee Society (KS) pain score at 3 months, which increased to 91.7% at 1 year
 - 37.5% patients scored ≥ 90 on the KS function score at 3 months which increased to 95.8% at 1 year


Figure. Overall satisfaction for JOURNEY II BCS and THA patients at 1 year post-op
Boxes represent range of values, median value represented by dark grey line.

Conclusions

JOURNEY II BCS TKA patients reported similar levels of satisfaction and patient reported outcomes as THA patients at 3 months and 1 year post-op. The kinematic design advancements of JOURNEY II BCS demonstrate improvements in patient satisfaction following TKA compared to past TKA procedures, via comparison with THAs.

Citation

*Snyder MA, Sympson A, Gregg J, Levit A. A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY II bi-cruciate stabilizing knee system and total hip arthroplasty patients. *Orthop Trauma Prosth*. 2018; <http://doi.org/10.15674.0030-5987201835>.

Available at: [Orthopaedics, Traumatology and Prosthetics](#)

+ Evidence in focus

Publication summary: Smith LA, et al. *J Arthroplasty* (2020)*

Femoral geometry and the cruciate ligaments may play a significant role in achieving more normal-like kinematic patterns post-total knee arthroplasty (TKA)

+ Plus points

JOURNEY[®] II XR[®] knees experienced more **normal-like kinematic patterns** compared to JOURNEY II CR, demonstrating the importance of the ACL

The authors stated JOURNEY II CR knees demonstrated an **improvement in lateral femoral rollback and axial rotation** compared to previous studies on CR knees

Overview

- First study to assess the impact of the ACL, PCL and femoral condylar geometry on kinematic patterns following TKA and comparing this to the normal knee
 - 40 JOURNEY II CR knees (average follow up, 24.8 months; mean age, 68.8 years)
 - 10 JOURNEY II XR knees (average follow up, 16 months; mean age, 62.3 years)
 - 10 normal knees (mean age, 57.4 years)
- All TKA subjects were selected as a result of having a well-functioning TKA, with a Knee Society Score (KSS) ≥ 80
- Each patient was asked to perform a weight-bearing deep knee bend (full extension to full flexion; Figure)
- Kinematics were measured at full extension and at 30° increments to full flexion

Results

Compared to the normal knee:

Early flexion (0-30°; ACL function)

- JOURNEY II XR demonstrated a similar magnitude of posterior femoral rollback (PFR)
- JOURNEY II CR demonstrated PFR, but significantly less than the normal knee (lateral, $p=0.004$; medial, $p=0.002$)
- Both JOURNEY II XR and CR exhibited external femoral rotation, but to a lower magnitude

Mid flexion (30-60°; ACL/PCL translation)


- JOURNEY II XR demonstrated no statistical difference in anterior-posterior translation or axial rotation
- JOURNEY II CR demonstrated significantly less anterior-posterior translation (lateral, $p<0.004$; medial, $p<0.03$) and axial rotation $p=0.01$
 - The authors stated that JOURNEY II CR did not experience the significant amount of paradoxical sliding seen in other PCR TKAs

Deeper flexion (60-90°; PCL dominant)

- No significant difference with JOURNEY II XR or CR in anterior-posterior translation
- JOURNEY II XR and CR experienced similar axial rotation patterns where the femur externally rotated relative to the tibia

Late flexion (90°+)

- JOURNEY II XR and CR demonstrated no significant differences in anterior-posterior translation or axial rotation

Conclusions

The retention of the ACL did benefit the JOURNEY II XR subjects compared to the CR subjects. However, the authors stated that JOURNEY II CR subjects did achieve increased amounts of lateral condyle PFR and axial rotation compared to previously studied CR TKAs. Anterior-cruciate ligament retention and subtle changes in femoral geometry may play an important role in kinematic outcomes.

Citation

*Smith LA, Nachtrab J, LaCour M, Cates H, Freeman MG, Komistek RD. In Vivo Knee Kinematics: How Important Are the Roles of Femoral Geometry and the Cruciate Ligaments? *J Arthroplasty*. 2020: doi: <https://doi.org/10.1016/j.arth.2020.10.020>.

Available at: [Journal of Arthroplasty](https://doi.org/10.1016/j.arth.2020.10.020)

+ Evidence in focus

Publication summary: Iriuchishima T and Ryu K. *Indian J Orthop* (2019)*

JOURNEY[◇] II BCS patients demonstrate significantly better stair climbing ability when compared with cruciate retaining (CR) and posterior stabilised (PS) total knee arthroplasty (TKA) patients

+ Plus points

89%

JOURNEY II BCS patients were able to climb stairs "one step at a time"

JOURNEY II BCS patients had **significantly better stair climbing ability** compared to CR and PS TKA patients

Overview


Case-controlled, single surgeon study comparing stair climbing and descending ability in JOURNEY II BCS patients compared to CR and PS TKA patients.

Patients were administered a questionnaire pre-TKA and 12 months post-TKA to assess their stair climbing and descending ability and necessity of a handrail.

- 63 JOURNEY II BCS patients (TKAs, 66; mean age, 75.2 years)
- 47 CR patients (FINE[®] Total Knee System [Teijin Nakashima Medical, Okayama, Japan]; number of TKAs, 56; mean age, 76.4 years)
- 38 PS patients (FNK[®] Total Knee System [Teijin Nakashima Medical, Okayama, Japan]; number of TKAs, 42; mean age, 78.6 years)

Results

- At 12 months post-TKA, a significantly higher percentage of JOURNEY II BCS patients (89%) were able to climb stairs "one step at a time" compared to CR (72%) and PS TKA patients (58%; $p < 0.05$; Figure)
- No significant differences in the percentage of patients able to descend stairs "one step at a time" or among handrail necessity between JOURNEY II BCS (59 and 75%), CR (49 and 77%) and PS (45 and 89%) TKA patients

JOURNEY II BCS: 89%
CR TKA: 72%
PS TKA: 58%

Figure: Percentage of patients able to climb stairs "one step at a time" 1-year post-TKA

Conclusions

JOURNEY II BCS results in significantly better "one step at a time" stair climbing ability compared to CR and PS TKAs.

Citation

*Iriuchishima T, Ryu K. Bicruciate substituting total knee arthroplasty improves stair climbing ability when compared with cruciate-retain or posterior stabilizing total knee arthroplasty. *Indian J Orthop*. 2019;53:641-645. Available at: [Indian Journal of Orthopaedics](https://doi.org/10.4103/0013-1226.268888)

+ Evidence in focus

Publication summary: Kono K, et al. *J Knee Surg* (2019)*

Smith+Nephew

JOURNEY[◇] II XR[◇] total knee arthroplasty (TKA) patients exhibit normal-like kinematics during high-flexion activities

+ Plus points

JOURNEY II XR kinematics were similar to the normal knee

Kinematics with JOURNEY II XR differed depending on high-flexion activity

Overview

- Evaluation of in vivo kinematics of JOURNEY II XR TKA patients during high-flexion activities
- 15 patients (17 TKAs; mean age, 72 years; mean follow-up, 7.6 months)
- Patients were examined under fluoroscopy during squatting and cross-legged sitting motions, starting with legs fully extended

Results

- Rotation angle
 - During squatting: steep femoral external rotation from extension to early-flexion (Figure), which may be indicative of a screw-home mechanism[‡]
 - During cross-legged sitting: slight femoral internal rotation from early-flexion to mid-flexion, followed by femoral external rotation beyond mid-flexion (Figure)
 - Significantly larger external rotation from 80 to 110° of flexion during squatting compared to cross-legged sitting ($p < 0.05$; Figure)
- Varus-valgus angle
 - During squatting: no significant movement
 - During cross-legged sitting: varus movement (6.1°) beyond 110° of flexion
- Anteroposterior translation
 - No significant differences in anterior movement during squatting or cross-legged sitting
 - Paradoxical motion was not observed[‡]
 - During squatting: medial and lateral contact points moved posteriorly from extension to mid-flexion
 - During cross-legged sitting: some posterior movement of the medial contact point, but no significant movement of the lateral contact point until movements beyond mid-flexion
- Kinematic pathway
 - During squatting: medial pivot pattern in early flexion up to 10°
 - During cross-legged sitting: medial pivot pattern from 60° of flexion


Figure. Rotation angle during squatting and cross-legged sitting. * $p < 0.05$

Conclusions

This is the first study to evaluate the in vivo kinematics of JOURNEY II XR patients during high-flexion activities of daily living. JOURNEY II XR patients exhibited kinematics which were consistent with those reported in the literature for normal knees and which differed depending on activity.

Citation

*Kono K, Inui H, Tomita T, Yamazaki T, Taketomi S, Tanaka S. In vivo kinematics of bicruciate-retaining total knee arthroplasty with anatomical articular surface under high-flexion conditions. *J Knee Surg*. 2019. [Epub ahead of print] Available at: [The Journal of Knee Surgery](https://doi.org/10.1177/1098314619851111)

[‡] Normal knee kinematics were not assessed in this study, comparison are drawn from existing literature

+ Evidence in focus



Poster summary: Mayman DJ, et al. ISPOR Symposium (2018)*

Smith+Nephew

JOURNEY[◇] II BCS is associated with significantly reduced total hospital cost and hospital stays, and patients are less likely to be readmitted within 30 days compared with other total knee arthroplasty (TKA) systems

+ Plus points

Compared to all other TKA systems, JOURNEY II BCS resulted in:

Overview

- A retrospective cohort study from the Premier Perspective Database (2014 to 2016) compared patients with osteoarthritis undergoing primary TKA with JOURNEY II BCS or other TKA systems
 - 1,692 JOURNEY II BCS patients (mean age, 64.1 years)
 - 1,692 other TKA patients (mean age, 63.9 years)
- Hospital-related clinical and economic outcomes (reported as 2016 US dollars) were compared between the groups
- 1:1 propensity score matching was used to control patients and provider characteristics

Results

Compared with patients receiving other TKA systems, those receiving JOURNEY II BCS:

- Experienced significantly reduced mean patient hospital costs (p<0.0001; Figure)
- Were 51% less likely to be readmitted to hospital within 30 days (p=0.0037; calculated value)
- Were 35% more likely to be discharged to home (p=0.0008; calculated value)
- Were 41% less likely to be discharged to a skilled nursing facility (p<0.0001)
- Experienced significantly reduced mean length of hospital stay (2.45 days vs 2.66 days; p<0.0001)


Figure. Mean total hospital costs (2016 US dollars; adjusted to nearest dollar)

Conclusions

JOURNEY II BCS led to significant reductions in total hospital cost, significant reductions in hospital stay, and patients were less likely to be readmitted within 30 days. Patients receiving this device were also more likely to be discharged to home when compared with patients undergoing primary TKA with other posterior-stabilized TKA systems.

Citation

*Mayman DJ, Patel AR, Carroll KM. Hospital related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patients. Poster presented at: ISPOR Symposium; May 19-23, 2018; Baltimore, Maryland, USA.

+ Evidence in focus

Publication summary: Di Benedetto P, et al. *Acta Biomed* (2019)*

Smith+Nephew

JOURNEY[◇] II Cruciate Retaining (CR) Knee System is associated with better early postoperative results in kinematics and muscle strength during walking compared to Attune[®] CR total knee arthroplasty (TKA)

+ Plus points

Compared to Attune CR, JOURNEY II CR patients demonstrated:

Significantly improved knee flexion during walking at 3 months post-TKA ($p < 0.01$)

Significantly improved external rotation during walking at 3 months post-TKA ($p < 0.01$)


Physiological activation timing of investigated muscles

Overview

- An independent pilot study evaluating pre- and 3 months post-TKA kinematics during walking in patients randomly assigned to receive JOURNEY II CR or Attune CR (Depuy, Warsaw, IN)
 - 12 JOURNEY II CR patients (mean age, 70.25 years; number of females, 8)
 - 14 Attune CR patients (mean age, 71.75 years; number of females, 8)
 - Two patients did not make the 3 month follow-up
- The kinematic parameters assessed were flexion-extension, internal-external rotation and abduction-adduction for both the operated and contralateral side, pre- and post-TKA
- Muscle strength and activation was evaluated using electromyography on both operated and contralateral side, pre- and post-TKA
- Knee Society Score (KSS) and Knee Injury and Osteoarthritis Outcomes Score (KOOS) were also assessed pre- and 3 months post-TKA

Results

- At 3 months post-TKA, compared with the Attune CR operated knee group, the JOURNEY II CR operated knee group demonstrated:
 - Significantly greater increase in flexion whilst walking compared to pre-TKA (3.76 vs -3.58° ; $p < 0.01$; Figure)
 - Significantly greater increase in external rotation whilst walking compared to pre-TKA (-1.35 vs 0.86° ; $p < 0.01$)
 - No significant difference in extension, internal rotation, abduction or adduction
 - Significantly lower reduction of muscle strength ($p < 0.01$) in most muscle types (4 out of 6 vs 1 out of 6 muscle groups for Attune CR)
 - Significantly lower reduction of muscle strength ($p < 0.01$) in 1 out of 6 muscle groups in the non-operated side (compared to 2 out of 6 for Attune CR)
 - Better muscle activation timing, comparable to the physiological knee
 - Attune CR knee showed incorrect activation timing of all muscles investigated
 - A numerically higher KSS satisfaction score (37.5 vs 35.4 ; maximum score, 40)
 - A numerically higher KOOS quality of life score (84.9 vs 76.2 ; maximum score, 100)

Increase in knee flexion:

+3.76°
JOURNEY II CR

-3.58°
Attune CR

Figure. Change in knee flexion angle during walking at 3 months post-TKA, compared to pre-TKA for JOURNEY II CR and Attune CR

Conclusions

Compared to Attune CR, JOURNEY II CR patients experienced significantly better flexion, external rotation and muscle strength during walking at 3 months post-TKA. In addition, JOURNEY II CR patients demonstrated physiological activation timing of most muscles.

Citation

*Di Benedetto P, Vidi D, Colombo, Buttironi MM, Cainero V, Causero A. Pre-operative and post-operative kinematic analysis in total knee arthroplasty. A pilot study. *Acta Biomed*. 2019;90:91-97.

Available at: [Acta Biomedica](https://doi.org/10.2165/0000131221990091097)

+ Evidence in focus


Publication summary: Hyodo K, et al. *Arthroplast Today* (2020)*

Smith+Nephew

JOURNEY[◇] II BCS patients demonstrate a more normal-like gait pattern compared to conventional posterior stabilised (PS) total knee arthroplasty (TKA)

+ Plus points

Compared to conventional PS TKA, JOURNEY II BCS demonstrated:


Overview

- Gait analysis study comparing JOURNEY II BCS and a conventional PS TKA (LEGION[◇] PS TKA)
 - JOURNEY II BCS group: n=12 patients (mean age, 69.4 years; percentage of males, 16.7%; mean BMI, 23.5kg/m²)
 - LEGION TKA group: n=12 patients (mean age, 70.0 years; percentage of males, 33.3%; mean BMI, 26.5kg/m²)
- Motion capture analysis with a force plate and electromyogram of the lower limb muscles were conducted at 6 months post-TKA

Results

Compared to LEGION PS TKA, JOURNEY II BCS resulted in:

- Significantly faster gait speed (1.1 vs 0.9m/s; p=0.03; Figure)
- Significantly longer step length (0.54 vs 0.45m; p=0.03; Figure)
- Significantly greater maximum knee flexion angle during initial stance phase (10.9 vs 8.0°; p=0.04)
- Significantly greater maximum knee extension moment (0.55 vs 0.34Nm/kg; p=0.04)
- No significant difference in quadriceps, hamstring or gluteus medius activity
 - Less quadriceps muscle activity in both stance and swing phases (p=ns)
 - Less hamstring muscle activity in the swing phase (p=ns)

Conclusions

JOURNEY II BCS patients demonstrate a faster gait speed, longer step length and greater knee extension moment, compared to conventional PS TKA. This suggests that the quadriceps muscles are more effectively activated and that anterior stability function of the anterior cruciate ligament can be achieved with JOURNEY II BCS TKA compared to conventional PS TKA.

Citation

*Hyodo K, Kanamori A, Kadone H, Takahashi T, Kajiwara M, Yamazaki M. Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty. *Arthroplast Today*. 2020;6:338-342.

Available at: [Arthroplast Today](https://doi.org/10.1016/j.artpl.2020.08.001)



+ Evidence in focus

Study summary: Iriuchishima T, Ryu K. *J Knee Surg* (2017)*

Smith+Nephew

JOURNEY[◇] II BCS total knee arthroplasty (TKA) results in comparable rollback ratio and knee flexion angle to unicompartmental knee arthroplasty (UKA) and the normal knee

+ Plus points


Overview

- Single-surgeon retrospective study comparing the rollback ratio in JOURNEY II BCS and Oxford™ UKA (Zimmer Biomet, Warsaw, Indiana, USA)
 - 64 JOURNEY II BCS patients (mean age, 71.3 years)
 - 50 Oxford UKA patients (mean age, 73.8 years)
 - 10 contralateral asymptomatic knees of subjects with UKA
- When subjects had recovered their range of knee flexion (6-9 months post-op), lateral radiographs in active flexion were taken to measure rollback ratio and flexion angle

Results

- Compared to Oxford UKA and the normal knee, JOURNEY II BCS demonstrated no significant difference in rollback ratio (Figure) or knee flexion angle
- There was a significant correlation between rollback ratio and knee flexion angle among the three groups ($p=0.002$)

Conclusions

JOURNEY II BCS showed no significant difference in rollback ratio when compared with UKA or asymptomatic control knees. The results suggest that JOURNEY II BCS reproduces native anterior and posterior cruciate function and native knee rollback.

Citation

*Iriuchishima T, Ryu K. A comparison of rollback ratio between bicruciate substituting total knee arthroplasty and oxford unicompartmental knee arthroplasty. *J Knee Surg*. 2018;31:568-572.

Available at: [The Journal of Knee Surgery](https://doi.org/10.1177/1098512217711111)

+ Evidence in focus

Study summary: Lutes W and Fitch D, SICOT Orthopaedic World Congress (2018)*

Smith+Nephew


JOURNEY[◇] II Cruciate Retaining (CR) Knee System demonstrates significantly improved functional outcomes compared to P.F.C. Sigma[™] CR total knee arthroplasty (TKA)

+ Plus points

Compared to P.F.C. Sigma, JOURNEY II CR demonstrated:

Significantly improved KSS at 3, 6, 12 and 24 months post-TKA ($p < 0.017$)

Significantly improved WOMAC scores at 6 and 12 months post-TKA ($p < 0.019$)

Significantly greater change in ROM from baseline at 3, 6 and 12 months post-TKA ($p < 0.05$)

Overview

- Retrospective comparison of functional outcomes following JOURNEY II CR or conventional CR TKA carried out by a single surgeon between September 2012 and June 2014
 - 52 JOURNEY II CR patients (mean age, 67.3 years)
 - 60 P.F.C. Sigma[™] [DePuy Synthes Orthopaedics, Warsaw, IN, USA] CR patients (mean age, 70.2 years)
- Knee society scores (KSS) and Western Ontario and McMaster Universities Arthritis Index (WOMAC) scores were assessed and compared pre-TKA and at 3, 6, 12 and 24 months post-TKA
- Range of motion (ROM) was assessed and compared pre-TKA and at 3, 6 and 12 months post-TKA

Results

Compared to P.F.C. Sigma, JOURNEY II CR patients reported:

- Significantly improved KSS at all time points post-TKA (Figure):
 - 3 months; 69.5 vs 63.0 ($p = 0.016$)
 - 6 months; 84.4 vs 70.1 ($p = 0.043$)
 - 12 months; 93.0 vs 86.1 ($p < 0.001$)
 - 24 months; 96.4 vs 91.7 ($p = 0.006$)
- Significantly improved WOMAC scores at 6 (17.8 vs 24.6; $p = 0.018$) and 12 (12.4 vs 18.5; $p = 0.008$) months post-TKA, differences at other time points were not significant
- Significantly greater change in ROM from baseline at 3 (-4.4 vs -10.1; $p < 0.05$), 6 (5.8 vs -1.8; $p < 0.05$) and 12 (11.4 vs 4; $p < 0.05$) months post-TKA


Figure. Mean KSS scores of JOURNEY II CR and P.F.C. Sigma from pre-TKA to 24 months post-TKA

Conclusions

JOURNEY II CR demonstrated significant improvements in short-term functional outcomes compared to conventional CR TKA. The results of this study suggest that choosing a knee implant designed to more closely replicate normal knee kinematics may improve patient-reported outcomes, and in turn patient satisfaction, compared to conventional TKA.

Citation

*Lutes W and Fitch D. Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics. Presented at: 39th SICOT Orthopaedic World Congress; October 10-13, 2018; Montreal, Canada.

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©2021 Smith+Nephew