MICRORAPTOR[®] Knotless Suture Anchor

10601402-A

EC REP

European Representative

Smith & Nephew Medical Limited
101 Hessle Road
Hull HU3 2BN
United Kingdom
ec.rep@smith-nephew.com

Endoscopy

Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810
USA

09/2018 10601402 Rev. A

MICRORAPTOR Knotless Suture Anchor

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MICRORAPTOR[®] Knotless Suture Anchor

STERILE EO

CE mark on the package insert (IFU) is not valid unless there is a CE mark on the product (description) label.

Device Description

The Smith & Nephew MICRORAPTOR[®] Knotless Suture Anchor is a fixation device intended to provide secure reattachment of soft tissue to bone. The device consists of a knotless suture anchor assembled on to an insertion device.

This device is provided sterile, for single use only.

Material

	MICRORAPTOR Knotless REGENESORB	MICRORAPTOR Knotless PEEK
Anchor body	PLGA (poly(lactic-co-glycolic acid)), β -TCP (beta tricalcium phosphate), and calcium sulfate.	PEEK (Polyether ether ketone)
Anchor tip	PEEK (Polyether ether ketone)	PEEK (Polyether ether ketone)

Contents

The Smith & Nephew MICRORAPTOR Knotless Suture Anchor contains:

- 1 ea. Absorbable biocomposite suture anchor with non-absorbable PEEK tip
or
- 1 ea. Non-absorbable PEEK suture anchor with non-absorbable PEEK tip
with
- 1 ea. Insertion device – stainless steel shaft with ABS/polycarbonate handle
- 1 ea. Suture threader – nitinol wire crimped into a stainless steel tube with ABS/polycarbonate handle

Intended Use

The Smith & Nephew MICRORAPTOR Knotless Suture Anchor is intended for the reattachment of soft tissue to bone.

Indications for Use

Smith & Nephew MICRORAPTOR Knotless Suture Anchor is intended for use only for the reattachment of soft tissue to bone for the following indications:

Hip

Acetabular labrum repair/reconstruction

Shoulder

- Capsular stabilization
 - Bankart repair
 - Anterior shoulder instability
 - SLAP lesion repairs
 - Capsular shift or capsulolabral reconstructions
- Biceps tenodesis
- Rotator cuff tear repairs

Contraindications

- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Pathological conditions of bone, such as cystic changes or severe osteopenia, which would compromise secure anchor fixation.
- Pathological conditions in the soft tissues to be attached that would impair secure fixation by suture.
- Comminuted bone surface, which would compromise secure anchor fixation.
- Physical conditions which would eliminate, or tend to eliminate, adequate anchor support or retard healing.
- The presence of infection.
- Conditions which would limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- Contraindications may be relative or absolute and must be carefully weighed against the patient's entire evaluation.

Warnings

 Do not use if package is damaged. Do not use if the product sterilization barrier or its packaging is compromised.

- Contents are sterile unless package is opened or damaged. **DO NOT RESTERILIZE.** For single use only. Discard any open, unused product. Do not use after the expiration date.
- It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.
- Read these instructions completely prior to use.
- Product must be stored in the original sealed pouch.
- Incomplete anchor insertion may result in the anchor not functioning as intended.
- Breakage of the suture anchor can occur if insertion sites are not prepared with the appropriate instrumentation prior to implantation.
- Only use the recommended drills, MICRORAPTOR drill guides, and MICRORAPTOR obturators intended for use with the MICRORAPTOR Knotless Suture Anchor. Use of other instruments may injure the patient, damage the instruments, or compromise fixation.
- Maintaining drill guide alignment throughout drilling is required to ensure drill integrity.
- Do not attempt to implant this device within cartilage epiphyseal growth plates or non-osseous tissue.
- Do not resterilize or reuse anchors, sutures, insertion devices, or disposable hole preparation instrumentation.

Precautions

 U.S. Federal law restricts this device to sale by or on the order of a physician.

- Hazards associated with reuse of this device include, but are not limited to, patient infection and/or device malfunction.
- Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.
- Postoperative care is important. A patient should be instructed on the limitations of the implant and should be cautioned regarding weight bearing and body stresses on the appliance prior to secure bone healing. Postoperative range of motion is to be determined by the physician.
- Do not use sharp instruments to manage or control the suture.
- As in all suture anchor or suturing techniques, until biological attachment of tissue to bone is complete the fixation should be considered temporary and may not withstand weight bearing or other unsupported stresses. The suture anchor and suture are not intended to provide indefinite biomechanical integrity.
- Implantation of the MICRORAPTOR Knotless Suture Anchor requires preparation of the insertion site. Pre-drilling with the appropriate Smith & Nephew drill is the preferred method of site preparation.
- Ensure the anchor placement is aligned with the drilled hole. Proper alignment is essential for successful repair.
- Use of excessive force during insertion can cause failure of the suture anchor or insertion device.
- Bone quality must be adequate to allow proper placement of suture anchor. Inadequate bone quality could result in loss of fixation or pullout of the suture anchor.
- Do not alter implant or instrumentation or performance may be compromised.
- After use, this device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.
- For dense bone, it may be required to cycle the drill in and out of the implantation site to ensure complete removal of bone.

Adverse Reactions

- Loss of fixation or pullout of suture anchors can occur
- Mild inflammatory reaction
- Foreign body reaction
- Infection, both deep and superficial
- Allergic reaction

Instructions for Use

⚠ WARNING: Use of the appropriate Smith & Nephew instrumentation is necessary to prepare the insertion site for the MICRORAPTOR[°] Knotless Suture Anchor.

When using a MICRORAPTOR Knotless Suture Anchor, use Smith & Nephew drill guide(s), obturator(s), and drill specific to the MICRORAPTOR Knotless Suture Anchor to prepare the insertion site (each sold separately).

1. Pass suture through soft tissue in desired stitch configuration.
2. Thread the free ends of suture through the suture threader, remove the suture threader from the shaft of the insertion device, and pull to feed the sutures through the anchor eyelet.
3. Insert the appropriate drill guide and obturator through the operative cannula or percutaneously. Remove the obturator.
4. Place the distal tip of the guide onto the bone at the desired implantation site.
5. While firmly holding the guide in place, use the following recommended drills to prepare the insertion site:

Shoulder General Use	Shoulder Hard Bone
REF 72205022, 2.2 mm Drill	REF 72205169, 2.6 mm Drill

Hip General Use	Hip Soft Bone
REF 72205169, 2.6 mm Drill	REF 72205022, 2.2 mm Drill

- a. The depth stop on the drill will bottom out on the proximal end of the guide when proper hole depth is reached.
CAUTION: For dense bone, it may be required to cycle the drill in and out of the implantation site to ensure complete removal of bone.
6. While holding the drill guide steady, remove the drill from the insertion site by backing the drill axially out of the hole.
CAUTION: Ensure drill is always powered during removal. Do not stop the drill between steps 5 and 6.
 7. Remove guide.
 8. Advance the loaded suture anchor into the operative cannula. Remove excess suture slack and approximate the anchor to the prepared bone site.
 9. Align the eyelet to the repair site and ensure that the sutures are not twisted around the anchor.
 10. Establish and maintain axial alignment of the suture anchor with the prepared insertion site. Place the tip of the anchor into the prepared hole.
 11. Apply tension to suture limbs to reduce labrum to desired location. Use a mallet to tap the inserter handle until the laser mark is flush with the cortical bone. This places the suture anchor approximately 2 mm below the bone surface.
 12. While maintaining slight downward pressure on the inserter, locate the knob on the proximal end of the inserter handle and rotate it clockwise. Continue rotating the knob until a minimum of five (5) clicks are heard.
 13. Slowly disengage the suture anchor from the inserter by pulling straight back or back tapping the distal end of the handle with a mallet. Do not wiggle the handle to disengage the anchor.
 14. Discard the insertion device and cut the excess suture.

Warranty

For single use only. This product is warranted to be free from defects in material and workmanship. Do not reuse.

For Further Information

If further information on this product is needed, contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized representative.

[°]Trademark of Smith & Nephew, registered U.S. Patent & Trademark Office.

©2018 Smith & Nephew, Inc. All rights reserved.

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[®]Varumärke som tillhör Smith & Nephew, och har registrerats vid U.S. Patent and Trademark Office (USA:s Patent- och registreringsverk).
°Varemærke, der tilhører Smith & Nephew og er registreret i det amerikanske patent- og varemærkedirektorat (U.S. Patent and Trademark Office).