

Smith+Nephew

JOURNEY[◇] II
Total Knee Arthroplasty

System
Specification
Guide



Table of contents

Please use this guide to learn more about the JOURNEY[◇] II TKA System technical specifications.

Implant Dimensions and Insert Compatibility

Femorals	1
Tibial Inserts	4
Tibial Baseplates	11
Patella	13
Tibial Stems	15

Catalog Item Numbers

Femoral Implants	18
Tibial Inserts	19
Revision Insert Options	24
Tibial Baseplates	27
Patella	28
Tibial Wedges	31
Tibial Stems	32


SYNC Instrumentation

Alignment Instruments	35
General Instruments	41
Femoral Finishing Instruments	42
Patella Instruments	43
Femoral Trials and Outliers	46
Tibial Baseplate Trials	48
Insert Trials	49

Femoral component dimensions (mm)

JOURNEY[◇] II BCS



Size	Anterior Posterior	Medial Lateral	PS Box Width	Posterior Condylar Offset	Flange Height	Distal Medial Thickness	Distal Lateral Thickness	Posterior Medial Thickness	Posterior Lateral Thickness	Box Height
Size	A	B	C	D	E	F	G	H	J	K
1	51.7	59.0	16.5	1.7	49.5	9.5	7	9	7.4	16.0
2	53.7	60.0	16.5	1.7	50.7	9.5	7	9	7.4	17.0
3	56.7	61.5	16.5	1.7	52.5	9.5	7	9	7.4	17.0
4	59.7	64.5	16.5	1.7	54.3	9.5	7	9	7.4	20.5
5	62.7	67.5	16.5	1.7	56.0	9.5	7	9	7.4	20.5
6	65.7	70.5	16.5	1.8	57.7	9.5	7	9	7.4	22.0
7	68.8	73.5	16.5	1.8	59.5	9.5	7	9	7.4	22.0
8	71.8	76.0	16.5	1.8	61.2	9.5	7	9	7.4	22.0
9	75.8	80.0	16.5	1.8	63.5	11.5	9	11	9.4	22.8
10	79.8	82.0	16.5	1.8	65.7	11.5	9	11	9.4	22.8



Femoral component dimensions (mm)

JOURNEY[◇] II CR

Size	A	B	C	D	E	F	G	H	J
1	51.7	59.0	19	1.7	49.5	9.5	7	9	7.4
2	53.7	60.0	19	1.7	50.7	9.5	7	9	7.4
3	56.7	61.5	19	1.7	52.5	9.5	7	9	7.4
4	59.7	64.5	19	1.7	54.3	9.5	7	9	7.4
5	62.7	67.5	19	1.7	56.0	9.5	7	9	7.4
6	65.7	70.5	19	1.8	57.7	9.5	7	9	7.4
7	68.8	73.5	19	1.8	59.5	9.5	7	9	7.4
8	71.8	76.0	19	1.8	61.2	9.5	7	9	7.4
9	75.8	80.0	19	1.8	63.5	11.5	9	11	9.4
10	79.8	82.0	19	1.8	65.7	11.5	9	11	9.4




Femoral component dimensions (mm)

JOURNEY[◇] II vs LEGION[◇]


Femoral sizing guide A/P

JOURNEY II		LEGION	
Size	(mm)	Size	(mm)
1	51.7	2	50.0
2	53.7	3N	54.0
3	56.7	3	55.0
4	59.7	4N	58.0
5	62.7	4	59.0
6	65.7	5N	61.0
7	68.8	5	62.0
8	71.8	6N	65.0
9	75.8	6	66.0
10	79.8	7	70.0
		8	75.0



Femoral sizing guide M/L

JOURNEY II		LEGION	
Size	(mm)	Size	(mm)
1	59.0	2	58.0
2	60.0	3N	58.0
3	61.5	3	62.0
4	64.5	4N	62.0
5	67.5	4	66.0
6	70.5	5N	66.0
7	73.5	5	70.0
8	76.0	6N	70.0
9	80.0	6	73.0
10	82.0	7	77.0
		8	80.0



Tibial insert dimensions (mm)

JOURNEY[◇] II BCS

9mm Insert	Anterior Posterior	Medial Lateral	Medial Thickness*	Lateral Thickness*	Post Height*
	A	B	C	D	E
Sz 1-2	42	60	9.6	11.9	34.1
Sz 3-4	48	68	9.6	11.9	35.1
Sz 5-6	52	74	9.6	11.9	38.6
Sz 7-8	56	81	9.6	11.9	40.1



Minimum polyethylene thickness for a 9mm metal-backed component is 6.7mm on the medial side.

*Baseplate thickness included.

Insert offering/compatibility

Insert Size	Femoral Size									
	1	2	3	4	5	6	7	8	9	10
1-2	●	●	●	●						
3-4		●	●	●	●	●				
5-6				●	●	●	●	●	●	
7-8*						●	●	●	●	●

*Size 7-8 tibia inserts can also be used with size 9 tibia baseplates.



Tibial insert dimensions (mm)

JOURNEY[◇] II Articular Stabilized (AS)[†]

9mm Insert	Anterior Posterior A	Medial Lateral B	Medial Thickness* C	Lateral Thickness* D	Post Height E
Sz 1-2	42	60	9.6	12.1	34.1
Sz 3-4	48	68	9.6	12.1	35.1
Sz 5-6	52	74	9.6	12.1	38.6
Sz 7-8	56	81	9.6	12.1	40.1

Minimum polyethylene thickness for a 9mm metal-backed component is 6.7mm on the medial side.


*Baseplate thickness included.

Insert offering/compatibility

Insert Size	Femoral Size									
	1	2	3	4	5	6	7	8	9	10
1-2	●	●	●	●						
3-4		●	●	●	●	●				
5-6				●	●	●	●	●	●	
7-8*						●	●	●	●	●

*Size 7-8 tibia inserts can also be used with size 9 tibia baseplates.

†Not widely available



Tibial insert dimensions (mm)

JOURNEY[◇] II Constrained

	Anterior Posterior	Medial Lateral	Medial Thickness*	Lateral Thickness*	Post Height**	Post Width
9mm Insert	A	B	C	D	E	F
Sz 1-2	42	60	9.6	12.1	34.1	16.1
Sz 3-4	48	68	9.6	12.1	35.3	16.1
Sz 5-6	52	74	9.6	12.1	38.6	16.1
Sz 7-8	56	81	9.6	12.1	40.1	16.1


Minimum polyethylene thickness for a 9mm metal-backed component is 6.7mm on the medial side.

*Baseplate thickness included.

Insert offering/compatibility

Insert Size	Femoral Size									
	1	2	3	4	5	6	7	8	9	10
1-2	●	●	●	●						
3-4		●	●	●	●	●				
5-6				●	●	●	●	●	●	
7-8*						●	●	●	●	●

*Size 7-8 tibia inserts can also be used with size 9 tibia baseplates.



JOURNEY II Constrained Insert allows:

+/- 3° of Varus/Valgus freedom

+/- 10° of rotational freedom



Tibial insert dimensions (mm)

JOURNEY[◇] II CR

	Anterior Posterior	Medial Lateral	Medial Thickness*	Lateral Thickness*
9mm Insert	A	B	C	D
Sz 1-2	42	60	9.6	12.0
Sz 3-4	48	68	9.6	12.0
Sz 5-6	52	74	9.6	11.9
Sz 7-8	56	81	9.6	11.9



Minimum polyethylene thickness for a 9mm metal-backed component is 6.7mm on the medial side.

*Baseplate thickness included.

Insert offering/compatibility

JOURNEY II CR is completely interchangeable with all size femoral components.

*Size 7-8 tibia inserts can also be used with size 9 tibia baseplates.




Tibial insert dimensions (mm)

JOURNEY[◇] II Medial Dished (MD)

9mm Insert	Anterior Posterior	Medial Lateral	Medial Thickness*	Lateral Thickness*	Anterior Thickness
	A	B	C	D	E
Sz 1-2	42	60	9.6	12.1	16.2
Sz 3-4	48	68	9.6	12.1	16.9
Sz 5-6	52	74	9.6	12.1	18.8
Sz 7-8	56	81	9.6	12.1	19.5

Minimum polyethylene thickness for a 9mm metal-backed component is 6.7mm on the medial side.


*Baseplate thickness included.



Insert offering/compatibility

Insert Size	Femoral Size									
	1	2	3	4	5	6	7	8	9	10
1-2	●	●	●	●						
3-4		●	●	●	●	●				
5-6				●	●	●	●	●	●	
7-8*						●	●	●	●	●

*Size 7-8 tibia inserts can also be used with size 9 tibia baseplates.



Tibial insert dimensions (mm)


JOURNEY[◇] II Deep-Dished

	Anterior Posterior	Medial Lateral	Medial Thickness*	Lateral Thickness*	Anterior Thickness*
9mm Deep-Dished Insert	A	B	C	D	E
Sz 1-2	42	60	9.6	12.1	16.9
Sz 3-4	48	68	9.6	12.1	18.1
Sz 5-6	52	74	9.6	12.1	19.3
Sz 7-8	56	81	9.6	12.1	19.9

Minimum polyethylene thickness for a 9mm metal-backed component is 6.7mm on the medial side.

*Baseplate thickness included.


**Compatible with JOURNEY BCS, JOURNEY II BCS and CR femoral components.



Insert offering/compatibility

	Femoral Size									
Insert Size	1	2	3	4	5	6	7	8	9	10
1-2	●	●	●	●						
3-4		●	●	●	●	●				
5-6				●	●	●	●	●	●	
7-8*						●	●	●	●	●

*Size 7-8 tibia inserts can also be used with size 9 tibia baseplates.





Thickness offerings and tibial slopes

JOURNEY[◇] II tibial inserts

Design	Thickness (mm)							
CR/MD	9	10	11	12	13	15	18	
Deep-Dished	9	10	11	12	13	15	18	21
BCS/AS	9	10	11	12	13	15	18	21
Constrained	9	10	11	12	13	15	18	21 25

Average Slope for JOURNEY II tibial inserts

Design	Medial (X)	Lateral (Y)
CR	0°	5°
MD	8.5°	6.8°
DD	11°	6°
BCS	0°	4°
AS	0°	1.3°
Constrained	0°	0°




Tibia baseplate dimensions (mm)

JOURNEY[◇] Tibia Baseplate

Size	Anterior Posterior	Medial Lateral
	A/P	M/L
1	42	60
2	45	64
3	48	68
4	50	71
5	52	74
6	54	77
7	56	81
8	59	85
9	61	89

Note: Stem sloped 3° posteriorly. Stem length is 50mm on all nonporous sizes





Tibial baseplate dimensions (mm)

CONCELOC[◇] Tibia Baseplate


Size	A/P	M/L
1	42	60
2	45	64
3	48	68
4	50	71
5	52	74
6	54	77
7	56	81
8	59	85

All dimensions in millimeters (mm)



Distal view

Dimensions above are the same for all size tibias and are in millimeters (mm).



Dimensions above are the same for all size tibias and are in millimeters (mm).



Trials account for thickness.



Patellar dimensions (mm)


JOURNEY[◇] Patella

Round	Diameter	Thickness	Peg Height	Peg Diameter	Peg Circle Diameter
Size	A	B	C	D	E
26mm Std	26	9	4.4	5.1	15.3
29mm Std	29	9	4.4	5.1	15.3
32mm Std	32	9	4.4	5.1	15.3
35mm Std	35	9	4.4	5.1	17.9
38mm Std	38	9	4.4	5.1	17.9
41mm Std	41	9	4.4	5.1	17.9



Oval	Diameter	Thickness	Peg Height	Peg Diameter	Peg Circle Diameter	Major Diameter
Size	A	B	C	D	E	F
29mm Std	29	8.5	4.4	5.1	15.3	35
32mm Std	32	9.0	4.4	5.1	15.3	38
35mm Std	35	9.0	4.4	5.1	17.9	41
38mm Std	38	9.5	4.4	5.1	17.9	44
41mm Std	41	10.0	4.4	5.1	17.9	47

Biconvex	Diameter	Thickness	Peg Height	Peg Diameter
Size	A	B	C	D
23mm Std	23	13	4.1	4.7
26mm Std	26	13	4.1	4.7
29mm Std	29	13	3.1	4.7
32mm Std	32	13	3.1	4.7



CONCELOC[◇] Patella options




CONCELOC Patellar Dimensions

Round	Diameter	Thickness	Peg height	Peg diameter	Peg circle diameter	CONCELOC pad diameter
Size	A	B	C	D	E	F
26mm Std	26	9	4.4	6	15.3	22
29mm Std	29	9	4.4	6	15.3	25
32mm Std	32	9	4.4	6	15.3	26
35mm Std	35	9	4.4	6	17.9	29
38mm Std	38	10	4.4	6	17.9	32
41mm Std	41	10	4.4	6	17.9	35

All dimensions in millimeters (mm)

Oval	Diameter	Thickness	Peg height	Peg diameter	Peg circle diameter	CONCELOC pad diameter	Major diameter
Size	A	B	C	D	E	F	G
29mm Std	29	9	4.4	6	15.3	25	35
32mm Std	32	9.5	4.4	6	15.3	26	38
35mm Std	35	9.5	4.4	6	17.9	29	41
38mm Std	38	10.5	4.4	6	17.9	32	44
41mm Std	41	10.5	4.4	6	17.9	35	47

All dimensions in millimeters (mm)




GENESIS[◇] II stems

These stems could be used with a male tapered baseplate during a complex primary case or for a revision case.

- Non-slotted stem
- ◐ Slotted stems


Size	Length (mm)	
	100	150
10	●	●
12	●	●
14	●	●
16	●	◐
18	●	
20	●	
22	●	
24	●	

Note: the length of the stem is measured from the undersurface of the tibial baseplate.




LEGION[◇] stems


120mm stem



160mm stem



220mm stem




Size	Cemented (mm)		
9			
10	80	120	160
11			
12	80	120	160
13			
14	80	120	160
15			
16		120	160
18			
20		120	160
22			
24			




LEGION[◇] stems offset


120mm stem offset



160mm stem offset



220mm stem offset



Size	Press-fit Straight (mm)		
9	120	160	
10	120	160	220
11	120	160	220
12	120	160	220
13	120	160	220
14	120	160	220
15	120	160	220
16	120	160	220
18	120	160	220
20	120	160	220
22	120	160	220
24	120	160	220



Femoral implants



JOURNEY^o II BCS Femoral Component

Size	OXINIUM ^o		CoCr	
	Left	Right	Left	Right
1	74022121	74022111	74024921	74024211
2	74022122	74022112	74024922	74024212
3	74022123	74022113	74024923	74024213
4	74022124	74022114	74024924	74024214
5	74022125	74022115	74024925	74024215
6	74022126	74022116	74024926	74024216
7	74022127	74022117	74024927	74024217
8	74022128	74022118	74024928	74024218
9	74022129	74022119	74024929	74024219
10	74022120	74022110		

JOURNEY II CR Femoral Component

Size	OXINIUM		CoCr	
	Left	Right	Left	Right
1	74021161	74021151	74021261	74021251
2	74021162	74021152	74021262	74021252
3	74021163	74021153	74021263	74021253
4	74021164	74021154	74021264	74021254
5	74021165	74021155	74021265	74021255
6	74021166	74021156	74021266	74021256
7	74021167	74021157	74021267	74021257
8	74021168	74021158	74021268	74021258
9	74021169	74021159	74021269	74021259
10	74021160	74021150		



Inserts



JOURNEY[®] II BCS Insert

Size 1-2	Left	Right
9mm	74027221	74027211
10mm	74027222	74027212
11mm	74027223	74027213
12mm	74027224	74027214
13mm	74027225	74027215
15mm	74027226	74027216
18mm	74027227	74027217
21mm	74027228	74027218

Size 3-4	Left	Right
9mm	74027241	74027231
10mm	74027242	74027232
11mm	74027243	74027233
12mm	74027244	74027234
13mm	74027245	74027235
15mm	74027246	74027236
18mm	74027247	74027237
21mm	74027248	74027238

Size 5-6	Left	Right
9mm	74027261	74027251
10mm	74027262	74027252
11mm	74027263	74027253
12mm	74027264	74027254
13mm	74027265	74027255
15mm	74027266	74027256
18mm	74027267	74027257
21mm	74027268	74027258

Size 7-8*	Left	Right
9mm	74027281	74027271
10mm	74027282	74027272
11mm	74027283	74027273
12mm	74027284	74027274
13mm	74027285	74027275
15mm	74027286	74027276
18mm	74027287	74027277
21mm	74027288	74027278

*Size 7-8 tibia inserts can also be used with size 9 tibia baseplates.



Inserts *continued*



JOURNEY^o II Constrained Insert

Size 1-2	Left	Right
9mm	74029221	74029211
10mm	74029222	74029212
11mm	74029223	74029213
12mm	74029224	74029214
13mm	74029225	74029215
15mm	74029226	74029216
18mm	74029227	74029217
21mm	74029228	74029218
25mm	74029229	74029219

Size 3-4	Left	Right
9mm	74029241	74029231
10mm	74029242	74029232
11mm	74029243	74029233
12mm	74029244	74029234
13mm	74029245	74029235
15mm	74029246	74029236
18mm	74029247	74029237
21mm	74029248	74029238
25mm	74029249	74029239

Size 5-6	Left	Right
9mm	74029261	74029251
10mm	74029262	74029252
11mm	74029263	74029253
12mm	74029264	74029254
13mm	74029265	74029255
15mm	74029266	74029256
18mm	74029267	74029257
21mm	74029268	74029258
25mm	74029269	74029259

Size 7-8*	Left	Right
9mm	74029281	74029271
10mm	74029282	74029272
11mm	74029283	74029273
12mm	74029284	74029274
13mm	74029285	74029275
15mm	74029286	74029276
18mm	74029287	74029277
21mm	74029288	74029278
25mm	74029289	74029279

*Size 7-8 tibia inserts can also be used with size 9 tibia baseplates.



Inserts *continued*



JOURNEY[®] II CR Insert

Size 1-2	Left	Right
9mm	74025621	74025611
10mm	74025622	74025612
11mm	74025623	74025613
12mm	74025624	74025614
13mm	74025625	74025615
15mm	74025626	74025616
18mm	74025627	74025617

Size 3-4	Left	Right
9mm	74025641	74025631
10mm	74025642	74025632
11mm	74025643	74025633
12mm	74025644	74025634
13mm	74025645	74025635
15mm	74025646	74025636
18mm	74025647	74025637

Size 5-6	Left	Right
9mm	74025661	74025651
10mm	74025662	74025652
11mm	74025663	74025653
12mm	74025664	74025654
13mm	74025665	74025655
15mm	74025666	74025656
18mm	74025667	74025657

Size 7-8*	Left	Right
9mm	74025681	74025671
10mm	74025682	74025672
11mm	74025683	74025673
12mm	74025684	74025674
13mm	74025685	74025675
15mm	74025686	74025676
18mm	74025687	74025677

*Size 7-8 tibia inserts can also be used with size 9 tibia baseplates.





JOURNEY^o II MD Insert

Size 1-2	Left	Right
9mm	74021431	74021421
10mm	74021432	74021422
11mm	74021433	74021423
12mm	74021434	74021424
13mm	74021435	74021425
15mm	74021436	74021426
18mm	74021437	74021427

Size 3-4	Left	Right
9mm	74021451	74021441
10mm	74021452	74021442
11mm	74021453	74021443
12mm	74021454	74021444
13mm	74021455	74021445
15mm	74021456	74021446
18mm	74021457	74021447

Size 5-6	Left	Right
9mm	74021471	74021461
10mm	74021472	74021462
11mm	74021473	74021463
12mm	74021474	74021464
13mm	74021475	74021465
15mm	74021476	74021466
18mm	74021477	74021467

Size 7-8*	Left	Right
9mm	74021491	74021481
10mm	74021492	74021482
11mm	74021493	74021483
12mm	74021494	74021484
13mm	74021495	74021485
15mm	74021496	74021486
18mm	74021497	74021487

*Size 7-8 tibia inserts can also be used with size 9 tibia baseplates.





JOURNEY^o II Deep-Dished Insert**

Size 1-2	Left	Right
9mm	74025721	74025711
10mm	74025722	74025712
11mm	74025723	74025713
12mm	74025724	74025714
13mm	74025725	74025715
15mm	74025726	74025716
18mm	74025727	74025717
21mm	74025728	74025718

Size 3-4	Left	Right
9mm	74025741	74025731
10mm	74025742	74025732
11mm	74025743	74025733
12mm	74025744	74025734
13mm	74025745	74025735
15mm	74025746	74025736
18mm	74025747	74025737
21mm	74025748	74025738

** Compatible with JOURNEY BCS, JOURNEY II BCS and CR femoral components.

Size 5-6	Left	Right
9mm	74025761	74025751
10mm	74025762	74025752
11mm	74025763	74025753
12mm	74025764	74025754
13mm	74025765	74025755
15mm	74025766	74025756
18mm	74025767	74025757
21mm	74025768	74025758

Size 7-8*	Left	Right
9mm	74025781	74025771
10mm	74025782	74025772
11mm	74025783	74025773
12mm	74025784	74025774
13mm	74025785	74025775
15mm	74025786	74025776
18mm	74025787	74025777
21mm	74025788	74025778

*Size 7-8 tibia inserts can also be used with size 9 tibia baseplates.



Revision insert options



Note Inserts needed when revising to a LEGION Revision femoral component

LEGION® PS High Flex XLPE Insert with JOURNEY® Lock Detail

Thickness	Size 1-2	Size 3-4
9mm	71425401	71425409
10mm	71425402	71425411
11mm	71425403	71425412
12mm	71425404	71425413
13mm	71425405	71425414
15mm	71425406	71425415
18mm	71425407	71425416
21mm	71425408	71425417

Thickness	Size 5-6	Size 7-8
9mm	71425418	71425427
10mm	71425419	71425428
11mm	71425421	71425429
12mm	71425422	71425431
13mm	71425423	71425432
15mm	71425424	71425433
18mm	71425425	71425434
21mm	71425426	71425435

LEGION PS High Flex Insert Trial with JOURNEY Lock Detail

Thickness	Size 1-2	Size 3-4
9mm	71433101	71433103
10mm	71433102	71433104

Thickness	Size 5-6	Size 7-8
9mm	71433105	71433107
10mm	71433106	71433108



Revision insert options *continued*

LEGION[®] Constrained Insert with JOURNEY[®] Lock Detail

Thickness	Size 1-2	Size 3-4
9mm	71425463	71425472
11mm	71425464	71425473
13mm	71425465	71425474
15mm	71425466	71425475
18mm	71425467	71425476
21mm	71425468	71425477
25mm	71425469	71425478
30mm	71425471	71425479

Thickness	Size 5-6	Size 7-8
9mm	71425481	71425489
11mm	71425482	71425491
13mm	71425483	71425492
15mm	71425484	71425493
18mm	71425485	71425494
21mm	71425486	71425495
25mm	71425487	71425496
30mm	71425488	71425497

LEGION Constrained Insert Trial with JOURNEY Lock Detail

Thickness	Size 1-2	Size 3-4
9mm	71433114	71433115

Thickness	Size 5-6	Size 7-8
9mm	71433116	71433117



Revision insert options *continued*

JOURNEY BCS Revision Insert (Tall Post)**

Thickness	Size 1-2 Left	Size 1-2 Right
9mm	71934506	71934500
10mm	71934507	71934501
11mm	71934508	71934502
13mm	71934509	71934503
15mm	71934510	71934504
18mm	71934511	71934505

Thickness	Size 3-4 Left	Size 3-4 Right
9mm	71934518	71934512
10mm	71934519	71934513
11mm	71934520	71934514
13mm	71934521	71934515
15mm	71934522	71934516
18mm	71934523	71934517

Thickness	Size 5-6 Left	Size 5-6 Right
9mm	71934530	71934524
10mm	71934531	71934525
11mm	71934532	71934526
13mm	71934533	71934527
15mm	71934534	71934528
18mm	71934535	71934529

Thickness	Size 7-8 Left	Size 7-8 Right
9mm	71934542	71934536
10mm	71934543	71934537
11mm	71934544	71934538
13mm	71934545	71934539
15mm	71934546	71934540
18mm	71934547	71934541

**Only for use with first generation JOURNEY BCS femoral components when a poly exchange is needed.



Tibial baseplates



JOURNEY^o Tibia Baseplate

Size	Left	Right
1	74022221	74022211
2	74022222	74022212
3	74022223	74022213
4	74022224	74022214
5	74022225	74022215
6	74022226	74022216
7	74022227	74022217
8	74022228	74022218
9*	71933634	71933635
9 Trial*	71933636	71933637

LEGION^o Revision Baseplate w/ JOURNEY Lock Detail

Size	Left	Right
1*	71933567	71933575
2*	71933568	71933576
3*	71933569	71933577
4*	71933570	71933578
5*	71933571	71933579
6*	71933572	71933580
7*	71933573	71933581
8*	71933574	71933582

*Available through SPG

CONCELOC^o Tibia Baseplate

Size	Left	Right
1	71425321	71425301
2	71425322	71425302
3	71425323	71425303
4	71425324	71425304
5	71425325	71425305
6	71425326	71425306
7	71425327	71425307
8	71425328	71425308



Patellas



JOURNEY^o Patella, Round

Catalog Item	Size
74024826	26mm
74024829	29mm
74024832	32mm
74024835	35mm
74024838	38mm
74024841	41mm

JOURNEY Patella, Biconvex

Catalog Item	Size
74024623	23mm
74024626	26mm
74024629	29mm
74024632	32mm

GENESIS^o II Patella, Biconvex

Catalog Item	Size
71420566	23mm
71420568	26mm
71420570	29mm
71420572	32mm

GENESIS^o II Patella, Round

Catalog Item	Size
71420580	26mm
71420574	29mm
71420576	32mm
71420578	35mm

GENESIS II Patella, Oval

Catalog Item	Size	Thickness
71421029	29mm	8.5mm
71421032	32mm	9.0mm
71421035	35mm	9.0mm
71421038	38mm	9.5mm
71421041	41mm	10mm

GENESIS II Round Patella w/ JOURNEY^o Peg

Catalog Item	Size
71934488*	26
71934489*	29
71934490*	32
71934491*	35
71934492*	38
71934493*	41

GENESIS II Round Patella Trial with JOURNEY Peg

Catalog Item	Size
71934494*	26
71934495*	29
71934496*	32
71934497*	35
71934498*	38
71934499*	41

*Available through SPG



Patellas *continued*

GENESIS[®] II 7.5mm Round Patella

Catalog Item	Size
71932634*	26mm
71932635*	29mm
71932636*	32mm
71932637*	35mm

GENESIS II 7.5mm Round Patella Trial

Catalog Item	Size
71928372*	26mm
71928373*	29mm
71928374*	32mm
71928375*	35mm

JOURNEY[®] 7.5mm Round Patella

Catalog Item	Size
71932910*	26mm
71932911*	29mm
71932912*	32mm
71932913*	35mm

JOURNEY 7.5mm Round Patella Trial

Catalog Item	Size
71932924*	26mm
71932925*	29mm
71932926*	32mm
71932927*	35mm

*Available through SPG



CONCELOC[◇] Patellas



Porous Patella Oval Implants Set #74027420

Catalog Item	Size
74027429	29mm
74027432	32mm
74027435	35mm
74027438	38mm
74027441	41mm

Porous Patella Round Implants Set #74025460

Catalog Item	Size
74025462	26mm
74025463	29mm
74025464	32mm
74025465	35mm
74025466	38mm
74025467	41mm



Tibial wedges

LEGION® 5mm Hemi-Stepped Tibial Wedge

Catalog item	Description	Catalog item	Description
71421141	Sizes 1-2, 34, 5-6, 7-8	71421145	Sizes 1-2, 34, 5-6, 7-8
71421142	LT-Medial/	71421146	LT-Lateral/
71421143	RT-Lateral	71421147	RT-Medial
71421144		71421148	

LEGION 10mm Hemi-Stepped Tibial Wedge

Catalog item	Description	Catalog item	Description
71423033	Sizes 1-2, 3-4, 5-6, 7-8	71423037	Sizes 1-2, 3-4, 5-6, 7-8
71423034	LT-Medial/	71423038	LT-Medial/
71423035	RT-Lateral	71423039	RT-Lateral
71423036		71423040	

LEGION 15mm Hemi-Stepped Tibial Wedge

Catalog item	Description	Catalog item	Description
71423041	Sizes 1-2, 3-4, 5-6, 7-8	71423045	Sizes 1-2, 3-4, 5-6, 7-8
71423042	LT-Medial/	71423046	LT-Medial/
71423043	RT-Lateral	71423047	RT-Lateral
71423044		71423048	

LEGION® Full-Stepped Tibial Wedge

Catalog item	Description	Catalog item	Description
71423057	Sizes 1-2, 3-4, 5-6, 7-8, 10mm	71423061	Sizes 1-2, 3-4, 5-6, 7-8, 15mm
71423058		71423062	
71423059		71423063	
71423060		71423064	

LEGION Screws

Catalog item	Description
71423097	Tibial Wedge Screw



Tibial Stems

GENESIS[®] II Stem

Stem Diameter	100mm Catalog Item	150mm Catalog Item
10mm	71420628	71420647
12mm	71420630	71420649
14mm	71420632	71420648
16mm	71420634	71420650
18mm	71420636	
20mm	71420638	
22mm	71420640	
24mm	71420642	

LEGION[®] Press-Fit Stem

Stem Diameter	120mm Catalog Item	160mm Catalog Item	220mm Catalog Item
9mm	71424022*	71424042*	N/A
10mm	71424023	71424043	71424063
11mm	71424024	71424044	71424064*
12mm	71424025	71424045	71424065
13mm	71424026	71424046	71424066*
14mm	71424027	71424047	71424067
15mm	71424028	71424048	71424068*
16mm	71424029	71424049	71424069
18mm	71424031	71424051	71424071
20mm	71424033	71424053	71424073
22mm	71424035*	71424055	71424075*
24mm	71424037*	71424057*	71424077*

*Not in standard set, limited availability

LEGION[®] Cemented Stem

Stem Diameter	120mm Catalog Item	160mm Catalog Item
10mm	71424182	71424202
12mm	71424184	71424204
14mm	71424186	71424206
16mm	71424188	71424208
18mm	71424190	71424210
20mm	71424192	71424212

LEGION Stem Extensions and Offset Coupler

Catalog Item	Size
71424161	Short Stem Extension 10 x 80
71424163	Short Stem Extension 12 x 80
71424165	Short Stem Extension 14 x 80
71424221*	Offset Coupler Angled
71424223	Offset Coupler + 2mm
71424225	Offset Coupler + 4mm
71424227	Offset Coupler + 6mm

LEGION Male-to-Male Coupler









Catalog Item	Size
71933693	LEGION Male-to-Male Mini Coupler 2mm
71933694	LEGION Male-to-Male Mini Coupler 4mm
71933695	LEGION Male-to-Male Mini Coupler 6mm
71933696	LEGION Male-to-Male Mini Coupler Trial 2mm
71933697	LEGION Male-to-Male Mini Coupler Trial 4mm
71933698	LEGION Male-to-Male Mini Coupler Trial 6mm
7193K178	LEGION Male-to-Male Mini Coupler Set (includes the Implants and trials)

*Not in standard set, limited availability





SYNC Instrumentation

JOURNEY[◇] II Overview



Cat. No.	Description
	SYNC Alignment Options
74039980	SYNC General Instruments
74034240	JOURNEY II Modular Femoral Finishing
74034100	SYNC Patella Resection
74010260	SYNC Patella Reaming Instruments
74011290	SYNC Biconvex Patella Instruments
74010360	JOURNEY II Femoral Trials
74031170	JOURNEY II Small (Sz 1-2) Femoral Outliers
74031220	JOURNEY II Large (Sz 9-10) Femoral Outliers
74010330	Modular Tibia Trial and Prep
74033850	JOURNEY II Insert Trials BCS
74033880	JOURNEY II Insert Trials CR
74031250	JOURNEY II Insert Trials Medial Dished*




*JOURNEY II MD insert trials are not available in the standard set.




SYNC Alignment Options

Note Where noted on the following pages, a femoral alignment set and tibial alignment set will need to be ordered


	Set Number	Set Description	Standard/Alternate	Notes
Femur Options	74039990	Adjustable Femoral Alignment	Standard	Valgus angle is adjustable on the guide
	74016350	Fixed Femoral Alignment	Alternate	Uses a dedicated valgus bushing (4,5,6 or 7°)
	74039940	Kinematic Alignment Instruments	Alternate	
Tibia Options	74016450	Macro Guide Tibial Alignment	Standard	Tibial resection adjustments made with button on alignment tube
	74016360	Micro Extramedullary Tibial Alignment	Alternate	Adds additional fine tune adjustment in cutting guide
Femur/Tibia	74039930	Spiked EM/IM Tibial and Adjustable Femoral Alignment	Alternate	



Adjustable Femur



Fixed Femur



Kinematic Alignment




Macro Tibia



Micro Tibia



Spiked EM Tibia




Intramedullary Tibia



Adjustable Femoral Alignment

Set No. 74039990

Note In order to complete the Alignment Set, one Femoral Alignment and one Tibia Alignment option must be ordered.



Adjustable Femoral Alignment (Standard Option)


Cat. No.	Description	Ref
71674576	TRIGEN ^o T-Handle	A
74016205	Distal Femur Cutting Block	B
74012111	Intramedullary Drill 9.5mm	C
74016211	IM Rod Long	D
74016212	IM Rod Short	E
74016203	Distal Cut Valgus Guide Bridge Adjustable	F
74016210	Variable Angle Valgus Guide	G
74034130	S+N Adjustable Femoral Alignment Module	
74034286	72mm Modular Tray	
74034289	Modular Tray Metal Lid	



Fixed Femoral Alignment

Set No. 74016350

Note In order to complete the Alignment Set, one Femoral Alignment and one Tibia Alignment option must be ordered.



Fixed Femoral Alignment (Alternate Option)


Cat. No.	Description	Ref
71674576	TRIGEN T-Handle	A
74016205	Distal Femur Cutting Block	B
74012111	Intramedullary Drill 9.5mm	C
74016211	IM Rod Long	D
74016212	IM Rod Short	E
74016200	Distal Cut Valgus Guide	F
74016202	Distal Cut Valgus Guide Bridge Fixed	G
74016206	Four Degree Valgus Bushing	H
74016207	Five Degree Valgus Bushing	I
74016208	Six Degree Valgus Bushing	J
74016209	Seven Degree Valgus Bushing	K
74034133	S+N Fixed Femoral Alignment Module	
74034286	72mm Modular Tray	
74034289	Modular Tray Metal Lid	



Kinematic Alignment Instruments

Set No. 74039940

Note In order to complete the Alignment Set, one Femoral Alignment and one Tibia Alignment option must be ordered.




Cat. No.	Description	Ref
71674576	TRIGEN ^o T-Handle	A
42407562	Distal Cartilage Shim - 2mm (x2)	B
42407354	Kinematic Alignment Caliper	C
74012111	Intramedullary Drill 9.5mm	D
74016212	IM Rod Short	E
74016202	Distal Cut Valgus Guide Bridge Fixed	F
42407624	Varus/Valgus Recut Block - 2deg	G
74016205	Distal Femur Cutting Block	H
42407560	LEGION ^o Kinematic Alignment Distal Valgus Guide	I
42407565	JOURNEY ^o Kinematic Alignment Distal Valgus Guide - Left	J
42407566	JOURNEY Kinematic Alignment Distal Valgus Guide - Right	K
71441622	Kinematic Alignment Module	
74034286	72mm Modular Tray	
74034289	Modular Tray Metal Lid	



Macro Guide Tibial Alignment

Set No. 74016450

Note In order to complete the Alignment Set, one Femoral Alignment and one Tibia Alignment option must be ordered.



Macro Guide Tibial Alignment (Standard Option)


Cat. No.	Description	Ref
	Finger Mat	A
74016216	Tibial Cutting Guide Three Deg Left	B
74016217	Tibial Cutting Guide Three Deg Right	C
74016241	Ankle Clamp	D
74016231	EM Posterior Slope Guide	E
74016229	Tibial Stylus Adjustable Footed	F
74016232	EM Tibial Alignment Tube Macro	G
74032502	S+N Macro Guide Tibial Module	



Micro Extramedullary Tibial Module

Set No. 74016360

Note In order to complete the Alignment Set, one Femoral Alignment and one Tibia Alignment option must be ordered.




Micro Extramedullary Tibial Module (Alternate Option)

Cat. No.	Description	Ref
74016237	Tibial Cutting Guide Head Three Degree Left	A
74016238	Tibial Cutting Guide Head Three Degree Right	B
74016241	Ankle Clamp	C
74016230	Tibial Stylus 2 and 10mm Pegged	D
74016236	EM Tibial Alignment Guide Micro Head	E
74016231	EM Posterior Slope Guide	F
74016233	EM Tibial Alignment Tube Micro	G
74032503	S+N Micro Guide Tibial Module	



Spiked EM/IM Tibial and Adjustable Femoral Alignment

Set No. 74039930




Cat. No.	Description	Ref
42407455	Spiked EM Slope Guide	A
74012111	Intramedullary Drill 9.5mm	B
74016229	Tibial Stylus Adjustable Footed	C
74016203	Distal Cut Valgus Guide Bridge Adjustable	D
74034197	Variable Angle Valgus Guide	E
74016211	IM Rod Long	F
74016212	IM Rod Short	G
42407453	Spiked EM Alignment Jig	H
42407454	Spiked EM/IM Alignment Rod	I
42407457	Intramedullary Alignment Jig	J
71674576	TRIGEN [®] T-Handle	K
42407616	Spiked EM/IM 3 Degree Symmetric Tibia Cutting Block	L
42407611	Spiked EM/IM 3 Degree Left Tibia Cutting Block	M
42407612	Spiked EM/IM 3 Degree Right Tibia Cutting Block	N
74016205	Distal Femur Cutting Block	O
74016241	Ankle Clamp	P
71441621	Spiked EM/IM - Adjustable Femur Alignment Module	
74034286	72mm Modular Tray	
74034289	Modular Tray Metal Lid	



SYNC General Instruments

Set No. 74039980




Cat. No.	Description	Ref
71513331	Universal Pin Driver	A
115035	Hex Screwdriver	B
71440302	VISIONAIRE ^o Lightweight Alignment Rod	C
74016213	EM Alignment Guide Handle x2	D
74016254	Extension Spacer Block	E
71440366	GENESIS ^o II Universal Extractor	F
71440303	VISIONAIRE Alignment Connector	G
74016255	Standard Flexion Spacer Block	H
74013723	18mm Spacer Plate Shim	I
74013722	15mm Spacer Plate Shim	J
74013721	13mm Spacer Plate Shim	K
74016258	12mm Spacer Plate Shim	L
74016257	11mm Spacer Plate Shim	M

Cat. No.	Description	Ref
74016242	Modular Impact Handle	N
74036148	Pin Caddy SCD Insert	O
74012825	JOURNEY ^o II Removal Tool	P
71440040	GENESIS II 11mm Tibial Drill	Q
74016461	MIS Pin Driver	R
74034196	Insert Trial Pitchfork	S
74018911	JOURNEY Articular Insert Assembly Tool	T
71440491	Universal Pin Puller	U
74016271	Resection Check	V
74016249	Modular Universal Tibial Impactor	W
74034136	SYNC General Instruments Module	
74034286	72mm Modular Tray	
74034289	Modular Tray Metal Lid	



JOURNEY[◇] II SYNC Modular Femoral Finishing Instruments

Set No. 74034240





Cat. No.	Description	Ref
74012457	JOURNEY II TKA Femoral Sizing Stylus	A
74012455	JOURNEY II TKA Femoral Sizing Guide Left	B
74012456	JOURNEY II TKA Femoral Sizing Guide Right	C
74012413	JOURNEY DCF AP Femoral Cutting Block Size 3	D
74012414	JOURNEY DCF AP Femoral Cutting Block Size 4	E
74012415	JOURNEY DCF AP Femoral Cutting Block Size 5	F
74012416	JOURNEY DCF AP Femoral Cutting Block Size 6	G
74012417	JOURNEY DCF AP Femoral Cutting Block Size 7	H
74012418	JOURNEY DCF AP Femoral Cutting Block Size 8	I
74016248	Modular Box Chisel	J
74016269	PS Box Reamer	K
74018903	JOURNEY II CR Locking Femoral Impactor Bumper	L
74011821	JOURNEY Femoral Implant Impactor Bumper Right	M
74012821	JOURNEY Femoral Implant Impactor Bumper Left	N

Cat. No.	Description	Ref
74012575	JOURNEY II BCS Femoral Box Prep Guide Size 3-5	O
74012576	JOURNEY II BCS Femoral Box Prep Guide Size 6-8	P
74018902	Outer Grip Locking Femoral Implant Impactor	Q
74012514	JOURNEY II TKA Femoral Trial Impactor Size 1-10	R
74016252	Modular JOURNEY Femoral Impactor	S
74011855	JOURNEY II CR Femoral Lug Drill	T
74012421	JOURNEY AP Cutting Block Impactor	U
74012812	JOURNEY II BCS Locking Femoral Implant Impactor	V
74034141	S+N JOURNEY II Modular Femoral Finishing Module	
74034290	107mm Modular Tray	
74034289	Modular Tray Metal Lid	



SYNC Patella Resection Instruments

Set No. 74034100



Saw Guide + 7.5mm Trials

Cat. No.	Description	Ref
42404211	ENZTEC Classic Patella Saw Guide 8-10mm	A
74035026	Cemented 7.5 Round Patella Trial w/ JRNY Pegs 26mm*	B
74035029	Cemented 7.5 Round Patella Trial w/ JRNY Pegs 29mm*	C
74035032	Cemented 7.5 Round Patella Trial w/ JRNY Pegs 32mm*	D
74035035	Cemented 7.5 Round Patella Trial w/ JRNY Pegs 35mm*	E
74034166	S+N Non Porous Saw Guide And 7.5mm Trial Module	

*7.5mm patellas may not be available in all regions

Sizing and Trials

Cat. No.	Description	Ref
74015462	Universal Caliper	A
42404321	ENZTEC Premium Patella Clamp Resection	B
74015461	Round/Oval Sizing Guides	C
74014748	Universal Patella Peg Drill	D
42404360	ENZTEC Patella Swivel Cement Clamp	E
42404366	ENZTEC Patella Drill Guide Round Oval 41mm	F
42404365	ENZTEC Patella Drill Guide Round Oval 38mm	G
42404364	ENZTEC Patella Drill Guide Round Oval 35mm	H
42404363	ENZTEC Patella Drill Guide Round Oval 32mm	I
42404362	ENZTEC Patella Drill Guide Round Oval 29mm	J
42404361	ENZTEC Patella Drill Guide Round 26mm	K

Note: Cemented Round Patella Trials with JRNY Pegs are compatible with both JOURNEY Round Patella implants and the new Cemented Round Patella implants

Drill Guides and Drill


Cat. No.	Description	Ref
74037026	Cemented 9.0 Round Patella Trial w/ JRNY Pegs 26mm	L
74037029	Cemented 9.0 Round Patella Trial w/ JRNY Pegs 29mm	M
74037032	Cemented 9.0 Round Patella Trial w/ JRNY Pegs 32mm	N
74037035	Cemented 9.0 Round Patella Trial w/ JRNY Pegs 35mm	O
74037038	Cemented 9.0 Round Patella Trial w/ JRNY Pegs 38mm	P
74037041	Cemented 9.0 Round Patella Trial w/ JRNY Pegs 41mm	Q
71430429	GENESIS ^o II Oval Patella Trial 29mm	R
71430432	GENESIS II Oval Patella Trial 32mm	S
71430435	GENESIS II Oval Patella Trial 35mm	T
71430438	GENESIS II Oval Patella Trial 38mm	U
71430441	GENESIS II Oval Patella Trial 41mm	V
74034165	S+N Non-Porous Resection Patella Module	
74034286	72mm Modular Tray	
74034289	Modular Tray Metal Lid	



SYNC Patella Reaming Instruments

Set No. 74010260

Note SYNC Patella Resection Set (74034100) must also be ordered to complete Patella Reaming Instrument Set.




Cat. No.	Description	Ref
71440427	Large Patella Depth Stop	A
42404322	ENZTEC Premium Patella Clamp – Reaming	B
42404370	ENZTEC Patella Reamer Collet - S+N – 26	C
42404371	ENZTEC Patella Reamer Collet - S+N – 29	D
42404372	ENZTEC Patella Reamer Collet - S+N – 32	E
42404373	ENZTEC Patella Reamer Collet - S+N – 35	F
42404374	ENZTEC Patella Reamer Collet - S+N – 38	G
42404375	ENZTEC Patella Reamer Collet - S+N – 41	H
42404376	ENZTEC Patella Reamer Collet - S+N – 44	I
42404377	ENZTEC Patella Reamer Collet - S+N – 46	J
42404378	ENZTEC Patella Reamer Collet - S+N – 51	K
42407380	ENZTEC Patella - Depth Gage - S+N	L

Cat. No.	Description	Ref
74016726	Patella Resurface Reamer 26mm	M
74016729	Patella Resurface Reamer 29mm	N
74016732	Patella Resurface Reamer 32mm	O
74016735	Patella Resurface Reamer 35mm	P
74016738	Patella Resurface Reamer 38mm	Q
74016741	Patella Resurface Reamer 41mm	R
74016744	Patella Resurface Reamer 44mm	S
74016746	Patella Resurface Reamer 46mm	T
74016751	Patella Resurface Reamer 51mm	U
74034167	S+N Reaming Patella Module	
74034287	90mm Modular Tray	
74034289	Modular Tray Metal Lid	



SYNC Biconvex Patella Instruments

Set No. 74011290




Cat. No.	Description	Ref
42404322	ENZTEC Premium Patella Clamp - Reaming	A
74015462	Universal Caliper	B
42404369	ENZTEC Patella Reamer Collet - S+N - 23	C
42404370	ENZTEC Patella Reamer Collet - S+N - 26	D
42404371	ENZTEC Patella Reamer Collet - S+N - 29	E
42404372	ENZTEC Patella Reamer Collet - S+N - 32	F
71440326	GENESIS [◇] II Patellar Depth Stop	G
42407380	ENZTEC Patella - Depth Gage - S+N	H
42404360	Patella Swivel Cement Clamp - S+N	I
74018723	Patella Biconvex Reamer 23mm	J
74018726	Patella Biconvex Reamer 26mm	K
74018729	Patella Biconvex Reamer 29mm	L
74018732	Patella Biconvex Reamer 32mm	M
71430566	GENESIS II Biconvex Trial 23mm	N
71430568	GENESIS II Biconvex Trial 26mm	O

Cat. No.	Description	Ref
71430570	GENESIS II Biconvex Trial 29mm	P
71430572	GENESIS II Biconvex Trial 32mm	Q
74034171	S+N Biconvex Patella Module	
74034286	72mm Modular Tray	
74034289	Modular Tray Metal Lid	




JOURNEY[◇] II Femoral Trials

Set No. 74010360



Left

Cat. No.	Description	Ref
74031593	JOURNEY II TKA Knee Femoral Trial Size 3 Left	A
74031224	JOURNEY II TKA Knee Femoral Trial Size 4 Left	B
74031225	JOURNEY II TKA Knee Femoral Trial Size 5 Left	C
74031226	JOURNEY II TKA Knee Femoral Trial Size 6 Left	D
74031227	JOURNEY II TKA Knee Femoral Trial Size 7 Left	E
74031228	JOURNEY II TKA Knee Femoral Trial Size 8 Left	F
74032143	JOURNEY II BCS Cam Trial Size 3 Left	G
74032144	JOURNEY II BCS Cam Trial Size 4 Left	H
74032145	JOURNEY II BCS Cam Trial Size 5 Left	I
74032146	JOURNEY II BCS Cam Trial Size 6 Left	J
74032147	JOURNEY II BCS Cam Trial Size 7 Left	K
74032148	JOURNEY II BCS Cam Trial Size 8 Left	L
74031363	JOURNEY II CR Femoral Notch Trial Size 3 Left	M
74031364	JOURNEY II CR Femoral Notch Trial Size 4 Left	N
74031365	JOURNEY II CR Femoral Notch Trial Size 5 Left	O
74031366	JOURNEY II CR Femoral Notch Trial Size 6 Left	P
74031367	JOURNEY II CR Femoral Notch Trial Size 7 Left	Q
74031368	JOURNEY II CR Femoral Notch Trial Size 8 Left	R
74032508	S+N JOURNEY II TKA Left Femur Trial Module	




Right

Cat. No.	Description	Ref
74031353	JOURNEY II CR Femoral Notch Trial Size 3 Right	A
74031354	JOURNEY II CR Femoral Notch Trial Size 4 Right	B
74031355	JOURNEY II CR Femoral Notch Trial Size 5 Right	C
74031356	JOURNEY II CR Femoral Notch Trial t Size 6 Right	D
74031357	JOURNEY II CR Femoral Notch Trial Size 7 Right	E
74031358	JOURNEY II CR Femoral Notch Trial Size 8 Right	F
74032133	JOURNEY II BCS Cam Trial Size 3 Right	G
74032134	JOURNEY II BCS Cam Trial Size 4 Right	H
74032135	JOURNEY II BCS Cam Trial Size 5 Right	I
74032136	JOURNEY II BCS Cam Trial Size 6 Right	J
74032137	JOURNEY II BCS Cam Trial Size 7 Right	K
74032138	JOURNEY II BCS Cam Trial Size 8 Right	L
74031583	JOURNEY II TKA Knee Femoral Trial Size 3 Right	M
74031214	JOURNEY II TKA Knee Femoral Trial Size 4 Right	N
74031215	JOURNEY II TKA Knee Femoral Trial Size 5 Right	O
74031216	JOURNEY II TKA Knee Femoral Trial Size 6 Right	P
74031217	JOURNEY II TKA Knee Femoral Trial Size 7 Right	Q
74031218	JOURNEY II TKA Knee Femoral Trial Size 8 Right	R
74032509	S+N JOURNEY II TKA Right Femur Trial Module	

Cat. No.	Description
74034290	107mm Modular Tray
74034289	Modular Tray Metal Lid




JOURNEY[◇] II Femoral Outliers



Small (Sz 1-2)
Set No. 74031170

Cat. No.	Description	Ref
74031161	JOURNEY II Total Knee Femoral Trial Left Size 1	A
74031222	JOURNEY II Total Knee Femoral Trial Left Size 2	B
74031151	JOURNEY II Total Knee Femoral Trial Right Size 1	C
74031212	JOURNEY II Total Knee Femoral Trial Right Size 2	D
74011441	JOURNEY DCF AP Femoral Cutting Block Size 1	E
74011442	JOURNEY DCF AP Femoral Cutting Block Size 2	F
74032131	JOURNEY II BCS Cam Trial Size 1 Right	G
74032132	JOURNEY II BCS Cam Trial Size 2 Right	H
74032141	JOURNEY II BCS Cam Trial Size 1 Left	I
74032142	JOURNEY II BCS Cam Trial Size 2 Left	J
74031351	JOURNEY II CR Femoral Notch Trial Right Size 1	K
74031352	JOURNEY II CR Femoral Notch Trial Right Size 2	L
74031361	JOURNEY II CR Femoral Notch Trial Left Size 1	M
74031362	JOURNEY II CR Femoral Notch Trial Left Size 2	N
74012574	JOURNEY II BCS Femoral Box Prep Guide Size 1-2	O
74032512	S+N JOURNEY II TKA Small Femoral Outlier Module	
74034176	S+N Half Outer Tray	
74034177	S+N Half Tray Lid	




Large (Sz 9-10)
Set No. 74031220

Cat. No.	Description	Ref
74012577	JOURNEY II BCS Femoral Box Prep Guide Size 9-10	A
74031359	JOURNEY II CR Femoral Notch Trial Right Size 9	B
74031360	JOURNEY II CR Femoral Notch Trial Right Size 10	C
74031369	JOURNEY II CR Femoral Notch Trial Left Size 9	D
74031370	JOURNEY II CR Femoral Notch Trial Left Size 10	E
74032139	JOURNEY II BCS Cam Trial Size 9 Right	F
74032130	JOURNEY II BCS Cam Trial Size 10 Right	G
74032149	JOURNEY II BCS Cam Trial Size 9 Left	H
74032140	JOURNEY II BCS Cam Trial Size 10 Left	I
74012419	JOURNEY DCF AP Femoral Cutting Block Size 9	J
74012410	JOURNEY DCF AP Femoral Cutting Block Size 10	K
74031150	JOURNEY II TKA Femoral Trial Right Size 10	L
74031159	JOURNEY II TKA Femoral Trial Right Size 9	M
74031160	JOURNEY II TKA Femoral Trial Left Size 10	N
74031169	JOURNEY II TKA Femoral Trial Left Size 9	O
74032513	S+N JOURNEY II TKA Large Femoral Outlier Module	
74034176	S+N Half Outer Tray	
74034177	S+N Half Tray Lid	



Modular Tibial Trial and Prep Instrument

Set No. 74010330




Cat. No.	Description	Ref	Cat. No.	Description	Ref	Cat. No.	Description	Ref
74033801	Universal Tibial Base Trial Size 1 Left	A	74016246	Size 5-6 Non-Porous Modular Fin Punch	S	74033655	Universal Insert Spacer Size 5-6 15mm	AJ
74033802	Universal Tibial Base Trial Size 2 Left	B	74016247	Size 7-8 Non-Porous Modular Fin Punch	T	74033654	Universal Insert Spacer Size 5-6 13mm	AK
74033803	Universal Tibial Base Trial Size 3 Left	C	71442196	Tibia Plug*	U	74016265	Universal Insert Spacer Size 5-6 12mm	AL
74033804	Universal Tibial Base Trial Size 4 Left	D	74033617	Universal Insert Spacer Size 1-2 21mm	V	74016264	Universal Insert Spacer Size 5-6 11mm	AM
74033805	Universal Tibial Base Trial Size 5 Left	E	74033616	Universal Insert Spacer Size 1-2 18mm	W	74033677	Universal Insert Spacer Size 7-8 21mm	AN
74033806	Universal Tibial Base Trial Size 6 Left	F	74033615	Universal Insert Spacer Size 1-2 15mm	X	74033676	Universal Insert Spacer Size 7-8 18mm	AO
74033807	Universal Tibial Base Trial Size 7 Left	G	74033614	Universal Insert Spacer Size 1-2 13mm	Y	74033675	Universal Insert Spacer Size 7-8 15mm	AP
74033808	Universal Tibial Base Trial Size 8 Left	H	74016261	Universal Insert Spacer Size 1-2 12mm	Z	74033674	Universal Insert Spacer Size 7-8 13mm	AQ
74033811	Universal Tibial Base Trial Size 1 Right	I	74016260	Universal Insert Spacer Size 1-2 11mm	AA	74016267	Universal Insert Spacer Size 7-8 12mm	AR
74033812	Universal Tibial Base Trial Size 2 Right	J	74033637	Universal Insert Spacer Size 3-4 21mm	AB	74016266	Universal Insert Spacer Size 7-8 11mm	AS
74033813	Universal Tibial Base Trial Size 3 Right	K	74033636	Universal Insert Spacer Size 3-4 18mm	AC	74034153	S+N Modular Tibial Trial And Prep Module	
74033814	Universal Tibial Base Trial Size 4 Right	L	74033635	Universal Insert Spacer Size 3-4 15mm	AD	74034287	90mm Modular Tray	
74033815	Universal Tibial Base Trial Size 5 Right	M	74033634	Universal Insert Spacer Size 3-4 13mm	AE	74034289	Modular Tray Metal Lid	
74033816	Universal Tibial Base Trial Size 6 Right	N	74016263	Universal Insert Spacer Size 3-4 12mm	AF			
74033817	Universal Tibial Base Trial Size 7 Right	O	74016262	Universal Insert Spacer Size 3-4 11mm	AG			
74033818	Universal Tibial Base Trial Size 8 Right	P	74033657	Universal Insert Spacer Size 5-6 21mm	AH			
74016244	Size 1-2 Non-Porous Modular Fin Punch	Q	74033656	Universal Insert Spacer Size 5-6 18mm	AI			
74016245	Size 3-4 Non-Porous Modular Fin Punch	R						

* Not included in tray




JOURNEY[◇] II Insert Trials: BCS

Set No. 74033850



Left

Cat. No.	Description	Ref
74033831	JOURNEY II BCS Insert Trial Size 1-2 9mm Left	A
74033832	JOURNEY II BCS Insert Trial Size 1-2 10mm Left	B
74033833	JOURNEY II BCS Insert Trial Size 3-4 9mm Left	C
74033834	JOURNEY II BCS Insert Trial Size 3-4 10mm Left	D
74033835	JOURNEY II BCS Insert Trial Size 5-6 9mm Left	E
74033836	JOURNEY II BCS Insert Trial Size 5-6 10mm Left	F
74033837	JOURNEY II BCS Insert Trial Size 7-8 9mm Left	G
74033838	JOURNEY II BCS Insert Trial Size 7-8 10mm Left	H
74033871	JOURNEY II Constrained Insert Trial Size 1-2 9mm Left	I
74033872	JOURNEY II Constrained Insert Trial Size 1-2 10mm Left	J
74033873	JOURNEY II Constrained Insert Trial Size 3-4 9mm Left	K
74033874	JOURNEY II Constrained Insert Trial Size 3-4 10mm Left	L
74033875	JOURNEY II Constrained Insert Trial Size 5-6 9mm Left	M
74033876	JOURNEY II Constrained Insert Trial Size 5-6 10mm Left	N
74033877	JOURNEY II Constrained Insert Trial Size 7-8 9mm Left	O
74033878	JOURNEY II Constrained Insert Trial Size 7-8 10mm Left	P
74034178	S+N JOURNEY II Insert Trials Module	




Right

Cat. No.	Description	Ref
74033839	JOURNEY II BCS Insert Trial Size 1-2 9mm Right	A
74033841	JOURNEY II BCS Insert Trial Size 1-2 10mm Right	B
74033842	JOURNEY II BCS Insert Trial Size 3-4 9mm Right	C
74033843	JOURNEY II BCS Insert Trial Size 3-4 10mm Right	D
74033844	JOURNEY II BCS Insert Trial Size 5-6 9mm Right	E
74033845	JOURNEY II BCS Insert Trial Size 5-6 10mm Right	F
74033846	JOURNEY II BCS Insert Trial Size 7-8 9mm Right	G
74033847	JOURNEY II BCS Insert Trial Size 7-8 10mm Right	H
74033881	JOURNEY II Constrained Insert Trial Size 1-2 9mm Right	I
74033882	JOURNEY II Constrained Insert Trial Size 1-2 10mm Right	J
74033883	JOURNEY II Constrained Insert Trial Size 3-4 9mm Right	K
74033884	JOURNEY II Constrained Insert Trial Size 3-4 10mm Right	L
74033885	JOURNEY II Constrained Insert Trial Size 5-6 9mm Right	M
74033886	JOURNEY II Constrained Insert Trial Size 5-6 10mm Right	N
74033887	JOURNEY II Constrained Insert Trial Size 7-8 9mm Right	O
74033888	JOURNEY II Constrained Insert Trial Size 7-8 10mm Right	P
74034178	S+N JOURNEY II Insert Trials Module	
74034290	107mm Modular Tray	
74034289	Modular Tray Metal Lid	




JOURNEY[◇] II Insert Trials: CR

Set No. 74033880



Left

Cat. No.	Description	Ref
74033851	JOURNEY II CR Insert Trial Size 1-2 9mm Left	A
74033852	JOURNEY II CR Insert Trial Size 1-2 10mm Left	B
74033853	JOURNEY II CR Insert Trial Size 3-4 9mm Left	C
74033854	JOURNEY II CR Insert Trial Size 3-4 10mm Left	D
74033855	JOURNEY II CR Insert Trial Size 5-6 9mm Left	E
74033856	JOURNEY II CR Insert Trial Size 5-6 10mm Left	F
74033857	JOURNEY II CR Insert Trial Size 7-8 9mm Left	G
74033858	JOURNEY II CR Insert Trial Size 7-8 10mm Left	H
74033891	JOURNEY II Deep Dish Insert Trial Size 1-2 9mm Left	I
74033892	JOURNEY II Deep Dish Insert Trial Size 1-2 10mm Left	J
74033893	JOURNEY II Deep Dish Insert Trial Size 3-4 9mm Left	K
74033894	JOURNEY II Deep Dish Insert Trial Size 3-4 10mm Left	L
74033895	JOURNEY II Deep Dish Insert Trial Size 5-6 9mm Left	M
74033896	JOURNEY II Deep Dish Insert Trial Size 5-6 10mm Left	N
74033897	JOURNEY II Deep Dish Insert Trial Size 7-8 9mm Left	O
74033898	JOURNEY II Deep Dish Insert Trial Size 7-8 10mm Left	P
74034178	S+N JOURNEY II Insert Trials Module	




Right

Cat. No.	Description	Ref
74033859	JOURNEY II CR Insert Trial Size 1-2 9mm Right	A
74033861	JOURNEY II CR Insert Trial Size 1-2 10mm Right	B
74033862	JOURNEY II CR Insert Trial Size 3-4 9mm Right	C
74033863	JOURNEY II CR Insert Trial Size 3-4 10mm Right	D
74033864	JOURNEY II CR Insert Trial Size 5-6 9mm Right	E
74033865	JOURNEY II CR Insert Trial Size 5-6 10mm Right	F
74033866	JOURNEY II CR Insert Trial Size 7-8 9mm Right	G
74033867	JOURNEY II CR Insert Trial Size 7-8 10mm Right	H
74033901	JOURNEY II Deep Dish Insert Trial Size 1-2 9mm Right	I
74033902	JOURNEY II Deep Dish Insert Trial Size 1-2 10mm Right	J
74033903	JOURNEY II Deep Dish Insert Trial Size 3-4 9mm Right	K
74033904	JOURNEY II Deep Dish Insert Trial Size 3-4 10mm Right	L
74033905	JOURNEY II Deep Dish Insert Trial Size 5-6 9mm Right	M
74033906	JOURNEY II Deep Dish Insert Trial Size 5-6 10mm Right	N
74033907	JOURNEY II Deep Dish Insert Trial Size 7-8 9mm Right	O
74033908	JOURNEY II Deep Dish Insert Trial Size 7-8 10mm Right	P
74034178	S+N JOURNEY II Insert Trials Module	
74034290	107mm Modular Tray	
74034289	Modular Tray Metal Lid	



JOURNEY[◇] II Insert Trials: Medial Dished

Set No. 74031250



Cat. No.	Description	Ref
74031253	JOURNEY II Medial Dished Insert Trial Size 1-2 Left 9mm	A
74031254	JOURNEY II Medial Dished Insert Trial Size 1-2 Left 10mm	B
74031257	JOURNEY II Medial Dished Insert Trial Size 3-4 Left 9mm	C
74031258	JOURNEY II Medial Dished Insert Trial Size 3-4 Left 10mm	D
74031263	JOURNEY II Medial Dished Insert Trial Size 5-6 Left 9mm	E
74031264	JOURNEY II Medial Dished Insert Trial Size 5-6 Left 10mm	F
74031267	JOURNEY II Medial Dished Insert Trial Size 7-8 Left 9mm	G
74031268	JOURNEY II Medial Dished Insert Trial Size 7-8 Left 10mm	H
74031251	JOURNEY II Medial Dished Insert Trial Size 1-2 Right 9mm	I
74031252	JOURNEY II Medial Dished Insert Trial Size 1-2 Right 10mm	J
74031255	JOURNEY II Medial Dished Insert Trial Size 3-4 Right 9mm	K
74031256	JOURNEY II Medial Dished Insert Trial Size 3-4 Right 10mm	L
74031261	JOURNEY II Medial Dished Insert Trial Size 5-6 Right 9mm	M
74031262	JOURNEY II Medial Dished Insert Trial Size 5-6 Right 10mm	N
74031265	JOURNEY II Medial Dished Insert Trial Size 7-8 Right 9mm	O
74031266	JOURNEY II Medial Dished Insert Trial Size 7-8 Right 10mm	P
74034178	S+N JOURNEY II Insert Trials Module	
74034176	S+N Half Outer Tray	
74034177	S+N Half Tray Lid	



Caution: US Federal law restricts these devices to the sale to or on the order of a physician.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

Smith & Nephew, Inc.
1450 Brooks Road
Memphis, TN 38116
USA

Telephone: 1-901-396-2121
Information: 1-800-821-5700
Orders/Inquiries: 1-800-238-7538

www.smith-nephew.com

°Trademark of Smith+Nephew
©2025 Smith+Nephew
00828 V9 02/25