

+ Precision
powered
performance

Smith+Nephew



WEREWOLF⁺
FASTSEAL 6.0
Hemostasis Wand

Introducing the new WEREWOLF[®] FASTSEAL 6.0 Hemostasis Wand. Use in your orthopedic procedures for hemostasis of soft tissue and bone. Feel the performance with:

- Fast priming and initiation features^{1,2}
- Operating temperatures equal to or less than 100 degrees Celsius based on testing a single device in triplicate
- A combination of direct and indirect thermal RF energy

Experience the power of the WEREWOLF FASTSEAL 6.0 Hemostasis Wand using the new WEREWOLF+ System, available to use for all your Sports Medicine, ENT and now your orthopedic procedures.

Your existing WEREWOLF Systems can be upgraded to the new WEREWOLF+ System version. Contact your US Sales Representative to find out how.



Ordering information

WEREWOLF FASTSEAL 6.0	
Reference	Description
72290146	WEREWOLF+ Generator
72290042	WEREWOLF FASTSEAL 6.0 Hemostasis Wand



Competitively priced for value. To find out more, contact either your Smith+Nephew Sports Medicine sales representative or your Orthopaedic Reconstruction sales representative.

References

1. Smith+Nephew 2021.Design Verification Report: Challenge Condition, Expected Use & Every Tissue Every Level. Internal Report. 110127-01 Rev A. 2. Smith+Nephew 2021. Priming time associated with the Aquamantys Pump Generator. Internal Report. EO.SPM.PCS.037.001.v1. 3. Marulanda GA, Ulrich SD, Seyler TM, Delanois RE, Mont MA. Reductions in blood loss with a bipolar sealer in total hip arthroplasty. Expert Rev Med Devices. 2008;5(2):125-131. 4. Derman PB, Kamath AF, Lee GC. Saline-coupled bipolar sealing in revision total knee arthroplasty for infection. Am J Orthop (Belle Mead NJ). 2013;42(9):407-411. 5. Smith+Nephew 2020.Report, Peak Temperature Comparative Study, FASTSEAL 6.0 and Aquamantys 6.0. Internal Report. 110180-01 Rev A.

Learn more at smith-nephew.com

Distributed by:
Sports Medicine
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810

www.smith-nephew.com
T +978 749 1000
US Customer Service:
+1 800 343 5717

Manufactured by:
ArthroCare Corporation
7000 West William
Cannon Drive
Austin, TX 78735, USA

[®]Trademark of Smith+Nephew.
All trademarks acknowledged.
©2021 Smith+Nephew. All rights reserved. Printed in USA. 29430 V2 06/21