



+ Redefining peak performance

Consistent deployment¹⁻³, fixation strength^{4,5}
An integral part of Instability Excellence

Smith+Nephew

Q-FIX[◇]
All-Suture Anchor

Q-FIX[◇] All-Suture Anchor

Consistent deployment¹⁻³

- Designed to provide consistent deployment due to the radially expanding implant design and delivery system
- The active deployment against the inserter shaft and the tensioning of the sutures to at least 140N contributes to the low displacement during cyclic loading^{*4-8}


*As demonstrated in benchtop testing



Fixation strength^{4,5}

- Shown to provide greater fixation strength when compared to other commercially available all-suture anchors^{**4,5}
- Provides the benefits of a small, soft anchor with the fixation characteristics similar to or greater than traditional solid body anchors^{***5,9,10}

Max Load of Labral Repair Anchors Tested (in 30 pcf)⁵



**As demonstrated in benchtop testing; Compared to BIORAPTOR[®] 2.3 PK Suture Anchor; Zimmer Biomet JuggerKnot[™] 1.45 (#1 and #2), Conmed Y-Knot[™] Flex 1.3 and Stryker Iconix 1

*** As demonstrated in benchtop testing



Improved access¹¹

- Improve access to challenging shoulder anatomy with the Q-FIX Curved Drill Guide, complete with intuitive visual and tactile cues¹¹
- Greater accuracy than leading competitor curved drill guides*, producing less than 1 degree of deviation¹²

*Measured as drill deviation vs. drill guide angle; Compared to Cayenne SureLock[™] Curved, Stryker Iconix Curved, Arthrex Fibretak Curved, Zimmer Biomet JuggerKnot[™] Long Curved (without centering sleeve)


Low displacement

- Ultra-low displacement during cyclic loading^{**4-8}
- Q-FIX 1.8mm anchor shown to provide leading performance, with the lowest displacement during cyclic loading among competitor all-suture anchors^{***4,5,8}

** As demonstrated in benchtop testing

*** As demonstrated in benchtop testing; Compared to Stryker Iconix 1, 2, 2.5 and 3, Parcus Medical Draw Tight[™] 1.8 and 3.2, Zimmer Biomet Juggerknot[™] 1.4, 1.45 (#1 and #2), 1.5 and 2.9 and Conmed Y-Knot[™] Flex 1.3, 1.8 and 2.8

Displacement at 400 cycles (in 30 pcf)⁵



Ordering information

1.8mm Q-FIX [®] MINI All-Suture Anchor	
Reference #	Description
Anchors	
72290123	1.8mm Q-FIX MINI All-Suture Anchor
Reusable instruments	
25-1812	1.8mm Reusable PATHFINDER [®] Obturator
Disposable kits	
72290125	1.8mm Q-FIX MINI Disposable Kit
72290126	1.8mm Q-FIX MINI XL Disposable Kit

1.8mm Q-FIX All-Suture Anchor	
Reference #	Description
Anchors	
25-1800	1.8mm Q-FIX All-Suture Anchor
Reusable instruments	
72290119	1.8mm Reusable Curved Drill Guide, XL
72205224	1.8mm Reusable Curved Drill Guide
25-1812	1.8mm Reusable PATHFINDER Obturator
72290032	1.8mm Q-FIX Drill Guide, reusable
Disposable kits	
25-1810	1.8mm Q-FIX Disposable Kit
25-1811	1.8mm Q-FIX Disposable Kit XL
Disposable accessories	
72290118	1.8mm Disposable Flex Drill
72290030	1.8mm Q-FIX Straight Drill, disposable

Learn more at smith-nephew.com

Sports Medicine
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810

www.smith-nephew.com
T +978 749 1000
US Customer Service:
+1 800 343 5717

[®]Trademark of Smith+Nephew.
All trademarks acknowledged.
©2021 Smith+Nephew. All rights reserved. Printed in USA. 06770 V4 06/21



References

1. ArthroCare Corporation 2016. Report: Design Verification, 1.8mm & 2.8 mm Q-Fix Suture Anchor P/N 49190-01 Rev.B.
2. ArthroCare Corporation 2017. Report: Design Verification, 1.8mm Q-FIX MINI Soft Suture Anchor P/N 49190-03 Rev. B.
3. Smith+Nephew 2020. Safety and Performance of Q-Fix All-Suture Anchor System 17-5010-11.
4. Barber FA, Herbert MA. All-Suture Anchors: Biomechanical Analysis of Pullout Strength, Displacement, and Failure Mode. *Arthroscopy*. 2017;33(6):1113-1121.
5. Douglass NP, Behn AW, Safran MR. Cyclic and Load to Failure Properties of All-Suture Anchors in Synthetic Acetabular and Glenoid Cancellous Bone. *Arthroscopy*. 2017;33(5):977-985 e975.
6. ArthroCare 2019. Comparative Testing of Bone Anchor Devices, 1.8mm Q-FIX Mini Soft Suture Anchor P/N 49193-02 Rev.B.
7. Nagra NS, Zargar N, Smith RD, Carr AJ. Mechanical properties of all-suture anchors for rotator cuff repair. *Bone Joint Res*. 2017;6(2):82-89.
8. Ruder JA, Dickinson EY, Peindl RD, Habet NA, Trofa DP, Fleischli JE. Cyclic and Load-to-Failure Properties of All-Suture Anchors in Human Cadaveric Shoulder Glenoid Bone. *Arthroscopy*. 2019;35(7):1954-1959 e1954.
9. Bernardoni E, Frank RM, Veera SS, et al. Biomechanical Analysis of All-Suture Anchor Fixation for Rotator Cuff Repair. *Orthop J Sports Med*. 2018;6(7 suppl4).
10. Saper MG, Meijer K, Winnier S, Popovich J, Jr., Andrews JR, Roth C. Biomechanical Evaluation of Classic Solid and All-Soft Suture Anchors for Medial Patellofemoral Ligament Reconstruction. *Am J Sports Med*. 2017;45(7):1622-1626.
11. ArthroCare Corporation 2017. Report: Design Validation, Q-FIX Reusable Curved Drill Guides and Disposable Flex Drill, Irvine CA P/N 92017-01 Rev. A.
12. Smith+Nephew 2017. Q-FIX Family Drill Deviation 15007031 A.