

Smith+Nephew

JOURNEY[◊] II UK
Unicompartmental Knee System

System Specification Guide




JOURNEY[◇] II UK Femoral Dimensions (mm)



Size	AP	SI
1	40	28.8
2	42	30.6
3	44	32.3
4	46	34.0
5	48	35.9
6	50	37.7
7	52	39.6
8	54	41.3
9	56	43.2
10	58	45.2

JOURNEY[◇] II UK Tibia Baseplate Dimensions (mm)



Size	AP	ML	
		Medial	Lateral
0	38	-	23.0
1	40	22.3	24.3
2	42	23.7	25.7
3	44	25.0	27.0
4	46	26.0	28.3
5	48	27.0	29.7
6	50	28.0	31.0
7	52	29.0	32.3
8	54	30.0	-
9	56	31.0	-
10	58	32.0	-

JOURNEY[◇] II UK Femoral Components

Offered in OXINIUM[◇] Technology

Size	LM/RL	RM/LL
1	74026021	74026031
2	74026022	74026032
3	74026023	74026033
4	74026024	74026034
5	74026025	74026035
6	74026026	74026036
7	74026027	74026037
8	74026028	74026038
9	74026029	74026039
10	74026030	74026040

JOURNEY^o II UK Medial Tibia Baseplates

Size	Left*	Right*
1	74026711	74026721
2	74026712	74026722
3	74026713	74026723
4	74026714	74026724
5	74026715	74026725
6	74026716	74026726
7	74026717	74026727
8	74026718	74026728
9	74026719	74026729
10	74026720	74026730

*Medial Baseplates have a secondary use in the contralateral compartment. For example, a Left Medial Baseplate can be used in the Right Lateral compartment and a Right Medial Baseplate can be used in the Left Lateral compartment.

JOURNEY[◇] II UK Lateral Tibia Baseplates

Size	Left*	Right*
0	74026130	74026140
1	74026131	74026141
2	74026132	74026142
3	74026133	74026143
4	74026134	74026144
5	74026135	74026145
6	74026136	74026146
7	74026137	74026147

*Lateral Baseplates have a secondary use in the contralateral compartment. For example, a Left Lateral Baseplate can be used in the Right Medial compartment and a Right Lateral Baseplate can be used in the Left Medial compartment.

JOURNEY^o II UK Medial Inserts

Thickness	Size 1-2	Size 3-4	Size 5-6
8mm	74026148	74026158	74026168
9mm	74026149	74026159	74026169
10mm	74026150	74026160	74026170
11mm	74026151	74026161	74026171
12mm	74026152	74026162	74026172
13mm	74026153	74026163	74026173
14mm*	74026154	74026164	74026174

Thickness	Size 7-8	Size 9-10
8mm	74026178	74026188
9mm	74026179	74026189
10mm	74026180	74026190
11mm	74026181	74026191
12mm	74026182	74026192
13mm	74026183	74026193
14mm*	74026184	74026194

*Special order only

JOURNEY[◇] II UK Lateral Inserts

Thickness	Size 0-1	Size 2-3
8mm	74026248	74026258
9mm	74026249	74026259
10mm	74026250	74026260
11mm	74026251	74026261
12mm	74026252	74026262
13mm	74026253	74026263
14mm*	74026254	74026264

Thickness	Size 4-5	Size 6-7
8mm	74026268	74026278
9mm	74026269	74026279
10mm	74026270	74026280
11mm	74026271	74026281
12mm	74026272	74026282
13mm	74026273	74026283
14mm*	74026274	74026284

*Special order only

Femoral Inlier Sets

74026910 JOURNEY^o II UK OXINIUM^o Core (Size 3-7) Set

74026023 JOURNEY II UK OXINIUM Femoral
3 LM/RL

74026024 JOURNEY II UK OXINIUM Femoral
4 LM/RL

74026025 JOURNEY II UK OXINIUM Femoral
5 LM/RL

74026026 JOURNEY II UK OXINIUM Femoral
6 LM/RL

74026027 JOURNEY II UK OXINIUM Femoral
7 LM/RL

74026033 JOURNEY II UK OXINIUM Femoral
3 RM/LL

74026034 JOURNEY II UK OXINIUM Femoral
4 RM/LL

74026035 JOURNEY II UK OXINIUM Femoral
5 RM/LL

74026036 JOURNEY II UK OXINIUM Femoral
6 RM/LL

74026037 JOURNEY II UK OXINIUM Femoral
7 RM/LL

Femoral Outlier Sets

74026920 JOURNEY[◊] II UK OXINIUM[◊] **Size 1 Set**

74026021	JOURNEY II UK OXINIUM Femoral 1 LM/RL
----------	--

74026031	JOURNEY II UK OXINIUM Femoral 1 RM/LL
----------	--

74026840 JOURNEY II UK OXINIUM **Size 2 Set**

74026022	JOURNEY II UK OXINIUM Femoral 2 LM/RL
----------	--

74026032	JOURNEY II UK OXINIUM Femoral 2 RM/LL
----------	--

74026850 JOURNEY II UK OXINIUM **Size 8 Set**

74026028	JOURNEY II UK OXINIUM Femoral 8 LM/RL
----------	--

74026038	JOURNEY II UK OXINIUM Femoral 8 RM/LL
----------	--

Femoral Outlier Sets

74026860 JOURNEY[◊] II UK OXINIUM[◊] Size 9-10 Set

74026029	JOURNEY II UK OXINIUM Femoral 9 LM/RL
----------	--

74026030	JOURNEY II UK OXINIUM Femoral 10 LM/RL
----------	---

74026039	JOURNEY II UK OXINIUM Femoral 9 RM/LL
----------	--

74026040	JOURNEY II UK OXINIUM Femoral 10 RM/LL
----------	---

Medial Insert Inlier Sets

74026870 JOURNEY^o II UK Medial Insert Core (All sizes; 8-10 Thick) Set

74026148 JOURNEY II UK Medial XLPE Insert
Size 1-2 8mm

74026149 JOURNEY II UK Medial XLPE Insert
Size 1-2 9mm

74026150 JOURNEY II UK Medial XLPE Insert
Size 1-2 10mm

74026158 JOURNEY II UK Medial XLPE Insert
Size 3-4 8mm

74026159 JOURNEY II UK Medial XLPE Insert
Size 3-4 9mm

74026160 JOURNEY II UK Medial XLPE Insert
Size 3-4 10mm

74026168 JOURNEY II UK Medial XLPE Insert
Size 5-6 8mm

74026169 JOURNEY II UK Medial XLPE Insert
Size 5-6 9mm

74026170 JOURNEY II UK Medial XLPE Insert
Size 5-6 10mm

74026178 JOURNEY II UK Medial XLPE Insert
Size 7-8 8mm

74026179 JOURNEY II UK Medial XLPE Insert
Size 7-8 9mm

74026180 JOURNEY II UK Medial XLPE Insert
Size 7-8 10mm

Medial Insert Inlier Sets

**74026870 JOURNEY^o II UK Medial Insert
Core (All sizes; 8-10 Thick) Set**
(continued)

74026188	JOURNEY II UK Medial XLPE Insert Size 9-10 8mm
----------	---

74026189	JOURNEY II UK Medial XLPE Insert Size 9-10 9mm
----------	---

74026190	JOURNEY II UK Medial XLPE Insert Size 9-10 10mm
----------	--

Medial Insert Outlier Sets

74026880 JOURNEY^o II UK Medial Insert Outlier (All sizes; 11-13 Thick) Set

74026151 JOURNEY II UK Medial XLPE Insert
Size 1-2 11mm

74026152 JOURNEY II UK Medial XLPE Insert
Size 1-2 12mm

74026153 JOURNEY II UK Medial XLPE Insert
Size 1-2 13mm

74026161 JOURNEY II UK Medial XLPE Insert
Size 3-4 11mm

74026162 JOURNEY II UK Medial XLPE Insert
Size 3-4 12mm

74026163 JOURNEY II UK Medial XLPE Insert
Size 3-4 13mm

74026171 JOURNEY II UK Medial XLPE Insert
Size 5-6 11mm

74026172 JOURNEY II UK Medial XLPE Insert
Size 5-6 12mm

74026173 JOURNEY II UK Medial XLPE Insert
Size 5-6 13mm

74026181 JOURNEY II UK Medial XLPE Insert
Size 7-8 11mm

74026182 JOURNEY II UK Medial XLPE Insert
Size 7-8 12mm

74026183 JOURNEY II UK Medial XLPE Insert
Size 7-8 13mm

Medial Insert Outlier Sets

**74026880 JOURNEY^o II UK Medial Insert
Outlier (All sizes; 11-13 Thick)
Set** *(continued)*

74026191 JOURNEY II UK Medial XLPE Insert
Size 9-10 11mm

74026192 JOURNEY II UK Medial XLPE Insert
Size 9-10 12mm

74026193 JOURNEY II UK Medial XLPE Insert
Size 9-10 13mm

Lateral Insert Inlier Sets

74026890 JOURNEY^o II UK Lateral Insert Core(All sizes; 8-10 Thick) Set

74026248 JOURNEY II UK Lateral XLPE Insert
0-1 8mm

74026249 JOURNEY II UK Lateral XLPE Insert
0-1 9mm

74026250 JOURNEY II UK Lateral XLPE Insert
0-1 10mm

74026258 JOURNEY II UK Lateral XLPE Insert
2-3 8mm

74026259 JOURNEY II UK Lateral XLPE Insert
2-3 9mm

74026260 JOURNEY II UK Lateral XLPE Insert
2-3 10mm

74026268 JOURNEY II UK Lateral XLPE Insert
4-5 8mm

74026269 JOURNEY II UK Lateral XLPE Insert
4-5 9mm

74026270 JOURNEY II UK Lateral XLPE Insert
4-5 10mm

74026278 JOURNEY II UK Lateral XLPE Insert
6-7 8mm

74026279 JOURNEY II UK Lateral XLPE Insert
6-7 9mm

74026280 JOURNEY II UK Lateral XLPE Insert
6-7 10mm

Lateral Insert Outlier Sets

74026899 JOURNEY^o II UK Lateral Insert Outlier (All sizes; 11-13 Thick) Set

74026251 JOURNEY II UK Lateral XLPE Insert
0-1 11mm

74026252 JOURNEY II UK Lateral XLPE Insert
0-1 12mm

74026253 JOURNEY II UK Lateral XLPE Insert
0-1 13mm

74026261 JOURNEY II UK Lateral XLPE Insert
2-3 11mm

74026262 JOURNEY II UK Lateral XLPE Insert
2-3 12mm

74026263 JOURNEY II UK Lateral XLPE Insert
2-3 13mm

74026271 JOURNEY II UK Lateral XLPE Insert
4-5 11mm

74026272 JOURNEY II UK Lateral XLPE Insert
4-5 12mm

74026273 JOURNEY II UK Lateral XLPE Insert
4-5 13mm

74026281 JOURNEY II UK Lateral XLPE Insert
6-7 11mm

74026282 JOURNEY II UK Lateral XLPE Insert
6-7 12mm

74026283 JOURNEY II UK Lateral XLPE Insert
6-7 13mm

Medial Tibia Inlier Sets

74026710 JOURNEY^o II UK Medial Tibia Core (Sizes 3-8) Set

74026713 JOURNEY II UK Medial Tibia Size 3 Left

74026714 JOURNEY II UK Medial Tibia Size 4 Left

74026715 JOURNEY II UK Medial Tibia Size 5 Left

74026716 JOURNEY II UK Medial Tibia Size 6 Left

74026717 JOURNEY II UK Medial Tibia Size 7 Left

74026718 JOURNEY II UK Medial Tibia Size 8 Left

74026723 JOURNEY II UK Medial Tibia Size 3 Right

74026724 JOURNEY II UK Medial Tibia Size 4 Right

74026725 JOURNEY II UK Medial Tibia Size 5 Right

74026726 JOURNEY II UK Medial Tibia Size 6 Right

74026727 JOURNEY II UK Medial Tibia Size 7 Right

74026728 JOURNEY II UK Medial Tibia Size 8 Right

Medial Tibia Outlier Sets

74026799 JOURNEY^o II UK Medial Tibia Size 1 Set

74026711 JOURNEY II UK Medial Tibia Size
1 Left

74026721 JOURNEY II UK Medial Tibia Size
1 Right

74026930 JOURNEY II UK Medial Tibia Size 2 Set

74026712 JOURNEY II UK Medial Tibia Size
2 Left

74026722 JOURNEY II UK Medial Tibia Size
2 Right

74026940 JOURNEY II UK Medial Tibia Sizes 9-10 Set

74026719 JOURNEY II UK Medial Tibia Size
9 Left

74026729 JOURNEY II UK Medial Tibia Size
9 Right

74026720 JOURNEY II UK Medial Tibia Size
10 Left

74026730 JOURNEY II UK Medial Tibia Size
10 Right

Lateral Tibia Inlier Sets

74026760 JOURNEY[®] II UK Lateral Tibia Core (Sizes 2-7) Set

74026132 JOURNEY II UK Lateral Tibia
Baseplate 2 Left

74026142 JOURNEY II UK Lateral Tibia
Baseplate 2 Right

74026133 JOURNEY II UK Lateral Tibia
Baseplate 3 Left

74026143 JOURNEY II UK Lateral Tibia
Baseplate 3 Right

74026134 JOURNEY II UK Lateral Tibia
Baseplate 4 Left

74026144 JOURNEY II UK Lateral Tibia
Baseplate 4 Right

74026135 JOURNEY II UK Lateral Tibia
Baseplate 5 Left

74026145 JOURNEY II UK Lateral Tibia
Baseplate 5 Right

74026136 JOURNEY II UK Lateral Tibia
Baseplate 6 Left

74026146 JOURNEY II UK Lateral Tibia
Baseplate 6 Right

74026137 JOURNEY II UK Lateral Tibia
Baseplate 7 Left

74026147 JOURNEY II UK Lateral Tibia
Baseplate 7 Right

Lateral Tibia Outlier Sets

74026820 JOURNEY^o II UK Lateral Tibia Sizes 0-1 Set

74026130 JOURNEY II UK Lateral Tibia
Baseplate 0 Left

74026140 JOURNEY II UK Lateral Tibia
Baseplate 0 Right

74026131 JOURNEY II UK Lateral Tibia
Baseplate 1 Left

74026141 JOURNEY II UK Lateral Tibia
Baseplate 1 Right

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets.

Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area.

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

www.smith-nephew.com

Smith & Nephew, Inc.
1450 Brooks Road
Memphis, TN 38116
USA

Telephone: 1-901-396-2121
Information: 1-800-821-5700
Orders and Inquiries: 1-800-238-7538

© 2023 Smith+Nephew

°Trademark of Smith+Nephew

17502 V4 10/23