

Smith+Nephew

JOURNEY[◊] II UK
Unicompartmental Knee System

System Specification Guide




JOURNEY[◇] II UK Femoral Dimensions



Size	AP	SI
1	40	28.8
2	42	30.6
3	44	32.3
4	46	34.0
5	48	35.9
6	50	37.7
7	52	39.6
8	54	41.3
9	56	43.2
10	58	45.2

JOURNEY[◇] II UK Tibia Baseplate Dimensions



Size	AP	ML	
		Medial	Lateral
0	38		23
1	40	22.3	24.3
2	42	23.7	25.7
3	44	25	27
4	46	26	28.3
5	48	27	29.7
6	50	28	31
7	52	29	32.3
8	54	30	
9	56	31	
10	58	32	

JOURNEY[◇] II UK Femoral Components

Offered in OXINIUM[◇], Oxidized Zirconium

Size	LM/RL	RM/LL
1	7402-6021	7402-6031
2	7402-6022	7402-6032
3	7402-6023	7402-6033
4	7402-6024	7402-6034
5	7402-6025	7402-6035
6	7402-6026	7402-6036
7	7402-6027	7402-6037
8	7402-6028	7402-6038
9	7402-6029	7402-6039
10	7402-6030	7402-6040

JOURNEY[®] II UK Medial Tibia Baseplates

Size	Left*	Right*
1	7402-6711	7402-6721
2	7402-6712	7402-6722
3	7402-6713	7402-6723
4	7402-6714	7402-6724
5	7402-6715	7402-6725
6	7402-6716	7402-6726
7	7402-6717	7402-6727
8	7402-6718	7402-6728
9	7402-6719	7402-6729
10	7402-6720	7402-6730

*Medial baseplates have a secondary use in the contralateral compartment. For example, a left medial baseplate can be used in the right lateral compartment and a right medial baseplate can be used in the left lateral compartment.

JOURNEY[◇] II UK Lateral Tibia Baseplates

Size	Left*	Right*
0	7402-6130	7402-6140
1	7402-6131	7402-6141
2	7402-6132	7402-6142
3	7402-6133	7402-6143
4	7402-6134	7402-6144
5	7402-6135	7402-6145
6	7402-6136	7402-6146
7	7402-6137	7402-6147

*Lateral baseplates have a secondary use in the contralateral compartment. For example, a left lateral baseplate can be used in the right medial compartment and a right lateral baseplate can be used in the left medial compartment.

JOURNEY^o II UK Medial Inserts

Thickness	Size 1-2	Size 3-4	Size 5-6
8mm	7402-6148	7402-6158	7402-6168
9mm	7402-6149	7402-6159	7402-6169
10mm	7402-6150	7402-6160	7402-6170
11mm	7402-6151	7402-6161	7402-6171
12mm	7402-6152	7402-6162	7402-6172
13mm	7402-6153	7402-6163	7402-6173
14mm*	7402-6154	7402-6164	7402-6174

Thickness	Size 7-8	Size 9-10
8mm	7402-6178	7402-6188
9mm	7402-6179	7402-6189
10mm	7402-6180	7402-6190
11mm	7402-6181	7402-6191
12mm	7402-6182	7402-6192
13mm	7402-6183	7402-6193
14mm*	7402-6184	7402-6194

*Special order only

JOURNEY[◇] II UK Lateral Inserts

Thickness	Size 0-1	Size 2-3
8mm	7402-6248	7402-6258
9mm	7402-6249	7402-6259
10mm	7402-6250	7402-6260
11mm	7402-6251	7402-6261
12mm	7402-6252	7402-6262
13mm	7402-6253	7402-6263
14mm*	7402-6254	7402-6264

Thickness	Size 4-5	Size 6-7
8mm	7402-6268	7402-6278
9mm	7402-6269	7402-6279
10mm	7402-6270	7402-6280
11mm	7402-6271	7402-6281
12mm	7402-6272	7402-6282
13mm	7402-6273	7402-6283
14mm*	7402-6274	7402-6284

*Special order only

Loaner Instrument Kits

008991 JOURNEY^o II UK Gen Instr

008993 JOURNEY II UK 14mm Inst

008994 JOURNEY II UK LT Medial Inst

008995 JOURNEY II UK RT Medial Inst

009102 JOURNEY II UK Lat Tib Instr

Loaner Implant Kits

009854 JOURNEY II UK Med Tib Left Implants

009855 JOURNEY II UK Med Tib
Right Implants

009856 JOURNEY II UK OX Fem Left Implants

009857 JOURNEY II UK OX Fem
Right Implants

009858 JOURNEY II UK Lat Tib Left Implants

009859 JOURNEY II UK Lat Tib Right Implants

009860 JOURNEY II UK Med Insert Implants

009861 JOURNEY II UK Lat Insert Implants

Inlier Sets

74026710 JOURNEY[◊] II UK Reduced Med Tibia Implant Set

74026713 JOURNEY II UK Medial Tibia
Baseplate 3 Left

74026723 JOURNEY II UK Medial Tibia
Baseplate 3 Right

74026714 JOURNEY II UK Medial Tibia
Baseplate 4 Left

74026724 JOURNEY II UK Medial Tibia
Baseplate 4 Right

74026715 JOURNEY II UK Medial Tibia
Baseplate 5 Left

74026725 JOURNEY II UK Medial Tibia
Baseplate 5 Right

74026716 JOURNEY II UK Medial Tibia
Baseplate 6 Left

74026726 JOURNEY II UK Medial Tibia
Baseplate 6 Right

74026717 JOURNEY II UK Medial Tibia
Baseplate 7 Left

74026727 JOURNEY II UK Medial Tibia
Baseplate 7 Right

74026718 JOURNEY II UK Medial Tibia
Baseplate 8 Left

74026728 JOURNEY II UK Medial Tibia
Baseplate 8 Right

Inlier Sets

74026740 JOURNEY[®] II UK Reduced Femur Resec. Ox Implant Set

74026023 JOURNEY II UK Resec OX Femur
3 LM/RL

74026024 JOURNEY II UK Resec OX Femur
4 LM/RL

74026025 JOURNEY II UK Resec OX Femur
5 LM/RL

74026026 JOURNEY II UK Resec OX Femur
6 LM/RL

74026027 JOURNEY II UK Resec OX Femur
7 LM/RL

74026028 JOURNEY II UK Resec OX Femur
8 LM/RL

74026033 JOURNEY II UK Resec OX Femur
3 RM/LL

74026034 JOURNEY II UK Resec OX Femur
4 RM/LL

74026035 JOURNEY II UK Resec OX Femur
5 RM/LL

74026036 JOURNEY II UK Resec OX Femur
6 RM/LL

74026037 JOURNEY II UK Resec OX Femur
7 RM/LL

74026038 JOURNEY II UK Resec OX Femur
8 RM/LL

Inlier Sets

74026750 JOURNEY[®] II UK Reduced Medial Insert Implant Set

74026158 JOURNEY II UK XLPE Tibia
Articular Med 3-4 8mm

74026159 JOURNEY II UK XLPE Tibia
Articular Med 3-4 9mm

74026160 JOURNEY II UK XLPE Tibia
Articular Med 3-4 10mm

74026161 JOURNEY II UK XLPE Tibia
Articular Med 3-4 11mm

74026168 JOURNEY II UK XLPE Tibia
Articular Med 5-6 8mm

74026169 JOURNEY II UK XLPE Tibia
Articular Med 5-6 9mm

74026170 JOURNEY II UK XLPE Tibia
Articular Med 5-6 10mm

74026171 JOURNEY II UK XLPE Tibia
Articular Med 5-6 11mm

74026178 JOURNEY II UK XLPE Tibia
Articular Med 7-8 8mm

74026179 JOURNEY II UK XLPE Tibia
Articular Med 7-8 9mm

74026180 JOURNEY II UK XLPE Tibia
Articular Med 7-8 10mm

74026181 JOURNEY II UK XLPE Tibia
Articular Med 7-8 11mm

Inlier Sets

74026760 JOURNEY[®] II UK Reduced Lateral Tibia Implant Set

74026132 JOURNEY II UK Lateral Tibia Baseplate 2 Left

74026142 JOURNEY II UK Lateral Tibia Baseplate 2 Right

74026133 JOURNEY II UK Lateral Tibia Baseplate 3 Left

74026143 JOURNEY II UK Lateral Tibia Baseplate 3 Right

74026134 JOURNEY II UK Lateral Tibia Baseplate 4 Left

74026144 JOURNEY II UK Lateral Tibia Baseplate 4 Right

74026135 JOURNEY II UK Lateral Tibia Baseplate 5 Left

74026145 JOURNEY II UK Lateral Tibia Baseplate 5 Right

74026136 JOURNEY II UK Lateral Tibia Baseplate 6 Left

74026146 JOURNEY II UK Lateral Tibia Baseplate 6 Right

74026137 JOURNEY II UK Lateral Tibia Baseplate 7 Left

74026147 JOURNEY II UK Lateral Tibia Baseplate 7 Right

Inlier Sets

74026255 JOURNEY[®] II UK Reduced Lateral Insert Set

74026258 JOURNEY[®] II UK Lat XLPE Insert
2-3 8mm

74026259 JOURNEY II UK Lat XLPE Insert
2-3 9mm

74026260 JOURNEY II UK Lat XLPE Insert
2-3 10mm

74026261 JOURNEY II UK Lat XLPE Insert
2-3 11mm

74026268 JOURNEY II UK Lat XLPE Insert
4-5 8mm

74026269 JOURNEY II UK Lat XLPE Insert
4-5 9mm

74026270 JOURNEY II UK Lat XLPE Insert
4-5 10mm

74026271 JOURNEY II UK Lat XLPE Insert
4-5 11mm

74026278 JOURNEY II UK Lat XLPE Insert
6-7 8mm

74026279 JOURNEY II UK Lat XLPE Insert
6-7 9mm

74026280 JOURNEY II UK Lat XLPE Insert
6-7 10mm

74026281 JOURNEY II UK Lat XLPE Insert
6-7 11mm

Outlier Sets

74026790 JOURNEY[®] II UK Medial Tibia Implants

74026711 JOURNEY II UK Medial Tibia
Baseplate 1 Left

74026721 JOURNEY II UK Medial Tibia
Baseplate 1 Right

74026712 JOURNEY II UK Medial Tibia
Baseplate 2 Left

74026722 JOURNEY II UK Medial Tibia
Baseplate 2 Right

74026719 JOURNEY II UK Medial Tibia
Baseplate 9 Left

74026729 JOURNEY II UK Medial Tibia
Baseplate 9 Right

74026720 JOURNEY II UK Medial Tibia
Baseplate 10 Left

74026730 JOURNEY II UK Medial Tibia
Baseplate 10 Right

Outlier Sets

74026800 JOURNEY[®] II UK Femur Resec. OX Implants

74026021	JOURNEY II UK Resec OX Femur 1 LM/RL
----------	---

74026022	JOURNEY II UK Resec OX Femur 2 LM/RL
----------	---

74026029	JOURNEY II UK Resec OX Femur 9 LM/RL
----------	---

74026030	JOURNEY II UK Resec OX Femur 10 LM/RL
----------	--

74026031	JOURNEY II UK Resec OX Femur 1 RM/LL
----------	---

74026032	JOURNEY II UK Resec OX Femur 2 RM/LL
----------	---

74026039	JOURNEY II UK Resec OX Femur 9 RM/LL
----------	---

74026040	JOURNEY II UK Resec OX Femur 10 RM/LL
----------	--

Outlier Sets

74026810 JOURNEY[®] II UK Medial Insert Implants

74026148 JOURNEY II UK XLPE Tibia
Articular Med 1-2 8mm

74026149 JOURNEY II UK XLPE Tibia
Articular Med 1-2 9mm

74026150 JOURNEY II UK XLPE Tibia
Articular Med 1-2 10mm

74026151 JOURNEY II UK XLPE Tibia
Articular Med 1-2 11mm

74026152 JOURNEY II UK XLPE Tibia
Articular Med 1-2 12mm

74026153 JOURNEY II UK XLPE Tibia
Articular Med 1-2 13mm

74026188 JOURNEY II UK XLPE Tibia
Articular Med 9-10 8mm

74026189 JOURNEY II UK XLPE Tibia
Articular Med 9-10 9mm

74026190 JOURNEY II UK XLPE Tibia
Articular Med 9-10 10mm

74026191 JOURNEY II UK XLPE Tibia
Articular Med 9-10 11mm

74026192 JOURNEY II UK XLPE Tibia
Articular Med 9-10 12mm

Outlier Sets

74026810 JOURNEY^o II UK Medial Insert Implants (continued)

74026193 JOURNEY II UK XLPE Tibia
Articular Med 9-10 13mm

74026162 JOURNEY II UK XLPE Tibia
Articular Med 3-4 12mm

74026163 JOURNEY II UK XLPE Tibia
Articular Med 3-4 13mm

74026172 JOURNEY II UK XLPE Tibia
Articular Med 5-6 12mm

74026173 JOURNEY II UK XLPE Tibia
Articular Med 5-6 13mm

74026182 JOURNEY II UK XLPE Tibia
Articular Med 7-8 12mm

74026183 JOURNEY II UK XLPE Tibia
Articular Med 7-8 13mm

74026820 JOURNEY II UK Lateral Tibia Implants

74026130 JOURNEY II UK Lateral Tibia
Baseplate 0 Left

74026140 JOURNEY II UK Lateral Tibia
Baseplate 0 Right

74026131 JOURNEY II UK Lateral Tibia
Baseplate 1 Left

74026141 JOURNEY II UK Lateral Tibia
Baseplate 1 Right

Outlier Sets

74026830 JOURNEY[°] II UK Lateral Inserts

74026248 JOURNEY[°] II UK Lat XLPE Insert
0-1 8mm

74026249 JOURNEY II UK Lat XLPE Insert
0-1 9mm

74026250 JOURNEY II UK Lat XLPE Insert
0-1 10mm

74026251 JOURNEY II UK Lat XLPE Insert
0-1 11mm

74026252 JOURNEY II UK Lat XLPE Insert
0-1 12mm

74026253 JOURNEY II UK Lat XLPE Insert
0-1 13mm

74026262 JOURNEY II UK Lat XLPE Insert
2-3 12mm

74026263 JOURNEY II UK Lat XLPE Insert
2-3 13mm

74026272 JOURNEY II UK Lat XLPE Insert
4-5 12mm

74026273 JOURNEY II UK Lat XLPE Insert
4-5 13mm

74026282 JOURNEY II UK Lat XLPE Insert
6-7 12mm

74026283 JOURNEY II UK Lat XLPE Insert
6-7 13mm

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

www.smith-nephew.com

Smith & Nephew, Inc.
1450 Brooks Road
Memphis, TN 38116
USA

Telephone: 1-901-396-2121
Information: 1-800-821-5700
Orders and Inquiries: 1-800-238-7538

© 2020 Smith+Nephew

°Trademark of Smith+Nephew

17502 V3 12/20